

Datasheet: 12004158

| | |
|----------------------|--|
| Description: | GOAT ANTI-MOUSE IgG StarBright™ Blue 700 |
| Specificity: | IgG |
| Format: | StarBright Blue 700 |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 400 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | |
| Western Blotting | ▪ | | | 1/2,500 - 1/5,000 |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | | | |
|-----------------------|---|----------------------------|--------------------------|
| Target Species | Mouse | | |
| Product Form | Purified IgG conjugated to StarBright Blue 700 - lyophilized | | |
| Reconstitution | Resuspend contents of the tube in the indicated volume of distilled or deionized water and leave on ice for at least 30 min prior to use. Brief centrifugation (pulse spin for 2-3 sec at max speed in a tabletop microcentrifuge) may be used to collect the contents at the bottom of the tube. Do not centrifuge more than 10 sec. | | |
| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
| | StarBright Blue 700 | 470 | 700 |

Antiserum Preparation Antiserum to Mouse IgG was raised by repeated immunisation of Goats with highly purified antigen.

Buffer Solution Phosphate buffered saline

Preservative >80% Sucrose
Stabilisers <1% Polyethylene Glycol

<0.04% 2-Methyl-2H-Isothiazol-3-One
<10% Bovine Serum Albumin

Product Information StarBright Blue 700 Secondary Antibodies are ideal for fluorescent western blotting — either for the detection of a single target protein or for multiplex detection of several proteins on one blot, without stripping and reprobing. The StarBright Fluorophore is composed of a condensed polymer made up of multiple light-absorbing and -emitting monomers which provides an exceptionally bright signal (high quantum yield) compared to most traditional fluorophores.

StarBright Blue 700 Fluorescent Secondary Antibodies can be used with traditional fluorophores like the DyLight 800 Fluorophore for multiplexing. In addition, StarBright Antibodies can be used with Bio-Rad's stain-free technology and/or the hFAB™ Rhodamine Anti-Housekeeping Primary Antibodies for protein normalization. These antibodies are optimized for use with the ChemiDoc™ MP Imaging System, permitting detection of up to three proteins in a single blot.

Features and Benefits

Easy multiplexing — StarBright Secondary Antibodies can be used with IRDye 800 or other common fluorophores to detect multiple proteins on the same blot

Short exposure time — the exceptional brightness of the StarBright Fluorophore leads to shorter exposure times compared to traditional fluorescent dyes: seconds vs. minutes

High signal-to-noise ratio — StarBright Secondary Antibodies emit in the red/far red, where there is minimal background ($Ex_{max}/Em_{max} = 470\text{ nm}/700\text{ nm}$)

Low nonspecific binding — StarBright Fluorophore is conjugated to highly cross-adsorbed antibodies

Goat Anti-Mouse IgG StarBright Blue 700 has been cross-adsorbed against bovine, goat, human, rabbit and rat.

Western Blotting StarBright Blue 700 Secondary Antibodies can be used for multiplex fluorescent western blotting with the right combination of fluorophores. These fluorescent secondary antibodies can also be used to detect post-translational modifications like phosphorylation of a target protein.

Note: StarBright Antibodies are not suitable for stripping and reprobing.

Instructions For Use Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/12004157.pdf

References

1. Dunbar, K. *et al.* (2021) IMiDs induce FAM83F degradation via an interaction with CK1 α to attenuate Wnt signalling. [Life Sci Alliance. 4 \(2\): e202000804.](#)
2. Dunbar, K. *et al.* (2020) FAM83F regulates canonical Wnt signalling through an interaction with CK1 α . [Life Sci Alliance.4 \(2\): e202000805.](#)
3. Holoch, D. *et al.* (2021) A cis-acting mechanism mediates transcriptional memory at Polycomb target genes in mammals. [Nat Genet. 53 \(12\): 1686-97.](#)
4. Thomas, S. *et al.* (2021) Altered expression of fragile X mental retardation-1 (FMR1) in the thymus in autoimmune myasthenia gravis. [J Neuroinflammation. 18 \(1\): 270.](#)

5. Han, H. *et al.* (2022) Loss of aryl hydrocarbon receptor suppresses the response of colonic epithelial cells to IL22 signaling by upregulating SOCS3. [Am J Physiol Gastrointest Liver Physiol. 322 \(1\): G93-G106.](#)
6. Meissner, M.E. *et al.* (2021) Differential Activity of APOBEC3F, APOBEC3G, and APOBEC3H in the Restriction of HIV-2. [J Mol Biol. 434 \(2\): 167355.](#)
7. , G.D.W. *et al.* (2022) Acquired chemoresistance can lead to increased resistance of pancreatic cancer cells to oncolytic vesicular stomatitis virus [Molecular Therapy - Oncolytics. 24: 59-76.](#)
8. Baumann, C. *et al.* (2021) Changes in chromatin accessibility landscape and histone H3 core acetylation during valproic acid-induced differentiation of embryonic stem cells. [Epigenetics Chromatin. 14 \(1\): 58.](#)
9. Riera-Tur, I. *et al.* (2022) Amyloid-like aggregating proteins cause lysosomal defects in neurons via gain-of-function toxicity. [Life Sci Alliance. 5 \(3\): e202101185.](#)
10. Cejas, R.B. *et al.* (2022) Impact of DYRK1A Expression on TNNT2 Splicing and Daunorubicin Toxicity in Human iPSC-Derived Cardiomyocytes. [Cardiovasc Toxicol. 22 \(8\): 701-712.](#)
11. Lach, S.R. *et al.* (2023) Conformational rearrangements in the sensory RcsF/OMP complex mediate signal transduction across the bacterial cell envelope. [PLoS Genet. 19 \(1\): e1010601.](#)
12. Chae, S.A. *et al.* (2024) Exercise enhances placental labyrinth trophoblast development by activation of PGC-1 α and FNDC5/irisin†. [Biol Reprod. 110 \(2\): 355-364.](#)
13. Ferreira-Santos, L. *et al.* (2024) Neuraminidase-induced externalization of phosphatidylserine activates ADAM17 and impairs insulin signaling in endothelial cells. [Am J Physiol Heart Circ Physiol. 326 \(1\): H270-H277.](#)
14. Hibino, K. *et al.* (2024) Single-nucleosome imaging unveils that condensins and nucleosome-nucleosome interactions differentially constrain chromatin to organize mitotic chromosomes. [Nat Commun. 15 \(1\): 7152.](#)
15. Ramzan, F. *et al.* (2024) Fatty links between multisystem proteinopathy and small VCP-interacting protein. [Cell Death Discov. 10 \(1\): 358.](#)
16. Solís, C. *et al.* (2023) Cardiomyocyte external mechanical unloading activates modifications of α -actinin differently from sarcomere-originated unloading. [FEBS J. 290 \(22\): 5322-5339.](#)
17. Smith, S. *et al.* (2025) Placental ischemia during pregnancy induces hypertension, cerebral inflammation, and oxidative stress in dams postpartum. [Hypertens Pregnancy. 44 \(1\): 2454597.](#)
18. Morin, G. *et al.* (2025) Somatic PIK3R1 mutations in the iSH2 domain are accessible to PI3K α inhibition. [EMBO Mol Med. May 19 \[Epub ahead of print\].](#)

| | |
|--------------------------------------|---|
| Storage | 1 year at -20°C lyophilized; 6 months at +4°C after resuspension. This product is photosensitive and should be protected from light DO NOT FREEZE the solubilized material |
| Guarantee | Guaranteed until date of expiry. Please see product label. |
| Health And Safety Information | Material Safety Datasheet documentation #12004157 available at: https://www.bio-rad-antibodies.com/SDS/12004158 Material Safety Datasheet documentation #12004157 available at https://www.bio-rad-antibodies.com/uploads/MSDS/12004157.pdf |

Related Products

Recommended Useful Reagents

[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004161\)](#)

[GOAT ANTI-RABBIT IgG StarBright™ Blue 700 \(12004162\)](#)

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004163\)](#)

[ANTI-ACTIN hFAB™ RHODAMINE ANTIBODY \(12004164\)](#)

[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004165\)](#)

[ANTI-TUBULIN hFAB™ RHODAMINE ANTIBODY \(12004166\)](#)

[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004167\)](#)

[ANTI-GAPDH hFAB™ RHODAMINE ANTIBODY \(12004168\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M438163:250331'

Printed on 27 Jun 2025