

Datasheet: 1720-9007

| Description: | GOAT ANTI RAT CALCITONIN GENE-RELATED PEPTIDE |
|---------------|---|
| Specificity: | CGRP |
| Format: | Purified |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 0.1 ml |

Product Details

| Applications | This product has been re | ported to | work in th | ne following application | ns. This information is | |
|--|--|--|--|---|--|----------------|
| | derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further | | | | | |
| | | | | | | |
| | information. For general p | protocol re | ecommen | dations, please visit <u>w</u> | ww.bio- | |
| | rad-antibodies.com/proto | <u>cols</u> . | | | | |
| | | Yes | No | Not Determined | Suggested Dilution | |
| | Immunohistology - Frozen | | | | | |
| | Immunohistology - Paraffin | | | | | |
| | Immunohistology - Resin | | - | | | |
| | ELISA | • | | | 1/500 - 1/2500 | |
| | Western Blotting | | | | | |
| | Immunofluorescence | | | | | |
| | Where this product has n | ot been te | ested for | use in a particular tech | nnique this does not | |
| | necessarily exclude its us | se in such | procedu | res. Suggested workin | g dilutions are given as | |
| a guide only. It is recommended that the user titrates the product for use in the system using the appropriate negative/positive controls. | | | | | | |
| | | | | | | Target Species |
| Species Cross Reactivity | Reacts with: Mouse, Guir N.B. Antibody reactivity a reactivity is derived from personal communications further information. | nea Pig, E and workir testing wi s from the | Emu ng conditio thin our la originato | ons may vary between aboratories, peer-revie rs. Please refer to refe | a species. Cross wed publications or erences indicated for | |
| Product Form | Purified Ig - liquid | | | | | |
| Preparation | Purified Ig prepared by a | ffinity chro | omatogra | phy on Protein G | | |
| Buffer Solution | Rhaanhata huffarad aalin | <u>م</u> | | | | |
| | Phosphale bullered sailin | 6 | | | | |

Stabilisers

| Approx. Protein Concentrations | IgG concentration 5.0 mg/ml |
|-----------------------------------|---|
| Immunogen | Synthetic rat Tyr-CGRP (23-37) conjugated to gamma globulin. |
| External Database Links | UniProt: <u>P01256</u> <u>Related reagents</u> Entrez Gene: <u>24241</u> Calca <u>Related reagents</u> |
| Synonyms | Calc |
| RRID | AB_2290729 |
| Specificity | Goat anti Rat Calcitonin Gene-Related Peptide antibody recognizes Calcitonin gene-related peptide, also known as CGRP, a neuropeptide that acts as a vasodilator and plays a role in in the pathophysiology of migraine (<u>Recober & Russo 2009</u>). Goat anti Rat Calcitonin Gene-Related Peptide antibody reacts with both the whole molecule (amino acids 1-37) and the C-terminal fragment (23-37). |
| References | Collins, J.J. <i>et al.</i> (2000) Distribution and origin of secretoneurin-immunoreactive nerves in the female rat uterus. <u>Neuroscience. 95 (1): 255-64.</u> Chetty, R. <i>et al.</i> (2006) Pancreatic endocrine tumour with ductules: further observations of an unusual histological subtype. <u>Pathology. 38 (1): 5-9.</u> Weir, K.A. and Lunam, C.A. (2006) Immunohistochemical study of cutaneous nerves in the emu. <u>Cell Tissue Res. 326: 697-705.</u> Brock, J.A. <i>et al.</i> (2007) Postnatal androgen deprivation dissociates the development of smooth muscle innervation from functional neurotransmission in mouse vas deferens. J <u>Physiol. 581: 665-78.</u> Marchant, N.J. <i>et al.</i> (2007) Coexpression of prodynorphin and corticotrophin-releasing hormone in the rat central amygdala: evidence of two distinct endogenous opioid systems in the lateral division. J <u>Comp Neurol. 504: 702-15.</u> Pitzer, C. <i>et al.</i> (2008) Granulocyte-colony stimulating factor improves outcome in a mouse model of amyotrophic lateral sclerosis. <u>Brain. 131: 3335-47.</u> Ikeda, E. <i>et al.</i> (2009) Eputy functional bioengineered tooth replacement as an organ replacement therapy. <u>Proc Natl Acad Sci U S A. 106: 13475-80.</u> Bilff, J.J. <i>et al.</i> (2009) Disruption of laminin in the peripheral nervous system impedes nonmyelinating Schwann cell development and impairs nociceptive sensory function. <u>Glia.</u> <u>57: 850-9.</u> Golden, J.P. <i>et al.</i> (2010) RET signaling is required for survival and normal function of nonpeptidergic nociceptors. <u>J Neurosci. 30: 3983-94.</u> Iliff, J.J. <i>et al.</i> (2010) Epoxyeicosatrienoic acids are endogenous regulators of |

vasoactive neuropeptide release from trigeminal ganglion neurons. <u>J Neurochem. 115:</u> 1530-42.

12. Fan, W. *et al.* (2010) Structural and cellular features in metaphyseal and diaphyseal periosteum of osteoporotic rats. J Mol Histol. 41: 51-60.

13. Gnanamanickam, G.J. and Llewellyn-Smith, I.J. (2011) Innervation of the rat uterus at estrus: a study in full-thickness, immunoperoxidase-stained whole-mount preparations. J Comp Neurol. 519: 621-43.

14. Tague, S.E. and Smith, P.G. (2011) Vitamin D receptor and enzyme expression in dorsal root ganglia of adult female rats: modulation by ovarian hormones. <u>J Chem</u> <u>Neuroanat. 41: 1-12.</u>

15. Hamed, K. *et al.* (2011) Changes in cutaneous innervation in patients with chronic pain after burns. <u>Burns. 37: 631-7</u>

16. Zou, M. *et al.* (2012) Brn3a/Pou4f1 regulates dorsal root ganglion sensory neuron specification and axonal projection into the spinal cord. <u>Dev Biol. 364: 114-27.</u>

17. Zimmerman, A.L. et al. (2012) Monoaminergic Modulation of Spinal Viscero-

Sympathetic Function in the Neonatal Mouse Thoracic Spinal Cord PLoS One. 7: e47213.

18. Drummond, E.S. *et al.* (2014) Increased expression of cutaneous α1-adrenoceptors after chronic constriction injury in rats. <u>J Pain. 15 (2): 188-96.</u>

19. Drummond, P.D. *et al.* (2014) Upregulation of α 1-adrenoceptors on cutaneous nerve fibres after partial sciatic nerve ligation and in complex regional pain syndrome type II. <u>Pain. 155: 606-16.</u>

20. Li, Z. *et al.* (2014) Activation of MrgC receptor inhibits N-type calcium channels in small-diameter primary sensory neurons in mice. <u>Pain. 155 (8): 1613-21.</u>

21. Park, S.I. *et al.* (2015) Soft, stretchable, fully implantable miniaturized optoelectronic systems for wireless optogenetics. <u>Nat Biotechnol. 33 (12): 1280-1286.</u>

22. O'Brien, D.E. *et al.* (2015) ERK2 Alone Drives Inflammatory Pain But Cooperates with ERK1 in Sensory Neuron Survival. <u>J Neurosci. 35 (25): 9491-507.</u>

23. Vogelaar CF *et al.* (2015) Pharmacological Suppression of CNS Scarring by Deferoxamine Reduces Lesion Volume and Increases Regeneration in an *In Vitro* Model for Astroglial-Fibrotic Scarring and in Rat Spinal Cord Injury *In Vivo*. <u>PLoS One. 10 (7):</u> <u>e0134371.</u>

24. Sheahan, T.D. *et al.* (2015) Voluntary Exercise Training: Analysis of Mice in Uninjured, Inflammatory, and Nerve-Injured Pain States. <u>PLoS One. 10 (7): e0133191.</u>

25. Huang, A.Y. & Wu, S.Y. (2015) Calcitonin Gene-Related Peptide Reduces Taste-Evoked ATP Secretion from Mouse Taste Buds. <u>J Neurosci. 35 (37): 12714-24.</u>

26. Watanabe, M. *et al.* (2015) Expression and Regulation of Cav3.2 T-Type Calcium Channels during Inflammatory Hyperalgesia in Mouse Dorsal Root Ganglion Neurons. <u>PLoS One. 10 (5): e0127572.</u>

27. Chucair-Elliott, A.J. *et al.* (2015) Degeneration and regeneration of corneal nerves in response to HSV-1 infection. Invest Ophthalmol Vis Sci. 56 (2): 1097-107.

Van Steenwinckel, J. *et al.* (2015) Stromal cell-derived CCL2 drives neuropathic pain states through myeloid cell infiltration in injured nerve. <u>Brain Behav Immun. 45: 198-210.</u>
 Payne, S.C. *et al.* (2015) Regeneration of sensory but not motor axons following visceral nerve injury. <u>Exp Neurol. 266: 127-42.</u>

30. Valtcheva, M.V. *et al.* (2015) Enhanced nonpeptidergic intraepidermal fiber density and an expanded subset of chloroquine-responsive trigeminal neurons in a mouse model of dry skin itch. J Pain. 16 (4): 346-56.

31. Chucair-Elliott, A.J. *et al.* (2015) Degeneration and regeneration of corneal nerves in response to HSV-1 infection. <u>Invest Ophthalmol Vis Sci. 56 (2): 1097-107.</u>

32. Marvaldi L *et al.* (2015) Enhanced axon outgrowth and improved long-distance axon regeneration in sprouty2 deficient mice. <u>Dev Neurobiol. 75 (3): 217-31.</u>

33. Wong, A.W. *et al.* (2015) Neurite outgrowth in normal and injured primary sensory neurons reveals different regulation by nerve growth factor (NGF) and artemin. <u>Mol Cell</u> <u>Neurosci. 65: 125-34.</u>

34. Lin, S.H. *et al.* (2016) Evidence for the involvement of ASIC3 in sensory mechanotransduction in proprioceptors. <u>Nat Commun. 7: 11460.</u>

35. Kim, Y.S. *et al.* (2016) Coupled Activation of Primary Sensory Neurons Contributes to Chronic Pain. <u>Neuron. 91 (5): 1085-96.</u>

36. Takahashi, N. *et al.* (2016) Morphology of P2X3-immunoreactive nerve endings in the rat laryngeal mucosa. <u>Histochem Cell Biol. 145 (2): 131-46.</u>

37. Takiguchi, M. *et al.* (2018) Neonatal spinal injury induces *de novo* projections of primary afferents to the lumbosacral intermediolateral nucleus in rats. <u>IBRO Rep. 4: 1-6.</u>
38. Sun, Y. *et al.* (2020) Somatostatin neurons in the central amygdala mediate anxiety by disinhibition of the central sublenticular extended amygdala. <u>Mol Psychiatry. Oct 01 [Epub ahead of print].</u>

39. Rytel, L. & Całka, J. (2020) Aspirin Administration Affects Neurochemical Characterization of Substance P-Like Immunoreactive (SP-LI) Nodose Ganglia Neurons Supplying the Porcine Stomach. <u>Biomed Res Int. 2020: 1049179.</u>

40. Morellini, N. *et al.* (2020) Decreased neural expression of the noradrenaline transporter in the papillary dermis after partial sciatic nerve lesion. <u>J Chem Neuroanat.</u> <u>107: 101806.</u>

41. Lin, J.H. *et al.* (2021) ATF3-Expressing Large-Diameter Sensory Afferents at Acute Stage as Bio-Signatures of Persistent Pain Associated with Lumbar Radiculopathy. <u>Cells.</u> <u>10 (5): 992</u>

42. Godinho, M.J. *et al.* (2020) Regeneration of adult rat sensory and motor neuron axons through chimeric peroneal nerve grafts containing donor Schwann cells engineered to express different neurotrophic factors. <u>Exp Neurol. 330: 113355.</u>

43. Bertrand, M.M. *et al.* (2020) Functional segregation within the pelvic nerve of male rats: a meso- and microscopic analysis. <u>J Anat. 237 (4): 757-773.</u>

44. Morgan, M. *et al.* (2024) BDNF sensitizes bone and joint afferent neurons at different stages of MIA-induced osteoarthritis <u>Bone. : 117260.</u>

45. Zhang, F.M. *et al.* (2023) Transcriptional profiles of TGF-β superfamily members in the lumbar DRGs and the effects of activins A and C on inflammatory pain in rats. <u>J Physiol</u> <u>Biochem. 79 (2): 313-325.</u>

46. Alber, S. *et al.* (2023) PTBP1 regulates injury responses and sensory pathways in adult peripheral neurons. <u>Sci Adv. 9 (30): eadi0286.</u>

StorageThis product is shipped at ambient temperature. It is recommended to aliquot and store at
-20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for
short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

| Guarantee | 12 months from date of despatch | | | | |
|----------------------------------|---|--|--|--|--|
| Health And Safety Information | Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/1720-9007 10040 | | | | |
| Regulatory | For research purposes only | | | | |
| Related Products | | | | | |
| Recommended Secondary Antibodies | | | | | |

Rabbit Anti Goat IgG (Fc) (STAR122...) FITC, HRP

Recommended Useful Reagents

ANTIGEN RETRIEVAL BUFFER, pH8.0 (BUF025A)

| North & South | Tel: +1 800 265 7376 | Wor |
|---------------|-----------------------|-----------------|
| America | Fax: +1 919 878 3751 | |
| | Email: antibody_sales | _us@bio-rad.com |

rldwide

Tel: +44 (0)1865 852 700 Europe Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M438082:250321'

Printed on 21 Mar 2025

© 2025 Bio-Rad Laboratories Inc | Legal | Imprint