

Datasheet: 6120-0004 **BATCH NUMBER 150742** 

Description:	SHEEP ANTI METHOTREXATE
Specificity:	METHOTREXATE
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 ml

## **Product Details**

### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
ELISA	-			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Chemical
Product Form	Serum - liquid
Preservative Stabilisers	0.09% Sodium Azide (NaN <sub>3</sub> )
Immunogen	Methotrexate (4-amino-10-methylpteroyl-L-glutamic acid) conjugated via the carboxyl groups of the glutamic acid residue to KLH.
RRID	AB_619727
Specificity	Sheep anti methotrexate antibody recognizes Methotrexate, a drug originally

Sheep anti methotrexate antibody recognizes Methotrexate, a drug originally manufactured as a chemotherapeutic agent for the treatment or a number of cancers.

Methotrexate also acts as a disease-modifying anti-rheumatic drug used for the treatment of autoimmune conditions such as rheumatoid arthritis, Chron's disease and vasculitis.

Sheep anti methotrexate antibody exhibits the following cross reactivities with related compounds:-

Aminopterin	18%			
Folic acid	<0.01%			
Dyhydrofolic acid	<0.01%			
Tetrahydrofolic acid	<0.01%			
N-5-Methyltetrahydrofolic acid	<0.01%			
Folinic acid	<0.01%			
N-(p-Aminobenzoyl)-L-glutamic acid <0.01%				

#### References

- 1. Suzuki, Y. *et al.* (2003) The immunoassay of methotrexate by capillary electrophoresis with laser-induced fluorescence detection. Anal Sci. 19 (1): 111-5.
- 2. Tesfalidet, S. *et al.* (2016) Detection of methotrexate in a flow system using electrochemical impedance spectroscopy and multivariate data analysis. <u>Anal Chim Acta.</u> 914: 1-6.
- 3. Sereilakena, Pet al. (2017) Determination of methotrexate in spiked human blood serum using multi-frequency electrochemical immittance spectroscopy and multivariate data analysis <u>Analytica Chimica Acta. 987: 15-24.</u>
- 4. Phal, S. *et al.* (2018) Covalently electrografted carboxyphenyl layers onto gold surface serving as a platform for the construction of an immunosensor for detection of methotrexate Journal of Electroanalytical Chemistry. Jan 3 [Epub ahead of print].

### **Storage**

Guarantee

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

# Health And Safety

18 months from date of despatch.

Material Safety Datasheet documentation #10081 available at: <a href="https://www.bio-rad-antibodies.com/SDS/6120-0004">https://www.bio-rad-antibodies.com/SDS/6120-0004</a>

10081

**Regulatory** For research purposes only

## Related Products

## **Recommended Secondary Antibodies**

Rabbit Anti Sheep IgG (H/L) (5184-2304...) Biotin

North & South Tel: +1 800 265 7376 America Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody\_sales\_us@bio-rad.com

Email: antibody\_sales\_uk@bio-rad.com

Email: antibody\_sales\_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M354527:190724'

### Printed on 15 Mar 2024