

Datasheet: AAI19

BATCH NUMBER 156397

Description: SHEEP ANTI BOVINE IgM		
Specificity:	lgM	
Format:	Purified	
Product Type:	Polyclonal Antibody	
Isotype:	Polyclonal IgG	
Quantity:	1 mg	

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			•	
Immunohistology - Frozen				
Immunohistology - Paraffin				
ELISA	•			1/100 - 1/1000
Immunoprecipitation				
Western Blotting			•	
Immunodiffusion				

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Bovine	
Product Form	Purified IgG - liquid	
Antiserum Preparation	Antisera to bovine IgM were raised by repeated immunisation antigen. Purified IgG prepared by affinity chromatography.	of sheep with highly purified
Buffer Solution	Phosphate buffered saline	
Preservative Stabilisers	0.09% Sodium Azide	
Approx. Protein	IaG concentration 1.0 mg/ml	

IgG concentration 1.0 mg/ml

Concentrations

Concentrations	
Immunogen	Purified bovine IgM.
RRID	AB_323077
Specificity	Sheep anti bovine IgM recognizes bovine IgM and shows no cross-reactivity with other bovine immunoglobulin classes in immunoelectrophoresis. Sheep anti bovine IgM may cross-react with IgM from other species.
References	 Nebl, T. et al. (2002) Proteomic analysis of a detergent-resistant membrane skeleton from neutrophil plasma membranes. J Biol Chem. 277 (45): 43399-409. Assad, A. et al. (2012) Immunophenotyping and characterization of BNP colostra revealed pathogenic alloantibodies of IgG1 subclass with specifity to platelets, granulocytes and monocytes of all maturation stages. Vet Immunol Immunopathol. 147: 25-34. Hamsten, C. et al. (2009) Recombinant surface proteomics as a tool to analyze humoral immune responses in bovines infected by Mycoplasma mycoides subsp. mycoides small colony type. Mol Cell Proteomics. 8: 2544-54. Mansilla FC et al. (2015) Safety and immunogenicity of a soluble native Neospora caninum tachyzoite-extract vaccine formulated with a soy lecithin/β-glucan adjuvant in pregnant cattle. Vet Immunol Immunopathol. 165 (1-2): 75-80. Hossain, M.M. et al. (2016) Multiplex Detection of IgG and IgM to Rift Valley Fever Virus Nucleoprotein, Nonstructural Proteins, and Glycoprotein in Ovine and Bovine. Vector Borne Zoonotic Dis. Jul 5. [Epub ahead of print] Van Meulder, F. et al. (2015) Analysis of the protective immune response following intramuscular vaccination of calves against the intestinal parasite Cooperia oncophora. Int J Parasitol. 45 (9-10): 637-46. Jankowska, A. et al. (2016) Humoral and cellular immune response to Histophilus somni recombinant heat shock protein 60 kDa in farm animals Veterinární Medicína. 60 (No. 11): 603-12. Trotta, M. et al. (2015) Simultaneous immunization of cattle with foot-and-mouth disease (FMD) and live anthrax vaccines do not interfere with FMD booster responses Trials in Vaccinology. 4: 38-42.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: https://www.bio-rad-antibodies.com/SDS/AAI19 10040
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) Biotin

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M363594:200528'

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint