

Datasheet: AAI19B

Description:	SHEEP ANTI BOVINE IgM:Biotin
Specificity:	lgM
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			•	
Immunohistology - Frozen			•	
Immunohistology - Paraffin				
ELISA	•			1:10000 - 1:100000
Western Blotting	•			1:10000 - 1:100000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified IgG conjugated to Biotin - liquid

Antiserum Preparation Antisera to bovine IgM were raised by repeated immunisation of sheep with highly purified antigen. Purified IgG was prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified bovine IgM.

RRID	AB_10672760
Specificity	Sheep anti bovine IgM rec bovine immunoglobulin clas cross-react with IgM from ot
References	 Nebl, T. et al. (2002) Prot from neutrophil plasma men Assad, A. et al. (2012) Im revealed pathogenic alloant granulocytes and monocyte 3. Hamsten, C. et al. (2009) humoral immune responses mycoides small colony type.

cognizes bovine IgM and shows no cross-reactivity with other sses in immunoelectrophoresis. Sheep anti bovine IgM may ther species.

- teomic analysis of a detergent-resistant membrane skeleton mbranes. J Biol Chem. 277 (45): 43399-409.
- nmunophenotyping and characterization of BNP colostra tibodies of IgG1 subclass with specifity to platelets, es of all maturation stages. Vet Immunol Immunopathol. 147:
-) Recombinant surface proteomics as a tool to analyze s in bovines infected by Mycoplasma mycoides subsp. e. Mol Cell Proteomics. 8: 2544-54.
- 4. Mansilla FC et al. (2015) Safety and immunogenicity of a soluble native Neospora caninum tachyzoite-extract vaccine formulated with a soy lecithin/β-glucan adjuvant in pregnant cattle. Vet Immunol Immunopathol. 165 (1-2): 75-80.
- 5. Hossain, M.M. et al. (2016) Multiplex Detection of IgG and IgM to Rift Valley Fever Virus Nucleoprotein, Nonstructural Proteins, and Glycoprotein in Ovine and Bovine. Vector Borne Zoonotic Dis. 16 (8): 550-7.
- 6. Van Meulder, F. et al. (2015) Analysis of the protective immune response following intramuscular vaccination of calves against the intestinal parasite Cooperia oncophora. Int J Parasitol. 45 (9-10): 637-46.
- 7. Jankowska, A. et al. (2016) Humoral and cellular immune response to Histophilus somni recombinant heat shock protein 60 kDa in farm animals Veterinární Medicína. 60 (No. 11): 603-12.
- 8. Di Giacomo, S et al. (2022) Assessment on Different Vaccine Formulation Parameters in the Protection against Heterologous Challenge with FMDV in Cattle Viruses. 14 (8): 1781.
- 9. Springer, A. et al. (2022) Immunization Trials with Recombinant Major Sperm Protein of the Bovine Lungworm Dictyocaulus viviparus.. Pathogens 2022, 11, 55.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/AAI19B 10041
Regulatory	For research purposes only

America

North & South Tel: +1 800 265 7376 Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody sales us@bio-rad.com

Email: antibody sales uk@bio-rad.com

Email: antibody sales de@bio-rad.com

Printed on 01 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint