

Datasheet: AAI19F

Description:	SHEEP ANTI BOVINE IgM:FITC
Specificity:	IgM
Format:	FITC
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			1/20 - 1/100
Immunohistology - Frozen	▪			1/50 - 1/500
Immunohistology - Paraffin			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Bovine		
Product Form	Purified IgG fraction conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525

Antiserum Preparation Antisera to bovine IgM were raised by repeated immunisation of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	0.2% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	Purified bovine IgM.
RRID	AB_323076
Specificity	Sheep anti bovine IgM recognizes bovine IgM and shows no cross-reactivity with other bovine immunoglobulin classes in immunoelectrophoresis. Sheep anti bovine IgM may cross-react with IgM from other species.
References	<ol style="list-style-type: none"> 1. Nebl, T. <i>et al.</i> (2002) Proteomic analysis of a detergent-resistant membrane skeleton from neutrophil plasma membranes. J Biol Chem. 277 (45): 43399-409. 2. Assad, A. <i>et al.</i> (2012) Immunophenotyping and characterization of BNP colostrum revealed pathogenic alloantibodies of IgG1 subclass with specificity to platelets, granulocytes and monocytes of all maturation stages. Vet Immunol Immunopathol. 147: 25-34. 3. Hamsten, C. <i>et al.</i> (2009) Recombinant surface proteomics as a tool to analyze humoral immune responses in bovines infected by <i>Mycoplasma mycoides</i> subsp. <i>mycoides</i> small colony type. Mol Cell Proteomics. 8: 2544-54. 4. Mansilla FC <i>et al.</i> (2015) Safety and immunogenicity of a soluble native <i>Neospora caninum</i> tachyzoite-extract vaccine formulated with a soy lecithin/β-glucan adjuvant in pregnant cattle. Vet Immunol Immunopathol. 165 (1-2): 75-80. 5. Hossain, M.M. <i>et al.</i> (2016) Multiplex Detection of IgG and IgM to Rift Valley Fever Virus Nucleoprotein, Nonstructural Proteins, and Glycoprotein in Ovine and Bovine. Vector Borne Zoonotic Dis. Jul 5. [Epub ahead of print] 6. Van Meulder, F. <i>et al.</i> (2015) Analysis of the protective immune response following intramuscular vaccination of calves against the intestinal parasite <i>Cooperia oncophora</i>. Int J Parasitol. 45 (9-10): 637-46. 7. Jankowska, A. <i>et al.</i> (2016) Humoral and cellular immune response to <i>Histophilus somni</i> recombinant heat shock protein 60 kDa in farm animals Veterinárni Medicína. 60 (No. 11): 603-12. 8. Trotta, M. <i>et al.</i> (2015) Simultaneous immunization of cattle with foot-and-mouth disease (FMD) and live anthrax vaccines do not interfere with FMD booster responses Trials in Vaccinology. 4: 38-42.
Storage	<p>Store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Regulatory	For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M363597:200528'

Printed on 11 Feb 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)