

Datasheet: ANNEX200PE

BATCH NUMBER 161079

| | |
|----------------------|------------------------|
| Description: | ANNEXIN V:PE ASSAY KIT |
| Name: | ANNEXIN V KIT |
| Format: | RPE |
| Product Type: | Kits |
| Quantity: | 200 TESTS |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | Neat |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
|-----------|-----------------|---------------------|-------------------|
| | RPE 488nm laser | 496 | 578 |

| | |
|---------------------|--|
| Preservative | 0.09% Sodium Azide (NaN ₃) |
| Stabilisers | 1% Bovine Serum Albumin |

External Database Links

UniProt:

[P08758](#) [Related reagents](#)

Entrez Gene:

[308](#) ANXA5 [Related reagents](#)

| | |
|-----------------|-----------------|
| Synonyms | ANX5, ENX2, PP4 |
|-----------------|-----------------|

Product Information This test employs the property of Annexin V to bind to the membrane phospholipid phosphatidylserine (PS) in the presence of Ca²⁺. PS is exposed at the cell surface during the early stages of apoptosis. Detection of PS is a very sensitive method for detecting cells entering apoptosis, at a time point considerably ahead of nuclear changes such as DNA degradation.

The conjugation protocol used to prepare this product has not changed the native phospholipid binding properties of Annexin V. This protocol is designed to measure apoptosis easily and quickly in a sample of suspended cells.

[View our complete list of formats and sizes of Annexin V kits](#)

| | |
|----------------------------|---|
| Reagents In The Kit | Annexin V:PE 200 tests 7-AAD Viability Staining Solution 2X 100 tests 10X Binding Buffer 100 ml |
|----------------------------|---|

| | |
|-----------------------------|---|
| Instructions For Use | <ol style="list-style-type: none">1) Dilute the 10X Binding Buffer to 1X in distilled water (1 ml Binding Buffer + 9 ml distilled water).2) Wash cells once in PBS by gentle shaking or pipetting up and down. Then wash once in 1X Binding Buffer.3) Resuspend cells in 1X Binding Buffer, adjusting to a cell density of $1-5 \times 10^6$ cells/ml.4) Add 5 μl of PE conjugated Annexin V to 100 μl of the cell suspension.5) Mix and incubate for 10 to 15 minutes at room temperature.6) Wash cells in 1X Binding Buffer.7) Resuspend cells in 200 μl of 1X Binding Buffer.8) Add 5 μl of 7-AAD Viability Staining Solution and incubate for 5 minutes at room temperature.9) Analyse by flow cytometry within 4 hours. Store at 2-8°C in the dark. |
|-----------------------------|---|

| | |
|----------------|--|
| Storage | Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use. |
|----------------|--|

| | |
|------------------|--|
| Guarantee | Guaranteed until date of expiry. Please see product label. |
|------------------|--|

| | |
|--------------------------------------|---|
| Health And Safety Information | Material Safety Datasheet documentation #10229 #10230 #10582 available at: Recombinant Annexin V (10229): https://www.bio-rad-antibodies.com/uploads/MSDS/10229.pdf Binding Buffer (10230): https://www.bio-rad-antibodies.com/uploads/MSDS/10230.pdf 7-AAD (10582): https://www.bio-rad-antibodies.com/uploads/MSDS/10582.pdf |
|--------------------------------------|---|

| | |
|-------------------|----------------------------|
| Regulatory | For research purposes only |
|-------------------|----------------------------|

Related Products

Recommended Useful Reagents

[ANNEXIN V:PE ASSAY KIT \(ANNEX50PE\)](#)

[ANNEXIN V:APC ASSAY KIT \(ANNEX50APC\)](#)

[ANNEXIN V:APC ASSAY KIT \(ANNEX200APC\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M371229:200608'

Printed on 08 Feb 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)