

Datasheet: APO006

Description:	FLISP™ FAM-Phe-CMK SERINE PROTEASE ASSAY KIT
Name:	SERINE PROTEASE ASSAY KIT
Format:	FAM (Green Fluorescence)
Product Type:	Kits
Quantity:	25 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FAM	494	520

Product Information

Fluorescent Labeled Inhibitors of Serine Protease (FLISP) kits allow the analysis of intracellular chymotrypsin-like serine protease activity in live cells. Serine protease activation plays major roles in several different cellular functions including apoptosis, tumor malignancy and as diagnostic and prognostic indicators of breast, head and neck carcinomas. Serine proteases are also thought to play a role in infections and transplant rejection.

Test Principle

Active chymotrypsin-like enzymes within the cell bind covalently to the FLISP FAM-Phe-CMK inhibitor probe, retaining the green fluorescent signal within the cell. The higher the concentration of active chymotrypsin-like enzymes within the cell, the higher the level of fluorescence observed. FLISP kits are suitable for use with suspension or adherent cells and can be used with a flow cytometer, fluorescence microscope or a fluorescence plate reader.

FLISP inhibitors are cell-permeable and non-cytotoxic fluorochrome-labeled analogs of the first serine protease inhibitor, tosyl-phenylalanylchloromethyl ketone (TPCK), labeled with

carboxyfluorescein (FAM) which fluoresces green. They are available with either a chloromethyl ketone (CMK) or diphenyl 1-(N-peptidylamino) alkanephosphonate (DAP) reactive group containing compound. Specificity for chymotrypsin-like enzymes is determined by the amino acid targeting peptide, the FAM-Phe-CMK probe inhibits chymotrypsin.

Reagents In The Kit FAM-Phe-CMK FLISP Reagent - lyophilized
10X Wash Buffer, 15 mL
Fixative, 6 mL
Propidium Iodide, 1 mL
Hoechst Stain, 1 mL

Instructions For Use Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/APO006.pdf

Storage MULTIPLE STORAGE CONDITIONS APPLY ON ARRIVAL. Store the unopened kit (and each unopened component) according to the storage instructions on each component label.
Store the FLISP reagent at -20°C. Once reconstituted in DMSO, the FLISP stock should be stored at -20°C for 6 months, protected from light. Avoid repeated freezing and thawing.

Guarantee Guaranteed until date of expiry. Please see product label.

Acknowledgements FLISP is a trademark of Immunochemistry Technologies, LLC.

Health And Safety Information Material Safety Datasheet documentation #20399 #10471 #10498 #10476 #10477 available at:
FAM-Phe-CMK FLISP Reagent (20399): <https://www.bio-rad-antibodies.com/uploads/MSDS/20399.pdf>
10X Wash Buffer (10471): <https://www.bio-rad-antibodies.com/uploads/MSDS/10471.pdf>
Fixative (10498): <https://www.bio-rad-antibodies.com/uploads/MSDS/10498.pdf>
Hoechst Stain (10476): <https://www.bio-rad-antibodies.com/uploads/MSDS/10476.pdf>
Propidium Iodide Reagent (10477): <https://www.bio-rad-antibodies.com/uploads/MSDS/10477.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[SR FLICA™ POLY CASPASE KIT \(ICT916\)](#)

[SR FLICA™ POLY CASPASE KIT \(ICT917\)](#)

[FAM FLICA™ POLY CASPASE KIT \(ICT091\)](#)

[FAM FLICA™ POLY CASPASE KIT \(ICT092\)](#)

[FLICA™ 660 POLY CASPASE KIT \(ICT9120\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M374916:201208'

Printed on 29 Aug 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)