

Datasheet: APO009

Description:	FLISP™ FAM-Phe-DAP SERINE PROTEASE ASSAY KIT
Name:	SERINE PROTEASE ASSAY KIT
Format:	FAM (Green Fluorescence)
Product Type:	Kits
Quantity:	25 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FAM	494	520

Product Information

Fluorescent Labeled Inhibitors of Serine Protease (FLISP) kits allow the analysis of intracellular chymotrypsin-like serine protease activity in live cells. Serine protease activation plays major roles in several different cellular functions including apoptosis, tumor malignancy and as diagnostic and prognostic indicators of breast, head and neck carcinomas. Serine proteases are also thought to play a role in infections and transplant rejection.

Test Principle

Active chymotrypsin-like enzymes within the cell bind covalently to the FLISP FAM-Phe-DAP inhibitor probe, retaining the green fluorescent signal within the cell. The higher the concentration of active chymotrypsin-like enzymes within the cell, the higher the level of fluorescence observed. FLISP kits are suitable for use with suspension or adherent cells and can be used with a flow cytometer, fluorescence microscope or a fluorescence plate reader.

FLISP inhibitors are cell-permeable and non-cytotoxic fluorochrome-labeled analogs of the first serine protease inhibitor, tosyl-phenylalanylchloromethyl ketone (TPCK), labeled with

carboxyfluorescein (FAM) which fluoresces green. They are available with either a chloromethyl ketone (CMK) or diphenyl 1-(N-peptidylamino) alkanephosphonate (DAP) reactive group containing compound. Specificity for chymotrypsin-like enzymes is determined by the amino acid targeting peptide, the FAM-Phe-DAP probe inhibits chymotrypsin.

Reagents In The Kit	FAM-Phe-DAP FLISP Reagent - lyophilized 10X Wash Buffer, 15 mL Fixative, 6 mL Propidium Iodide, 1 mL Hoechst Stain, 1 mL
Instructions For Use	Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/APO009.pdf
Storage	MULTIPLE STORAGE CONDITIONS APPLY ON ARRIVAL. Store the unopened kit (and each unopened component) according to the storage instructions on each component label. Store the FLISP reagent at -20°C. Once reconstituted in DMSO, the FLISP stock should be stored at -20°C for 6 months, protected from light. Avoid repeated freezing and thawing.
Guarantee	Guaranteed until date of expiry. Please see product label.
Acknowledgements	FLISP is a trademark of Immunochemistry Technologies, LLC.
Health And Safety Information	Material Safety Datasheet documentation #20400 #10471 #10498 #10476 #10477 available at: FAM-Phe-DAP FLISP Reagent (20400): https://www.bio-rad-antibodies.com/uploads/MSDS/20400.pdf 10X Apoptosis Wash Buffer (10471): https://www.bio-rad-antibodies.com/uploads/MSDS/10471.pdf Fixative (10498): https://www.bio-rad-antibodies.com/uploads/MSDS/10498.pdf Hoechst Stain (10476): https://www.bio-rad-antibodies.com/uploads/MSDS/10476.pdf Propidium Iodide Reagent (10477): https://www.bio-rad-antibodies.com/uploads/MSDS/10477.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[SR FLICA™ POLY CASPASE KIT \(ICT916\)](#)
[SR FLICA™ POLY CASPASE KIT \(ICT917\)](#)
[FLICA™ 660 POLY CASPASE KIT \(ICT9120\)](#)
[FAM FLICA™ POLY CASPASE KIT \(ICT091\)](#)
[FAM FLICA™ POLY CASPASE KIT \(ICT092\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M374919:201208'

Printed on 24 Aug 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)