

Datasheet: BUF017B **BATCH NUMBER 156453**

Name: PEROXIDE BLOCKING I	PEROXIDE BLOCKING REAGENT		
	PEROXIDE BLOCKING REAGENT		
Format: Reagent			
Product Type: Accessory Reagent	Accessory Reagent		
Quantity: 50 ml			

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.biorad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	•			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Intended Use

BUF017B is a 3% (v/v) solution of hydrogen peroxide. This solution may be used to block endogenous peroxidase activity in paraffin embedded tissue sections. Incubate re-hydrated sections in this buffer for 15 minutes, and rinse well with water before continuing with staining protocol.

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10233 available at: https://www.bio-rad-antibodies.com/SDS/BUF017B

For research purposes only

America

Regulatory

North & South Tel: +1 800 265 7376 Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

10233

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint