

## Datasheet: BUF036A BATCH NUMBER 154880

| Description:  | 10X PHOSPHATE BUFFERED SALINE |
|---------------|-------------------------------|
| Name:         | 10X PHOSPHATE BUFFERED SALINE |
| Format:       | 10 X Concentrate              |
| Product Type: | Accessory Reagent             |
| Quantity:     | 1 LITRE                       |

## **Product Details**

| Applications                     | This product has been reported to work in the following applications. This information derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further  |      |    |                |                    |  |  |  |
|----------------------------------|--|------|----|----------------|--------------------|--|--|--|
|                                  | information. For general protocol recommendations, please visit <u>www.bio-</u><br>rad-antibodies.com/protocols.   |      |    |                |                    |  |  |  |
|                                  |  |      |    |                |                    |  |  |  |
|                                  |  | Yes  | No | Not Determined | Suggested Dilution |  |  |  |
|                                  | ELISA  | •    |    |                |                    |  |  |  |
|                                  | Immunoassay  | -    |    |                | 1/10               |  |  |  |
|                                  | Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls. |      |    |                |                    |  |  |  |
| Product Form                     | 10 x concentrate - liquid pH7.3 +/-0.1   |      |    |                |                    |  |  |  |
| Intended Use                     | Dilute 1:9 with distilled water to prepare solution of 0.01M PBS for use in washing of sections in immunohistology staining procedures.  |      |    |                |                    |  |  |  |
| Storage                          | Store at room temperature  |      |    |                |                    |  |  |  |
| Guarantee                        | Guaranteed until date of expiry. Please see product label.   |      |    |                |                    |  |  |  |
| Health And Safety<br>Information | Material Safety Datasheet documentation #10216 available at:<br>https://www.bio-rad-antibodies.com/SDS/BUF036A<br>10216  |      |    |                |                    |  |  |  |
| Regulatory                       | For research purposes of   | only |    |                |                    |  |  |  |
|                                  |  |      |    |                |                    |  |  |  |

| North & South | Tel: +1 800 265 7376               | Vorldwide | Tel: +44 (0)1865 852 700        | Europe | Tel: +49 (0) 89 8090 95 21           |
|---------------|------------------------------------|-----------|---------------------------------|--------|--------------------------------------|
| America       | Fax: +1 919 878 3751               |           | Fax: +44 (0)1865 852 739        |        | Fax: +49 (0) 89 8090 95 50           |
|               | Email: antibody_sales_us@bio-rad.c | om        | Email: antibody_sales_uk@bio-ra | d.com  | Email: antibody_sales_de@bio-rad.com |

## Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint