

Datasheet: BUF036A BATCH NUMBER 161714

| Description: | 10X PHOSPHATE BUFFERED SALINE |
|---------------|-------------------------------|
| Name: | 10X PHOSPHATE BUFFERED SALINE |
| Format: | 10 X Concentrate |
| Product Type: | Accessory Reagent |
| Quantity: | 1 LITRE |

Product Details

| Applications | This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal | | | | | | | | | | | | | |
|----------------------------------|---|------|--|--|------|--|--|--|--|-------|-----|----|----------------|--------------------|
| | communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-</u> rad-antibodies.com/protocols. | | | | | | | | | | | | | |
| | | | | | | | | | | | Yes | No | Not Determined | Suggested Dilution |
| | | | | | | | | | | ELISA | | | | |
| | Immunoassay | • | | | 1/10 | | | | | | | | | |
| | Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls. | | | | | | | | | | | | | |
| Product Form | 10 x concentrate - liquid pH7.3 +/-0.1 | | | | | | | | | | | | | |
| Intended Use | Dilute 1:9 with distilled water to prepare solution of 0.01M PBS for use in washing of sections in immunohistology staining procedures. | | | | | | | | | | | | | |
| Storage | Store at room temperature | | | | | | | | | | | | | |
| Guarantee | Guaranteed until date of expiry. Please see product label. | | | | | | | | | | | | | |
| Health And Safety Information | Material Safety Datasheet documentation #10216 available at: https://www.bio-rad-antibodies.com/SDS/BUF036A 10216 | | | | | | | | | | | | | |
| Regulatory | For research purposes of | only | | | | | | | | | | | | |
| | | | | | | | | | | | | | | |

| North & South | Tel: +1 800 265 7376 | Worldwide | Tel: +44 (0)1865 852 700 | Europe | Tel: +49 (0) 89 8090 95 21 |
|---------------|---------------------------------|-----------|--------------------------------|--------|--------------------------------------|
| America | Fax: +1 919 878 3751 | | Fax: +44 (0)1865 852 739 | | Fax: +49 (0) 89 8090 95 50 |
| | Email: antibody_sales_us@bio-ra | id.com | Email: antibody_sales_uk@bio-r | ad.com | Email: antibody_sales_de@bio-rad.com |

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint