

Datasheet: BUF039A

Description:	ELISA NEPTUNE ASSAY DILUENT
Name:	ELISA NEPTUNE ASSAY DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Ready to use - liquid
---------------------	-----------------------

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative Stabilisers	0.1% Sodium Azide
---------------------------------	-------------------

Product Information **ELISA Neptune Assay Diluent** is an ELISA assay diluent designed to reduce background noise in problematic human, porcine and bovine serum or plasma samples. Specially formulated for use in antigen-down ELISA, it contains non-mammalian protein and proprietary detergent additives to reduce non-specific IgG adsorption to antigen coated and blocked ELISA well surfaces.

The difference in complexity, between serum or plasma samples and diluted standard curve proteins, is a potential cause of target analyte under-recovery. BUF039A equalizes the matrix, inhibits complement and thrombin activity and contains an antimicrobial agent.

Instructions For Use	Pipette buffer directly onto the plate just prior to adding the sample or standard (suggested volume 50 - 100ul/well).
-----------------------------	--

References	1. Defresne, F. <i>et al.</i> (2010) Differential influence of anticancer treatments and
-------------------	--

angiogenesis on the seric titer of autoantibody used as tumor and metastasis biomarker.

[Neoplasia 12: 562-70.](#)

2. Massei, G. *et al.* (2015) Immunocontraception for managing feral cattle in Hong Kong.

[PLoS One. 10 \(4\): e0121598.](#)

Storage	Store at +4°C. DO NOT FREEZE
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10377 available at: https://www.bio-rad-antibodies.com/SDS/BUF039A 10377
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[ELISA ULTRABLOCK \(BUF033A\)](#)

[ELISA ULTRABLOCK \(BUF033B\)](#)

[ELISA ULTRABLOCK \(BUF033C\)](#)

[5x ELISA COATING BUFFER \(BUF030A\)](#)

[10x ELISA WASH BUFFER \(BUF031A\)](#)

[ELISA BSA BLOCK \(BUF032A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350297:190307'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)