

Datasheet: BUF039B

BATCH NUMBER 156031

Description:	ELISA NEPTUNE ASSAY DILUENT
Name:	ELISA NEPTUNE ASSAY DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	500 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Ready to use - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide

Product Information **ELISA Neptune Assay Diluent** is an ELISA assay diluent designed to reduce background noise in problematic human, porcine and bovine serum or plasma samples. Specially formulated for use in antigen-down ELISA, it contains non-mammalian protein and proprietary detergent additives to reduce non-specific IgG adsorption to antigen coated and blocked ELISA well surfaces.

The difference in complexity, between serum or plasma samples and diluted standard curve proteins, is a potential cause of target analyte under-recovery. BUF039A equalizes the matrix, inhibits complement and thrombin activity and contains an antimicrobial agent.

Instructions For Use Pipette buffer directly onto the plate just prior to adding the sample or standard (suggested volume 50 - 100ul/well).

References	1. Defresne, F. <i>et al.</i> (2010) Differential influence of anticancer treatments and angiogenesis on the seric titer of autoantibody used as tumor and metastasis biomarker. Neoplasia 12: 562-70. 2. Massei, G. <i>et al.</i> (2015) Immunocontraception for managing feral cattle in Hong Kong. PLoS One. 10 (4): e0121598.
-------------------	--

Storage	Store at +4°C. DO NOT FREEZE
----------------	-------------------------------------

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10377 available at: https://www.bio-rad-antibodies.com/SDS/BUF039B 10377
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

- [5x ELISA COATING BUFFER \(BUF030A\)](#)
- [ELISA ULTRABLOCK \(BUF033B\)](#)
- [ELISA ULTRABLOCK \(BUF033C\)](#)
- [ELISA ULTRABLOCK \(BUF033A\)](#)
- [10x ELISA WASH BUFFER \(BUF031A\)](#)
- [ELISA BSA BLOCK \(BUF032A\)](#)

North & South America Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
---	---	--

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
 'M350298:190307'

Printed on 23 May 2025