

Datasheet: BUF071

BATCH NUMBER 155907

Description:	FIXATION BUFFER
Name:	FIXATION BUFFER
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	Paraformaldehyde

Product Information **Fixation Buffer** can be used for preparing cells which will undergo staining for intracellular cytokines or proteins. The Fixation Buffer should be used prior to permeabilization using ICS Perm/Wash Buffer ([BUF072](#)) and has been designed to produce the greatest signal to noise ratio.

Instructions For Use

1. Resuspend cells in 0.5 ml Fixation Buffer
2. Leave in dark for 20 min at room temperature

Note: It is recommended that the reagent be titrated for optimal performance for each application

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #20367 available at:
<https://www.bio-rad-antibodies.com/SDS/BUF071>
20367

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350325:190307'

Printed on 29 Feb 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)