

Datasheet: BUF071

BATCH NUMBER 159802

Description:	FIXATION BUFFER
Name:	FIXATION BUFFER
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	Paraformaldehyde

Product Information **Fixation Buffer** can be used for preparing cells which will undergo staining for intracellular cytokines or proteins. The Fixation Buffer should be used prior to permeabilization using ICS Perm/Wash Buffer ([BUF072](#)) and has been designed to produce the greatest signal to noise ratio.

- Instructions For Use**
1. Resuspend cells in 0.5 ml Fixation Buffer
 2. Leave in dark for 20 min at room temperature

Note: It is recommended that the reagent be titrated for optimal performance for each application

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #20367 available at:
<https://www.bio-rad-antibodies.com/SDS/BUF071>
20367

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350325:190307'

Printed on 29 Feb 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)