

Datasheet: C12CA BATCH NUMBER 165327

Description:	BABY RABBIT COMPLEMENT
Name:	BABY RABBIT COMPLEMENT
Format:	Serum
Product Type:	Serum
Quantity:	2 ml

Product Details

Applications	This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-</u> rad-antibodies.com/protocols.						
	Yes No Not Determined Suggested Dilution						
	Functional Assays (1)	•					
	Immunoassay	-					
	Where this product has not been tested for use in a particular technique this does not						
	necessarily exclude its u a guide only. It is recomm system using appropriate (1) This product is not It is preferable to dilute filtration in order to min	nended th e negative sold as s e the com	at the use /positive c terile but plement t	er titrates the product fo controls. can be sterilized by t to a final working cor	or use in their own filtration if necessary.		
Product Form	Baby rabbit serum - lyophilized						
Reconstitution	Reconstitute with 2ml ice cold distilled water						
Preservative Stabilisers	None present						
Product Information	Baby rabbit complement serum preparation is intended for use as a source of rabbit complement for cytotoxicity assays.						
Instructions For Use	Use within one hour of reconstitution, keeping on ice throughout.						
References	1. De clercq, L. <i>et al.</i> (19 preadipocytes <i>in vitro</i> an <u>75 (7): 1791-7.</u> 2. Anderson, L.D. Jr <i>et a</i>	d depress	es the dev	velopment of pig adipc	ose tissue. <u>J Anim Sci.</u>		

versus-host disease by pretransplant immunization of allogeneic bone marrow donors with a recipient-derived tumor cell vaccine. <u>Cancer Res. 59 (7): 1525-30.</u>

 Lidington, E.A. *et al.* (2000) Induction of decay-accelerating factor by thrombin through a protease-activated receptor 1 and protein kinase C-dependent pathway protects vascular endothelial cells from complement-mediated injury. <u>Blood. 96 (8): 2784-92.</u>
 Mason, J.C. *et al.* (2002) bFGF and VEGF synergistically enhance endothelial

cytoprotection via decay-accelerating factor induction. <u>Am J Physiol Cell Physiol. 282:</u> <u>C578-87.</u>

5. Mason, J.C. *et al.* (2002) Statin-induced expression of decay-accelerating factor protects vascular endothelium against complement-mediated injury. <u>Circ Res. 91 (8):</u> 696-703.

6. Li, S.H. *et al.* (2004) C-reactive protein upregulates complement-inhibitory factors in endothelial cells. <u>Circulation. 109: 833-6.</u>

7. Newcombe, J. *et al.* (2004) Infection with an avirulent phoP mutant of *Neisseria meningitidis* confers broad cross-reactive immunity. Infect Immun. 72: 338-44.

8. Sancho, D. *et al.* (2006) CD69 targeting differentially affects the course of collageninduced arthritis. <u>J Leukoc Biol. 80: 1233-41.</u>

9. Hyams, C. *et al.* (2010) *Streptococcus pneumoniae* resistance to complement-mediated immunity is dependent on the capsular serotype. <u>Infect Immun. 78: 716-25.</u>

10. Hung, M.C. *et al.* (2011) The *Neisseria meningitidis* Macrophage Infectivity Potentiator Protein Induces Cross-Strain Serum Bactericidal Activity and Is a Potential Serogroup B Vaccine Candidate. <u>Infect Immun. 79: 3784-91.</u>

11. Lee, S.J. *et al.* (2012) Identification of a common immune signature in murine and human systemic Salmonellosis. <u>Proc Natl Acad Sci U S A. 109 (13): 4998-5003.</u>

12. Hung MC *et al.* (2013) The adhesin complex protein (ACP) of *Neisseria meningitidis* is a new adhesin with vaccine potential. <u>MBio. 4 (2): pii: e00041-13.</u>

13. Goh, Y.S. & MacLennan, C.A. (2013) Invasive African nontyphoidal Salmonella requires high levels of complement for cell-free antibody-dependent killing. <u>J Immunol</u> <u>Methods. 387 (1-2): 121-9.</u>

14. Goh YS *et al.* (2016) Bactericidal Immunity to *Salmonella* in Africans and Mechanisms Causing Its Failure in HIV Infection. <u>PLoS Negl Trop Dis. 10 (4): e0004604.</u>

15. Humbert MV *et al.* (2016) Vaccine Potential and Diversity of the Putative Cell Binding Factor (CBF, NMB0345/NEIS1825) Protein of *Neisseria meningitidis*. <u>PLoS One. 11 (8):</u> <u>e0160403.</u>

16. Dierckx de Casterlé I *et al.* (2018) Reduction of myeloid-derived suppressor cells reinforces the anti-solid tumor effect of recipient leukocyte infusion in murine neuroblastoma-bearing allogeneic bone marrow chimeras. <u>Cancer Immunol Immunother.</u> <u>67 (4): 589-603.</u>

17. Valton, J. *et al.* (2018) A Versatile Safeguard for Chimeric Antigen Receptor T-Cell Immunotherapies. <u>Sci Rep. 8 (1): 8972.</u>

18. Dierckx de Casterlé, I. *et al.* (2018) Reduction of myeloid-derived suppressor cells reinforces the anti-solid tumor effect of recipient leukocyte infusion in murine neuroblastoma-bearing allogeneic bone marrow chimeras. <u>Cancer Immunol Immunother.</u> <u>67 (4): 589-603.</u>

19. Nganje, C.N. *et al.* (2019) PepN is a non-essential, cell wall-localized protein that contributes to neutrophil elastase-mediated killing of *Streptococcus pneumoniae*. <u>PLoS</u> <u>One. 14 (2): e0211632</u>.

	20. Cuesta-Mateos, C. <i>et al.</i> (2020) CCR7 as a novel therapeut PROLYMPHOCYTIC leukemia <u>Biomarker Research.8, 54.</u>					
	21. Mosti, L. <i>et al.</i> (21. Mosti, L. <i>et al.</i> (2021) Targeted multi-epitope switching enables straightforward				
	positive/negative se	election of CAR T cells. Ge	ene Ther. 28 (9) <u>: 602-12.</u>		
Storage	Prior to reconstitution aliquot and store at	n store at +4ºC for 1 hour or				
	•	ct. Should this product cor	•	ezing and thawing as this may ate we recommend		
Guarantee						
Health And SafetyMaterial Safety Datasheet documentation #10288 avaInformationhttps://www.bio-rad-antibodies.com/SDS/C12CA			able at:			
	10288					
Regulatory	Regulatory For research purposes only					
North & South Tel: +1 800 2 America Fax: +1 919 8		 Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50		
Email: antibo	dy_sales_us@bio-rad.com	Email: antibody_sales_uk@bi	o-rad.com	Email: antibody_sales_de@bio-rad.com		
To find a batch/lot spe	ecific datasheet for this pro	oduct, please use our online 'M391022:211008'	search tool at: b	bio-rad-antibodies.com/datasheets		
Printed on 18 Jan 2024						

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint