

Datasheet: HIS003B

BATCH NUMBER 155241

Description:	3x ANTIGEN RETRIEVAL BUFFER, UNIVERSAL
Name:	ANTIGEN RETRIEVAL BUFFER, Universal
Format:	3 X Concentrate
Product Type:	Accessory Reagent
Quantity:	250 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Buffer Solution	Bio-Rad's Antigen unmasking fluid is supplied as 3x concentrated fluid. Crystals in the fluid will dissolve during heating before use.
------------------------	--

Intended Use	Bio-Rad's Antigen unmasking fluid is a heat induced, epitope retrieval system. The 'universal' formulation is suitable for use with most antibodies and gives excellent preservation of tissue morphology.
---------------------	--

Instructions For Use	<ol style="list-style-type: none"> 1. Deparaffinize sections and re-hydrate. 2. Dilute concentrated Antigen unmasking fluid 1:3 with de-ionised water (1 part Antigen unmasking fluid:2 parts water). 3. Pre-heat retrieval buffer to 92°C *. 4. Place slides in buffer and incubate 10 minutes at 92°C *. 5. Allow to cool for 5-10 minutes. 6. Rinse in distilled water or PBS, proceed with staining.
-----------------------------	--

* N.B. Optimal unmasking is a function of time and temperature. Temperature (range 90°C-100°C) and time (5-30 minutes) should be determined by each investigator.

Storage	Store at room temperature. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10270 available at: https://www.bio-rad-antibodies.com/SDS/HIS003B 10270
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M374958:201210'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)