

Datasheet: ICT091

Description:	FAM FLICA™ POLY CASPASE KIT
Name:	POLY CASPASE
Format:	FAM (Green Fluorescence)
Product Type:	Kits
Quantity:	25 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FAM	494	520

Product Information

The FAM FLICA™ Poly Caspase kit uses a quick and easy method to analyze active caspases in apoptotic cells.

This FLICA™ reagent binds the following activated caspases: Caspase 1, Caspase 3, Caspase 4, Caspase 5, Caspase 6, Caspase 7, Caspase 8 and Caspase 9. It can also be used as a generic probe for the detection of most caspases.

Test Principle

Caspase FLICA™ kits measure apoptosis by detecting active caspases in whole, living cells. These kits do not work by using antibodies or as an ELISA. Instead, their methodology is based on a unique cell-permeable and non-cytotoxic reagent called the Fluorochrome Inhibitor of Caspases (FLICA). The FLICA™ reagent contains a caspase inhibitor sequence linked to a green (Carboxyfluorescein, FAM) fluorescent probe.

The Caspase FLICA™ Kits are suitable for cells in suspension and adherent cells from many species including mammalian, insect and yeast. Different cell types, e.g. HeLa, primary neurons, macrophages and lymphocytes have also been successfully studied with these kits.

This kit can be used with a flow cytometer, fluorescence microscope or a fluorescence plate reader using black microtitre plates.

Reagents In The Kit

1 vial of FAM-VAD-FMK FLICA™ Reagent - lyophilized
 10X Wash Buffer, 15 mL
 Fixative, 6 mL

Propidium Iodide, 1 mL
Hoechst Stain, 1 mL

Instructions For Use Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/ICT091.pdf

References

1. Brugnon, F. *et al.* (2010) Apoptosis and meiotic segregation in ejaculated sperm from Robertsonian translocation carrier patients. [Hum Reprod. 25 \(7\): 1631-42.](#)
2. Peppas, D. *et al.* (2013) Up-regulation of a death receptor renders antiviral T cells susceptible to NK cell-mediated deletion. [J Exp Med. 210 \(1\): 99-114.](#)
3. Chua, M.L. *et al.* (2014) DNA double-strand break repair and induction of apoptosis in *ex vivo* irradiated blood lymphocytes in relation to late normal tissue reactions following breast radiotherapy. [Radiat Environ Biophys. 53 \(2\): 355-64.](#)

Storage Store the unopened kit (and each unopened component) at +4°C until the expiration date. Protect the FLICA™ reagent from light at all times.

Once reconstituted, the 150X FLICA™ stock should be stored at -20°C protected from light.

Guarantee Guaranteed until date of expiry. Please see product label.

Acknowledgements FLICA™ is a trademark of Immunochemistry Technologies, LLC.

Health And Safety Information Material Safety Datasheet documentation #20286 #10471 #10498 #10476 #10477 available at: FAM-VAD-FMK FLICA™ Reagent (20286): <https://www.bio-rad-antibodies.com/uploads/MSDS/20286.pdf>
10X Wash Buffer (10471): <https://www.bio-rad-antibodies.com/uploads/MSDS/10471.pdf>
Fixative (10498): <https://www.bio-rad-antibodies.com/uploads/MSDS/10498.pdf>
Hoechst Stain (10476): <https://www.bio-rad-antibodies.com/uploads/MSDS/10476.pdf>
Propidium Iodide Reagent (10477): <https://www.bio-rad-antibodies.com/uploads/MSDS/10477.pdf>

Licensed Use Please see PDF datasheet for full Health & Safety information.

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[FAM FLICA™ POLY CASPASE KIT \(ICT092\)](#)

[FAM FLICA™ CASPASE-3/7 KIT \(ICT093\)](#)

[FAM FLICA™ CASPASE-1 KIT \(ICT097\)](#)

[FAM FLICA™ CASPASE-2 KIT \(ICT918\)](#)

[FAM FLICA™ CASPASE-6 KIT \(ICT095\)](#)

[FAM FLICA™ CASPASE-8 KIT \(ICT099\)](#)

[FAM FLICA™ CASPASE-9 KIT \(ICT912\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M350345:190307'

Printed on 20 May 2019

