

Datasheet: ICT938

Description:	MAGIC RED™ CATHEPSIN B KIT
Name:	CATHEPSIN B
Format:	MAGIC RED™ (Red Fluorescence)
Product Type:	Kits
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Information

The Magic Red Cathepsin B kit uses a quick and easy method to analyze intracellular Cathepsin B protease activity in whole living cells.

Test Principle

Cathepsin Magic Red kits measure cathepsin B protease activity by detecting active cathepsins in whole, living cells. These kits do not work by using antibodies or as an ELISA. Instead, their methodology is based on a cell-permeable and non-cytotoxic reagent which is cleaved in the presence of cathepsins to produce a fluorescent product.

The Magic Red reagent contains a cathepsin B target sequence peptide (RR)₂ linked to a red (Cresyl Violet) fluorescent probe.

The Cathepsin Magic Red kits are suitable for cells in suspension and adherent cells from many species including mammalian, insect and yeast. Different cell types, e.g. Jurkat, HL-60, THP-1, fibroblasts, UMUC-3, MCF-7 and U937 cells have also been successfully studied with these kits.

This kit can be used in fluorescence microscopy or with fluorescence plate readers using black microtitre plates.

Reagents In The Kit

Magic Red Substrate MR-(RR)₂, 100 Tests - lyophilized
Hoechst 33342, 1 mL

Acridine Orange, 0.5 ml.

Instructions For Use	Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/ICT938.pdf .
References	1. Hämälistö, S. <i>et al.</i> (2020) Spatially and temporally defined lysosomal leakage facilitates mitotic chromosome segregation. Nat Commun. 11 (1): 229.
Storage	Store the unopened kit and each unopened component at +4°C until the expiration date. Once reconstituted with DMSO, use Magic Red reagent immediately, or store at -20°C for 6 months protected from light and thawed no more than twice during that time.
Guarantee	Guaranteed until date of expiry. Please see product label.
Acknowledgements	Magic Red is a trademark of Immunochemistry Technologies, LLC.
Health And Safety Information	Material Safety Datasheet documentation #20263 #10476 #10478 available at: https://www.bio-rad-antibodies.com/SDS/ICT938 Magic Red Substrate MR-(RR) ₂ (20263) Hoechst 33342 (10476) Acridine Orange (10478)
Regulatory	For research purposes only.

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M404833:220906'

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)