

Datasheet: ICT961

BATCH NUMBER 171571

Description:	SR FLICA™ CASPASE-9 KIT
Name:	CASPASE-9
Format:	Sulforhodamine (Red Fluorescence)
Product Type:	Kits
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunofluorescence	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Sulforhodamine	565	590

Product Information **The SR FLICA Caspase 9 kit** uses a quick and easy method to analyze active caspases in apoptotic cells.

Test Principle

Caspase FLICA kits measure apoptosis by detecting active caspases in whole, living cells. These kits do not work by using antibodies or as an ELISA. Instead, their methodology is based on a unique cell-permeable and non-cytotoxic reagent called the Fluorochrome Inhibitor of Caspases (FLICA). The FLICA reagent contains a caspase inhibitor sequence linked to a red (Sulforhodamine, SR) fluorescent probe.

The Caspase FLICA Kits are suitable for cells in suspension and adherent cells from many species including mammalian, insect and yeast. Different cell types, e.g. HeLa, primary neurons, macrophages and lymphocytes have also been successfully studied with these kits.

This kit can be used with a fluorescence microscope or a fluorescence plate reader using black microtitre plates.

Reagents In The Kit	4 vials of SR-LEHD-FMK FLICA Reagent - lyophilized 10X Wash Buffer, 60 mL Fixative, 6 mL Hoechst Stain, 1 ml.
Instructions For Use	Instructions for use can be found at www.bio-rad-antibodies.com/uploads/IFU/ICT961.pdf
Storage	This product is shipped at ambient temperature. Store the unopened kit (and each unopened component) at +4°C until the expiration date. Protect the FLICA reagent from light at all times. Once reconstituted, the 150X FLICA stock should be stored at -20°C protected from light.
Guarantee	Guaranteed until date of expiry. Please see product label.
Acknowledgements	FLICA is a trademark of Immunochemistry Technologies, LLC.
Health And Safety Information	Material Safety Datasheet documentation #20288 #10471 #10498 #10476 available at: https://www.bio-rad-antibodies.com/SDS/ICT961
Regulatory	For research purposes only.

Related Products

Recommended Useful Reagents

[SR FLICA™ CASPASE-9 KIT \(ICT960\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M439393:250523'

Printed on 29 Jan 2026