

Datasheet: LNK113PECY7 BATCH NUMBER 161249

Description:	LYNX RAPID RPE	E-Cy7 AN1		ONJUGATION KIT		
Name:	RPE-Cy7 CONJU	GATION K	KIT			
Format:	Kit					
Product Type:	Conjugation Kit					
Quantity:	1 CONJUGATION	l for 600μç	g antibody	,		
Product Details						
Applications	This product has been re derived from testing with communications from the information. For general rad-antibodies.com/proto	nin our labo e originato protocol re	oratories, ors. Please	peer-reviewed publica e refer to references in	tions or personal dicated for further	
	Conjugation	•				
Product Information	LYNX Rapid RPE-Cy7 A pre-prepared lyophilized antibody or protein. Activ results in directional cove	mixture co vation of p	ontaining roprietary	R-Phycoerythrin (RPE reagents within the ar)-Cy7 label to an htibody-label solution	
	The LYNX Rapid Conjug antibody/protein at near antibody recovery.			•		
Reagents In The Kit	1 Vial of 1mg LYNX lyop 1 Vial LYNX Modifier rea 1 Vial LYNX Quencher re	agent	PE-Cy7 mi	x		
Preparing The Antibody	The following buffer solutions are recommended for preparing the antibody:					
	10-50mM amine-free buffer (e.g HEPES, MES, MOPS and phosphate) pH range 6.5-8.5, although moderate concentrations of Tris buffer (<20mM) may be tolerated.					
	lf possible, avoid buffe	ers contai	ning nucl	eophilic component	s such as primary	

chemicals. Azide (0.02-0.1%), EDTA, up to 50% Glycerol and common non-buffering salts and sugars have little or no effect on conjugation efficiency.

	Due to the large size of the RPE-Cy7 label, it is recommended that 50-60ug of antibody be used for every 100ug RPE-Cy7, to ensure a slight RPE-Cy7 molar excess (50ug antibody gives a 1:1 Ab:RPE-Cy7 molar ratio). For optimal results the antibody should be at a concentration of 1mg/ml, with a maximum volume of 600ul and a maximum antibody amount of 600ug. Antibody at a concentration of greater than 1mg/ml requires dilution. Antibody below 1mg/ml can still be used as long as the maximum volume is not exceeded. Using less than the recommended amount of antibody may result in unbound label, but this will be removed during subsequent application wash steps. Antibody below 0.5mg/ml should be concentrated before use with the kit.
Instructions For Use	1.To the antibody sample add 1ul of the Modifier reagent for every 10ul of antibody and mix gently.
	2.Pipette the mixed antibody-modifier sample directly onto the LYNX lyophilized mix and gently pipette up and down twice to resuspend.
	3.Replace cap onto vial and incubate in the dark at room temperature (20-25 ^o C) for 3 hours, or overnight if preferred.
	4.After incubation, add 1ul of Quencher reagent for every 10ul of antibody used. Leave to stand for 30 minutes before use.
References	 Gawronska-Kozak, B. <i>et al.</i> (2021) Dermal White Adipose Tissue (dWAT) Is Regulated by Foxn1 and Hif-1α during the Early Phase of Skin Wound Healing. Int J Mol Sci. 23 (1)Dec 27 [Epub ahead of print]. Haach, V. <i>et al.</i> (2023) A polyvalent virosomal influenza vaccine induces broad cellular and humoral immunity in pigs. Virol J. 20 (1): 181. Rotolo, A. <i>et al.</i> (2023) Unedited allogeneic iNKT cells show extended persistence in MHC-mismatched canine recipients. <u>Cell Rep Med. 4 (10): 101241.</u>
Storage	This kit contains lyophilized hygroscopic components that are moisture-sensitive. This kit is shipped under ambient conditions with silica packets to avoid exposure to moisture. On receipt, Bio-Rad recommend that the kit is stored at -20°C and protected from moisture. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing. Before opening, allow the components to reach room temperature to minimize condensation.
Guarantee	12 months from date of despatch
Acknowledgements	This product or portions thereof is manufactured under license from Carnegie Mellon University under U.S. Patent Number 5,268,486 and related patents. Cy and CyDye are trademarks of GE Healthcare Limited.
Health And Safety Information	Material Safety Datasheet documentation #10551 #10546 #10549 available at: <u>https://www.bio-rad-antibodies.com/SDS/LNK113PECY7</u>

	M	rophilized RPE-Cy7 odifier Reagent (10 uencher Reagent (1	546)		
Licensed	Ki pu riç us pu da pu	ngdom patent numb irchase of this produ iht (without the right se the product to ma irchaser cannot sell ata made using this	ber 2446088 and associa uct conveys to the buyer t to resell repackage or f ake conjugates for resea or otherwise transfer th product, or its compone	ated internati r the limited, further sublic arch and deve is product, ou ints to a third	atent protected under United ional patent applications. The non exclusive non-transferable ense) under these patents to elopment purposes only. The r its components, or materials or party. Further information on obtained by contacting Bio-Rad,
					_ane, Kidlington, Oxon. OX5 antibodies@bio-rad.com
Regulato	10		OM. Tel: +44 1865 852		
Regulato North & South America	10 D ry Fo	GE UNITED KINGD or research purpose Worldwide	OM. Tel: +44 1865 852	700. E-mail: Europe	
North & South America	10 Pry Fo Tel: +1 800 265 7376 Fax: +1 919 878 375' Email: antibody_sales	GE UNITED KINGD or research purpose Worldwide	OM. Tel: +44 1865 852 s only Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio	700. E-mail: Europe	antibodies@bio-rad.com Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint