

Datasheet: LNK132APCCY7 BATCH NUMBER 166007

LYNX RAPID A					
			ONJUGATION KIT		
APC-Cy7 CON	JUGATION P	KIT			
Kit					
Conjugation Kit	t				
3 CONJUGATI					
This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-rad-antibodies.com/protocols</u> .					
Conjugation	Yes	NO	Not Determined	Suggested Dilution	
 LYNX Rapid APC-Cy7 Antibody Conjugation Kit® enables the rapid conjugation of a pre-prepared lyophilized mixture containing Allophycocyanin (APC)-Cy7 label to an antibody or protein. Activation of proprietary reagents within the antibody-label solution results in directional covalent bonding of APC-Cy7 to the antibody. The LYNX Rapid Conjugation kit® can be used to label small quantities of antibody/protein at near neutral pH, allowing a high conjugation efficiency with 100% antibody recovery. 					
The LYNX Rapid Con antibody/protein at ne	jugation kit®	can be us	sed to label small quar	ntities of	
The LYNX Rapid Con antibody/protein at ne	jugation kit® ear neutral p⊦ X lyophilized reagent	can be us I, allowing	sed to label small quar a high conjugation ef	ntities of	
The LYNX Rapid Con antibody/protein at ne antibody recovery. 3 Vials of 100ug LYNX 1 Vial LYNX Modifier 1 Vial LYNX Quenche	jugation kit® ar neutral pH X lyophilized reagent er reagent	can be us I, allowing APC-Cy7	sed to label small quar a high conjugation ef	ntities of ficiency with 100%	
The LYNX Rapid Con antibody/protein at ne antibody recovery. 3 Vials of 100ug LYNX 1 Vial LYNX Modifier 1 Vial LYNX Quenche The following buffer s 10-50mM amine-free	jugation kit® ar neutral pH X lyophilized reagent er reagent olutions are r buffer (e.g H	can be us I, allowing APC-Cy7 recommer EPES, MI	sed to label small quar a high conjugation ef mix nded for preparing the	ntities of ficiency with 100% antibody: nate) pH range 6.5-8.5,	
	Conjugation Kit 3 CONJUGATIO This product has been derived from testing w communications from information. For gene rad-antibodies.com/pr Conjugation We recommend that f ratio. LYNX Rapid APC-Cy pre-prepared lyophiliz antibody or protein. A	Conjugation Kit 3 CONJUGATIONS for 100p This product has been reported to derived from testing within our labor communications from the originator information. For general protocol re- rad-antibodies.com/protocols. Yes Conjugation • We recommend that for each conju- ratio. LYNX Rapid APC-Cy7 Antibody of pre-prepared lyophilized mixture con- antibody or protein. Activation of p	Conjugation Kit 3 CONJUGATIONS for 100µg antibod This product has been reported to work in the derived from testing within our laboratories, communications from the originators. Please information. For general protocol recomment rad-antibodies.com/protocols. Yes No Conjugation • We recommend that for each conjugation the ratio. LYNX Rapid APC-Cy7 Antibody Conjugat pre-prepared lyophilized mixture containing a	Conjugation Kit 3 CONJUGATIONS for 100µg antibody This product has been reported to work in the following application derived from testing within our laboratories, peer-reviewed publica communications from the originators. Please refer to references in information. For general protocol recommendations, please visit w rad-antibodies.com/protocols. Yes No Not Determined Conjugation • We recommend that for each conjugation the user determines the ratio. LYNX Rapid APC-Cy7 Antibody Conjugation Kit® enables the ratio.	

	and common non-buffering salts and sugars have little or no effect on conjugation efficiency.
	For optimal results the antibody should be at a concentration of 1mg/ml, with a maximum volume of 100ul and a maximum antibody amount of 100ug. Antibody at a concentration of greater than 1mg/ml requires dilution. Antibody below 1mg/ml can still be used as long as the maximum volume is not exceeded. Using less than the recommended amount of antibody may result in unbound label, but this will be removed during subsequent application wash steps. Antibody below 0.5mg/ml should be concentrated before use with the kit.
Instructions For Use	1.To the antibody sample add 1ul of the Modifier reagent for every 10ul of antibody and mix gently.
	2.Pipette the mixed antibody-modifier sample directly onto the LYNX lyophilized mix and gently pipette up and down twice to resuspend.
	3.Replace cap onto vial and incubate in the dark at room temperature (20-25 ^o C) for 3 hours, or overnight if preferred.
	4.After incubation, add 1ul of Quencher reagent for every 10ul of antibody used. Leave to stand for 30 minutes before use.
References	1. Jax, E. <i>et al.</i> (2023) Evaluating Effects of AIV Infection Status on Ducks Using a Flow Cytometry-Based Differential Blood Count <u>Microbiology Spectrum. jun 15: e0435122.</u>
Storage	This kit contains lyophilized hygroscopic components that are moisture-sensitive. This kit is shipped under ambient conditions with silica packets to avoid exposure to moisture. On receipt, Bio-Rad recommend that the kit is stored at -20°C and protected from moisture. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing. Before opening, allow the components to reach room temperature to minimize condensation.
Guarantee	12 months from date of despatch
Acknowledgements	This product or portions thereof is manufactured under license from Carnegie Mellon University under U.S. Patent Number 5,268,486 and related patents. Cy and CyDye are trademarks of GE Healthcare Limited.
Health And Safety Information	Material Safety Datasheet documentation #10553 #10546 #10549 available at: https://www.bio-rad-antibodies.com/SDS/LNK132APCCY7 Lyophilized APC-Cy7 Mix (10553) Modifier Reagent (10546) Quencher Reagent (10549)
Licensed Use	These products and the methodology of conjugation are patent protected under United Kingdom patent number 2446088 and associated international patent applications. The purchase of this product conveys to the buyer the limited, non exclusive non-transferable

right (without the right to resell repackage or further sublicense) under these patents to use the product to make conjugates for research and development purposes only. The purchaser cannot sell or otherwise transfer this product, or its components, or materials or data made using this product, or its components to a third party. Further information on purchasing licenses for diagnostic and other uses may be obtained by contacting Bio-Rad, at. Endeavour House, Langford Business Park, Langford Lane, Kidlington, Oxon. OX5 1GE UNITED KINGDOM. Tel: +44 1865 852 700. E-mail: <u>antibodies@bio-rad.com</u>

Regulatory For research purposes only

North & South	Tel: +1 800 265 7376	Worldwide	Tel: +44 (0)1865 852 700	Europe	Tel: +49 (0) 89 8090 95 21
America	Fax: +1 919 878 3751		Fax: +44 (0)1865 852 739		Fax: +49 (0) 89 8090 95 50
	Email: antibody_sales_us@bio-rad.com		Email: antibody_sales_uk@bio-rad.com		Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M390415:210910'

Printed on 26 Jun 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint