

Datasheet: MCA1037A488

## **BATCH NUMBER 161593**

Description:	RAT ANTI DOG CD5:Alexa Fluor®488		
Specificity:	CD5		
Format:	ALEXA FLUOR® 488		
Product Type:	Monoclonal Antibody		
Clone:	YKIX322.3		
Isotype:	IgG2a		
Quantity:	100 TESTS/1ml		

## **Product Details**

### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	•			Neat - 1/10

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Dog			
Purified IgG conjug	- liquid		
Fluorophore	Excitation Max (nm)	Emission Max (nm)	
Alexa Fluor®488	495	519	
supernatant		raphy on Protein G from tissue o	cultu
0.09% Sodium Azid	de (NaN <sub>3</sub> )		
IgG concentration (	0.05mg/ml		
	Purified IgG conjug  Fluorophore  Alexa Fluor®488  Purified IgG preparations supernatant  Phosphate buffered  0.09% Sodium Azid  1% Bovine Serum	Purified IgG conjugated to Alexa Fluor 488  Fluorophore Excitation Max (nm)  Alexa Fluor®488 495  Purified IgG prepared by affinity chromatoge	Purified IgG conjugated to Alexa Fluor 488 - liquid  Fluorophore Excitation Max (nm) Emission Max (nm)  Alexa Fluor®488 495 519  Purified IgG prepared by affinity chromatography on Protein G from tissue of supernatant  Phosphate buffered saline  0.09% Sodium Azide (NaN <sub>3</sub> ) 1% Bovine Serum Albumin

Immunogen	Concanavilin A activated canine peripheral blood cells
RRID	AB_10863322
Fusion Partners	Spleen cells from an immunised DA rat were fused with cells of the rat Y3/Ag1.2.3 myeloma cell line
Specificity	Rat anti Dog CD5 antibody, clone YKIX322.3 recognizes canine CD5, a 67 kDa cell surface type 1 transmembrane glycoprotein also known as lymphocyte antigen T1, Ly-1 or Leu-1. CD5 is expressed on the surface of T-cells and thymocytes, CD5 is also expressed by NK cells at low levels ( <u>Huang et al. 2008</u> ). Rat anti dog CD5, cloneYKIX322.3 was clustered as canine CD5 in the First Canine Leucocyte Antigen Workshop ( <u>Cobbold et al. 1994</u> ).
	In a study of 73 cases of canine chronic lymphocytic leukemia (CLL) CD5 expression was absent on all cases of B-cell CLL as defined by CD21 expression and lack of CD3 or other T cell antigen expression (Vernau and Moore 1999). Rat anti dog CD5 serves as a useful marker for the discrimination of canine leukemias of differing origins (Deravi et al. 2017).
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells or 100ul of whole blood.
References	<ol> <li>Cobbold, S/ &amp; Metcalfe, S. (1994) Monoclonal antibodies that define canine homologues of human CD antigens: summary of the First International Canine Leukocyte Antigen Workshop (CLAW). <u>Tissue Antigens. 43 (3): 137-54.</u></li> <li>Hewicker-Trautwein, M. <i>et al.</i> (1999) Immunocytochemical demonstration of lymphocyte subsets and MHC class II antigen expression in synovial membranes from dogs with rheumatoid arthritis and degenerative joint disease. <u>Vet Immunol Immunopathol. 67 (4): 341-57.</u></li> </ol>
	3. Huang, Y.C. (2008) CD5-low expression lymphocytes in canine peripheral blood show characteristics of natural killer cells. <u>J Leukoc Biol. 84: 1501-10.</u> 4. Araújo, M.S. <i>et al.</i> (2011) Immunological changes in canine peripheral blood leukocytes triggered by immunization with first or second generation vaccines against canine visceral leishmaniasis. <u>Vet Immunol Immunopathol. 141: 64-75.</u> 5. Burnett, R.C. <i>et al.</i> (2003) Diagnosis of canine lymphoid neoplasia using clonal rearrangements of antigen receptor genes. <u>Vet Pathol. 40: 32-41.</u> 6. Fosmire, S.P. <i>et al.</i> (2007) Inactivation of the p16 cyclin-dependent kinase inhibitor in high-grade canine non-Hodgkin's T-cell lymphoma. <u>Vet Pathol. 44: 467-78.</u> 7. Guarga, J.L. <i>et al.</i> (2002) Evaluation of a specific immunochemotherapy for the treatment of canine visceral leishmaniasis. <u>Vet Immunol Immunopathol. 88: 13-20.</u>
	0 V M M DF

- 8. Vernau, W. Moore, P.F. *et al.* (1999) An immunophenotypic study of canine leukemias and preliminary assessment of clonality by polymerase chain reaction. <u>Vet Immunol Immunopathol. 69: 145-64.</u>

  9. Lamerato-kozicki. A.R. *et al.* (2006) Canine hemangiosarcoma originates from
- 9. Lamerato-kozicki, A.R. *et al.* (2006) Canine hemangiosarcoma originates from hematopoietic precursors with potential for endothelial differentiation. <u>Exp Hematol. 34 (7):</u> 870-8.
- 10. Rütgen BC *et al.* (2012) Authentication of primordial characteristics of the CLBL-1 cell line prove the integrity of a canine B-cell lymphoma in a murine in vivo model. <u>PLoS One.</u> 7 (6): e40078.

- 11. Aresu, L. *et al.* (2014) VEGF and MMP-9: biomarkers for canine lymphoma. <u>Vet Comp</u> Oncol. 12: 29-36.
- 12. Moreira, M.L. *et al.* (2016) Vaccination against canine leishmaniosis increases the phagocytic activity, nitric oxide production and expression of cell activation/migration molecules in neutrophils and monocytes. Vet Parasitol. 220: 33-45.
- 13. Gelain, M.E. *et al.* (2014) CD44 in canine leukemia: analysis of mRNA and protein expression in peripheral blood. <u>Vet Immunol Immunopathol. 159 (1-2): 91-6.</u>
- 14. Michael, H.T. *et al.* (2013) Isolation and characterization of canine natural killer cells. Vet Immunol Immunopathol. 155 (3): 211-7.
- 15. Karayannopoulou, M. *et al.* (2017) Evaluation of blood T-lymphocyte subpopulations involved in host cellular immunity in dogs with mammary cancer. <u>Vet Immunol Immunopathol</u>. 186: 45-50.
- 16. Bonnefont-Rebeix, C. *et al.* (2016) Characterization of a novel canine T-cell line established from a spontaneously occurring aggressive T-cell lymphoma with large granular cell morphology. <a href="Immunobiology.221">Immunobiology.221</a> (1): 12-22.
- 17. Gibbons, N. *et al.* (2017) Phenotypic heterogeneity of peripheral monocytes in healthy dogs. Vet Immunol Immunopathol. 190: 26-30.
- 18. Deravi, N. *et al.* (2017) Specific immunotypes of canine T cell lymphoma are associated with different outcomes. <u>Vet Immunol Immunopathol. 191: 5-13.</u>
- 19. GomesMde, O. *et al.* (2011) Old beagle dogs have lower faecal concentrations of some fermentation products and lower peripheral lymphocyte counts than young adult beagles. <u>Br J Nutr. 106 Suppl 1: S187-90.</u>
- 20. MariaA, P.J. *et al.* (2017) The effect of age and carbohydrate and protein sources on digestibility, fecal microbiota, fermentation products, fecal IgA, and immunological blood parameters in dogs. <u>J Anim Sci. 95 (6): 2452-66.</u>
- 21. Stokol, T. *et al.* (2015) Alkaline phosphatase is a useful cytochemical marker for the diagnosis of acute myelomonocytic and monocytic leukemia in the dog. <u>Vet Clin Pathol.</u> 44 (1): 79-93.
- 22. Martini, V. *et al.* (2019) Prognostic role of non-neoplastic lymphocytes in lymph node aspirates from dogs with diffuse large B-cell lymphoma treated with chemo-immunotherapy. Res Vet Sci. 125: 130-5.
- 23. Lin, C.S. *et al.* (2018) Activating natural killer (NK) cytotoxicity of canine CD5<sup>-</sup>CD21<sup>-</sup> cells requires low surface CD5 density NK cells. <u>Iran J Vet Res. 19 (2): 87-95.</u>
- 24. Roatt, B.M. *et al.* (2017) A Vaccine Therapy for Canine Visceral Leishmaniasis Promoted Significant Improvement of Clinical and Immune Status with Reduction in Parasite Burden. <u>Front Immunol. 8: 217.</u>
- 25. Aricò, A. *et al.* (2013) The role of vascular endothelial growth factor and matrix metalloproteinases in canine lymphoma: *in vivo* and *in vitro* study. <u>BMC Vet Res. 9: 94.</u>
  26. Aguiar-Soares, R.D.O. *et al.* (2020) Phase I and II Clinical Trial Comparing the LBSap,
- Leishmune<sup>®</sup>, and Leish-Tec<sup>®</sup> Vaccines against Canine Visceral Leishmaniasis. <u>Vaccines</u> (Basel). 8 (4)Nov 17 [Epub ahead of print].
- 27. Graves, S.S. *et al.* (2019) Development and characterization of a canine-specific anti-CD94 (KLRD-1) monoclonal antibody. <u>Vet Immunol Immunopathol. 211: 10-8.</u>
  28. Ito, D. *et al.* (2015) A double blinded, placebo-controlled pilot study to examine reduction of CD34 <sup>+</sup>/CD117 <sup>+</sup>/CD133 <sup>+</sup> lymphoma progenitor cells and duration of remission induced by neoadjuvant valspodar in dogs with large B-cell lymphoma. <u>F1000Res. 4: 42.</u>

- 29. Wolf-Ringwall, A. *et al.* (2020) Prospective evaluation of flow cytometric characteristics, histopathologic diagnosis and clinical outcome in dogs with naïve B-cell lymphoma treated with a 19-week CHOP protocol. <u>Vet Comp Oncol. 18 (3): 342-52.</u>
  30. Sayag, D. *et al.* (2020) Proof-concept study: Evaluation of plasma and urinary
- 30. Sayag, D. *et al.* (2020) Proof-of-concept study: Evaluation of plasma and urinary electrolytes as markers of response to L-asparaginase therapy in dogs with high-grade lymphoma. Vet Clin Pathol. 49 (3): 476-83.
- 31. Lee, J. *et al.* (2021) Canine Natural Killer Cell-Derived Exosomes Exhibit Antitumor Activity in a Mouse Model of Canine Mammary Tumor. <u>Biomed Res Int. 2021: 6690704</u>.
- 32. Grudzien, M. *et al.* (2021) A newly established canine NK-type cell line and its cytotoxic properties. Vet Comp Oncol. 19 (3): 567-77.
- 33. Lee, S.H. *et al.* (2021) Safety and immunological effects of recombinant canine IL-15 in dogs. Cytokine. 148: 155599.
- 34. Karayannopoulou, M. *et al.* (2022) Effect of major versus minor mastectomy on host immunity in canine mammary cancer <u>Vet Immunol Immunopathol. 24 Feb: 110403.</u>
- 35. Riccardo, F. *et al.* (2022) Antigen mimicry as an effective strategy to induce CSPG4-targeted immunity in dogs with oral melanoma: a veterinary trial. <u>J Immunother Cancer. 10</u> (5): e004007. [Epub ahead of print].
- 36. Jaensch, S. *et al.* (2022) Clinicopathologic and immunophenotypic features in dogs with presumptive large granular lymphocyte leukaemia <u>Australian Veterinary Journal.</u> [Epub ahead of print].

#### Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

## Guarantee

12 months from date of despatch

### Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

# Health And Safety Information

Material Safety Datasheet documentation #10041 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA1037A488">https://www.bio-rad-antibodies.com/SDS/MCA1037A488</a> 10041

# Regulatory

For research purposes only

# **Related Products**

# **Recommended Negative Controls**

## RAT IgG2a NEGATIVE CONTROL:Alexa Fluor® 488 (MCA1212A488)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M385295:210513'

### Printed on 05 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint