

Datasheet: MCA1076SBUV605

BATCH NUMBER 100006783

| Description: | MOUSE ANTI HUMAN CD62L:StarBright UltraViolet 605 |
|---------------|---|
| Specificity: | CD62L |
| Other names: | LECAM-1, L-SELECTIN |
| Format: | StarBright UltraViolet 605 |
| Product Type: | Monoclonal Antibody |
| Clone: | FMC46 |
| Isotype: | IgG2b |
| Quantity: | 100 TESTS/0.5ml |
| | |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | • | | | Neat |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| Species Cross | |
|---------------|--|
| B | |
| Reactivity | |

Target Species

Reacts with: Bovine, Cynomolgus monkey, Rhesus Monkey, Dog

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

turther information

Human

Product Form

Purified IgG conjugated to StarBright UltraViolet 605 - liquid

| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
|-----------|-------------------------------|---------------------|-------------------|
| | StarBright UltraViolet 605 | 340 | 609 |

Preparation

Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant

| Buffer Solution | Phosphate buffered saline | |
|-----------------------------|--|---|
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 0.1% Pluronic F68 0.1% PEG 3350 0.05% Tween 20 | |
| lmmunogen | PHA stimulated lymphoblasts | |
| External Database Links | UniProt: P14151 Related reagents Entrez Gene: 6402 SELL Related reagents | |
| Synonyms | LNHR, LYAM1 | |
| Fusion Partners | Spleen cells from immunized BALB/c mice were fused with ce myeloma cell line | lls of the mouse NS1 |
| Specificity | Mouse anti Human CD62L antibody, clone FMC46 recognizione known a L-selectin, a 74-95 kDa member of the selectin family which acts as a ligand for both CD62P (P-selectin) and CD62E CD62L is constitutively expressed on most leucocytes includin lymphocytes, NK cells, bone marrow myeloid progenitor cells at thymocytes. | of adhesion receptors, E (E-selectin). Human g monocytes, granulocytes, |
| | CD62L plays an important role in leucocyte tethering and rolling surface and for the homing of naïve lymphocytes to lymph noon HEV. Neutrophils require a constant supply of this molecule or migration into peripheral tissues and adhesion to activated end inflammation, where CD62L is rapidly shed as soluble L-select still remains. | les and Peyers patches via the cell surface for dothelium at sites of |
| | The expression of CD62L is down regulated on lymphocytes a stimulation. | and neutrophils by PMA |
| Flow Cytometry | Use 5ul of the suggested working dilution to label 10 ⁶ cells in suggest a 5 minutes centrifugation at 6,000g prior to sample a | • |
| References | 1. Pilarski, L.M. <i>et al.</i> (1991) FMC46, a cell protrusion-associal molecule-1 epitope on human lymphocytes and thymocytes. <u>J</u> 2. Zola, H. <i>et al.</i> (1991) The expression of sub-population mar re-evaluation using high-sensitivity fluorescence flow cytometric 103-18. 3. Sopp, P. & Howard, C.J. (1997) Cross-reactivity of monoclo | Immunol. 147 (1): 136-43. kers on B cells: a y. <u>Dis Markers</u> . 9 (2): |

human leucocyte differentiation antigens with bovine cells. <u>Vet Immunol Immunopathol. 56</u> (1-2): 11-25.

- 4. Haanstra, K.G. *et al.* (2008) Characterization of naturally occurring CD4+CD25+ regulatory T cells in rhesus monkeys. <u>Transplantation 85:1185-92.</u>
- 5. Dalli, J. *et al.* (2008) Annexin 1 mediates the rapid anti-inflammatory effects of neutrophil-derived microparticles. <u>Blood. 112 (6): 2512-9.</u>
- 6. Raposo, R.A. *et al.* (2011) Protein Kinase C and NF-{kappa}B-Dependent CD4 Downregulation in Macrophages Induced by T Cell-Derived Soluble Factors: Consequences for HIV-1 Infection. <u>J Immunol.</u> 187: 748-59.
- 7. Hughes, S.F. *et al.* (2010) Total hip and knee replacement surgery results in changes in leukocyte and endothelial markers. J Inflamm (Lond). 7:2.
- 8. Bismarck, D. *et al.* (2012) Canine CD4+CD8+ double positive T cells in peripheral blood have features of activated T cells. Vet Immunol Immunopathol. 149: 157-66.
- 9. Hartley, A.N. & Tarleton, R.L. (2015) Chemokine receptor 7 (CCR7)-expression and IFNγ production define vaccine-specific canine T-cell subsets. <u>Vet Immunol Immunopathol.</u> 164 (3-4): 127-36.
- 10. Hayhoe, R.P. *et al.* (2006) Annexin 1 and its bioactive peptide inhibit neutrophilendothelium interactions under flow: indication of distinct receptor involvement. <u>Blood. 107</u> (5): 2123-30.
- 11. Urquhart, P. et al. (2007) Carbon monoxide-releasing molecules modulate leukocyte-endothelial interactions under flow. <u>J Pharmacol Exp Ther. 321 (2): 656-62.</u>
- 12. Aspinall, A.I. *et al.* (2010) CX(3)CR1 and vascular adhesion protein-1-dependent recruitment of CD16(+) monocytes across human liver sinusoidal endothelium. <u>Hepatology. 51 (6): 2030-9.</u>
- 13. Rothe, K. *et al.* (2017) Canine peripheral blood CD4⁺CD8⁺ double-positive Tcell subpopulations exhibit distinct Tcell phenotypes and effector functions. <u>Vet Immunol Immunopathol</u>. 185: 48-56.
- 14. Withers, S.S. *et al.* (2018) Multi-color flow cytometry for evaluating age-related changes in memory lymphocyte subsets in dogs. <u>Dev Comp Immunol. 87: 64-74.</u>
- 15. Hughes, S.F. *et al.* (2020) The role of phagocytic leukocytes following flexible ureterenoscopy, for the treatment of kidney stones: an observational, clinical pilots-study. <u>Eur J Med Res. 25 (1): 68.</u>
- 16. Svitek, N. *et al.* (2018) An Ad/MVA vectored *Theileria parva* antigen induces schizont-specific CD8⁺ central memory T cells and confers partial protection against a lethal challenge. NPJ Vaccines. 3: 35.

| Storage | Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. | |
|----------------------------------|--|--------------------|
| Guarantee | 12 months from date of despatch | |
| Acknowledgements | This product is covered by U.S. Patent No. 10,150,841 and related counterparts | d U.S. and foreign |
| Health And Safety Information | Material Safety Datasheet documentation #20471 available at: https://www.bio-rad-antibodies.com/SDS/MCA1076SBUV605 20471 | |

Related Products

Recommended Useful Reagents

<u>HUMAN SEROBLOCK (BUF070A)</u> <u>HUMAN SEROBLOCK (BUF070B)</u>

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M403683:220728'

Printed on 08 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint