

Datasheet: MCA1133A647

Description:	MOUSE ANTI HUMAN CD49a:Alexa Fluor® 647
Specificity:	CD49a
Other names:	INTEGRIN ALPHA 1 CHAIN, VLA-1
Format:	ALEXA FLUOR® 647
Product Type:	Monoclonal Antibody
Clone:	TS2/7
Isotype:	IgG1
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Alexa Fluor® 647 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®647	650	665
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% sodium azide (NaN ₃)		
Stabilisers	1% bovine serum albumin		
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml		

Immunogen Human CTL line directed against HLA-DR.

External Database

Links

UniProt:

[P56199](#) [Related reagents](#)

Entrez Gene:

[3672](#) ITGA1 [Related reagents](#)

RRID

AB_566788

Fusion Partners

Spleen cells from immunized Balb.c mice were fused with cells of the NS1 mouse myeloma.

Specificity

Mouse anti Human CD49a antibody, clone TS2/7 recognizes the human alpha 1 integrin sub-unit, which forms the VLA-1 heterodimer in association with the beta 1 integrin. VLA-1 is a receptor for collagen and laminin, and is expressed by chronically activated T cells, melanoma cells and smooth muscle cells.

Mouse anti Human CD49a antibody, clone TS2/7 has been reported as being suitable for use in Western Blotting.

Mouse anti Human CD49a antibody, clone TS2/7 is routinely tested in flow cytometry using HeLa cells.

Flow Cytometry

Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

References

1. Hemler, M.E. *et al.* (1985) Biochemical characterization of VLA-1 and VLA-2. Cell surface heterodimers on activated T cells. [J Biol Chem. 260 \(28\): 15246-52.](#)
2. Hafler, D.A. *et al.* (1985) Investigation of in vivo activated T cells in multiple sclerosis and inflammatory central nervous system diseases. [Clin Immunol Immunopathol. 37 \(2\): 163-71.](#)
3. Klemmt, P.A. *et al.* (2007) Endometrial cells from women with endometriosis have increased adhesion and proliferative capacity in response to extracellular matrix components: towards a mechanistic model for endometriosis progression. [Hum Reprod. 22: 3139-47.](#)
4. Berry, P.A. *et al.* (2011) The calcium sensor STIM1 is regulated by androgens in prostate stromal cells. [Prostate. 71: 1646-55.](#)
5. Elshaw, S.R. *et al.* (2001) A comparison of ocular melanocyte and uveal melanoma cell invasion and the implication of alpha1beta1, alpha4beta1 and alpha6beta1 integrins. [Br J Ophthalmol. 85: 732-8.](#)
6. Cervar-Zivkovic, M. *et al.* (2011) Endothelin-1 Stimulates Proliferation of First-Trimester Trophoblasts via the A- and B-Type Receptor and Invasion via the B-Type Receptor. [J Clin Endocrinol Metab. 96: 3408-15.](#)
7. Kaneko, Y. *et al.* (2012) Integrin α1/β1 and α2/β1 as a receptor for IgA1 in human glomerular mesangial cells in IgA nephropathy. [Int Immunol. 24: 219-32.](#)
8. Majali-Martinez, A. *et al.* (2018) Temporal changes of the endothelin system in human cytotrophoblasts during the first trimester of pregnancy. [Physiol Res. 67 \(Suppl 1\):](#)

[S247-S255.](#)

9. Zhuang, L. *et al.* (2021) Dermal Pericytes Exhibit Declined Ability to Promote Human Skin Regeneration with Ageing in 3D Organotypic Culture Models [Cells. 10 \(11\): 3051.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee 12 months from date of despatch

Acknowledgements This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information Material Safety Datasheet documentation #10041 available at: <https://www.bio-rad-antibodies.com/SDS/MCA1133A647>
10041

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:Alexa Fluor® 647 \(MCA928A647\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M408265:221011'

Printed on 18 Jan 2024