

# Datasheet: MCA1396A BATCH NUMBER 152743

Description:	MOUSE ANTI HISTIDINE TAG:ALK.PHOS.		
Specificity:	HISTIDINE TAG		
Format:	Alk. Phos.		
Product Type:	Monoclonal Antibody		
Clone:	AD1.1.10		
Isotype:	lgG1		
Quantity:	0.1 mg		

## **Product Details**

#### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry				
Immunohistology - Frozen			•	
Immunohistology - Paraffin				
ELISA	•			1/20 - 1/1000
Immunoprecipitation				
Western Blotting	-			1/1000 - 1/5000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Synthetic Peptide	
Product Form	Purified IgG conjugated to Alkaline Phosphatase - liquid	
Preparation	Purified IgG prepared by affinity chromatography on Protein G supernatant	from tissue culture
Buffer Solution	Phosphate buffered saline	
Preservative Stabilisers	0.09% Sodium Azide (NaN <sub>3</sub> )	
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml	

Immunogen	PAX6 transcription factor linked to histidine tag.
RRID	AB_1605025
Fusion Partners	Spleen cells from immunised Balb/c mice were fused with cells of the mouse NS1 myeloma cell line.
Specificity	Mouse anti Histidine tag antibody, clone AD1.1.10 recognizes proteins and peptides containing the motif H-H-H-H-H and is therefore of value in detecting proteins containing histidine tags. Clone AD1.1.10 has been used to detect and purify histidine-tagged proteins expressed in mammalian (Hoffmann <i>et al.</i> 2007) and Hwang <i>et al.</i> 2008) and non-mammalian (Zheng <i>et al.</i> 2007; Gunnarsen <i>et al.</i> 2010; and de Vooght <i>et al.</i> 2012) cell lines.
	In Western blotting of bacterial extracts the antibody has been shown not to cross-react with any endogenous products, although some cross-reactivity may be seen with extracts of insect or mammalian cells.
	This antibody is routinely tested in Western blotting on histidine tagged recombinant proteins and reacts against all histidine-tagged proteins so far tested.

#### References

- 1. Els Conrath, K. *et al.* (2001) Camel single-domain antibodies as modular building units in bispecific and bivalent antibody constructs. <u>J Biol Chem. 276 (10): 7346-50.</u>
- 2. Suen, J.L. *et al.* (2001) Characterization of self-T-cell response and antigenic determinant of U1A protein with bone marrow-derived dendritic cells in NZB x NZW F<sub>1</sub> mice. Immunol. 103: 301-309.
- 3. Hoffmann, S.C. *et al.* (2007) Identification of CLEC12B, an inhibitory receptor on myeloid cells. <u>J Biol Chem. 282 (31): 22370-5.</u>
- 4. Zheng, J. *et al.* (2007) Serum from mice immunized in the context of Treg inhibition identifies DEK as a neuroblastoma tumor antigen. <u>BMC Immunol. 8: 4.</u>
- 5. Bahi, A. & Dreyer, J.L. (2008) Overexpression of plasminogen activators in the nucleus accumbens enhances cocaine-, amphetamine- and morphine-induced reward and behavioral sensitization. <u>Genes Brain Behav. 7 (2): 244-56.</u>
- 6. Wrighton, K.H. *et al.* (2009) Transforming Growth Factor {beta} Can Stimulate Smad1 Phosphorylation Independently of Bone Morphogenic Protein Receptors. <u>J Biol Chem. 284</u> (15): 9755-63.
- 7. Diefenbacher, M. *et al.* (2011) The Dsl1 Tethering Complex Actively Participates in Soluble NSF (N-Ethylmaleimide-sensitive Factor) Attachment Protein Receptor (SNARE) Complex Assembly at the Endoplasmic Reticulum in Saccharomyces cerevisiae. <u>J Biol Chem.</u> 286: 25027-38.
- 8. Alvarez, M.M. *et al.* (2010) Specific recognition of influenza A/H1N1/2009 antibodies in human serum: a simple virus-free ELISA method. <u>PLoS One. 5: e10176.</u>
- 9. Bahi, A. *et al.* (2008) The role of tissue-type plasminogen activator system in amphetamine-induced conditional place preference extinction and reinstatement. Neuropsychopharmacology. 33: 2726-34.
- 10. Gunnarsen, K.S. *et al.* (2010) Periplasmic expression of soluble single chain T cell receptors is rescued by the chaperone FkpA. <u>BMC Biotechnol.</u> 10: 8.

- 11. Hwang, H.Y. *et al.* (2008) Highly specific inhibition of C1q globular-head binding to human IgG: a novel approach to control and regulate the classical complement pathway using an engineered single chain antibody variable fragment. <u>Mol Immunol. 45: 2570-80.</u>
- 12. De Vooght, L. *et al.* (2012) Expression and extracellular release of a functional anti-trypanosome Nanobody® in Sodalis glossinidius, a bacterial symbiont of the tsetse fly. <u>Microb Cell Fact. 11: 23.</u>
- 13. Saerens, D. *et al.* (2004) Single domain antibodies derived from dromedary lymph node and peripheral blood lymphocytes sensing conformational variants of prostate-specific antigen. J Biol Chem. 279 (50): 51965-72.
- 14. Than, N.G. *et al.* (2014) Evolutionary origins of the placental expression of chromosome 19 cluster galectins and their complex dysregulation in preeclampsia. Placenta. 35: 855-65.
- 15. Elders RC *et al.* (2014) Recombinant canine IgE Fc and an IgE Fc-TRAIL fusion protein bind to neoplastic canine mast cells. <u>Vet Immunol Immunopathol. 159 (1-2): 29-40.</u>
  16. Chin, S.E. *et al.* (2015) Isolation of high-affinity, neutralizing anti-idiotype antibodies by phage and ribosome display for application in immunogenicity and pharmacokinetic analyses. J Immunol Methods. 416: 49-58.
- 17. Peyrassol, X. *et al.* (2016) Development by Genetic Immunization of Monovalent Antibodies (Nanobodies) Behaving as Antagonists of the Human ChemR23 Receptor. <u>J. Immunol.</u> 196 (6): 2893-901.
- 18. Kim H & Loparo JJ (2016) Multistep assembly of DNA condensation clusters by SMC. Nat Commun. 7: 10200.
- 19. Borg M *et al.* (2014) A novel interaction between Rab7b and actomyosin reveals a dual role in intracellular transport and cell migration. <u>J Cell Sci. 127 (Pt 22): 4927-39.</u>
- 20. De Meyer, T. *et al.* (2015) Comparison of VHH-Fc antibody production in Arabidopsis thaliana, Nicotiana benthamiana and Pichia pastoris. <u>Plant Biotechnol J. 13 (7): 938-47.</u>
- 21. Siddiqui AA *et al.* (2015) Humoral immune responses to a recombinant Plasmodium vivax tryptophan-rich antigen among *Plasmodium vivax*-infected patients and its localization in the parasite. Appl Biochem Biotechnol. 175 (4): 2166-77.
- 22. Warnecke, A. *et al.* (2017) Nitration of MOG diminishes its encephalitogenicity depending on MHC haplotype. J Neuroimmunol. 303: 1-12.
- 23. Bertucci, A. *et al.* (2011) A new coral carbonic anhydrase in *Stylophora pistillata*. Mar Biotechnol (NY). 13 (5): 992-1002.
- 24. Liu, C.C. *et al.* (2016) The Fimbrial Protein is a Virulence Factor and Potential Vaccine Antigen of *Avibacterium paragallinarum*. Avian Dis. 60 (3): 649-55.
- 25. Boujon, C.L. *et al.* (2017) Development and validation of an immunohistochemistry procedure for the detection of a neurotropic bovine astrovirus. <u>J Virol Methods. 239:</u> 26-33.
- 26. Cartwright, S.P. *et al.* (2017) Rapid expression and purification of the hepatitis delta virus antigen using the methylotropic yeast *Pichia pastoris*. <u>BMC Res Notes</u>. 10 (1): 340.
- 27. Thanongsaksrikul, J. *et al.* (2018) Identification and production of mouse scFv to specific epitope of enterovirus-71 virion protein-2 (VP2). <u>Arch Virol. 163 (5): 1141-1152.</u>
- 28. Gunnarsen, K.S. *et al.* (2018) Soluble T-cell receptor design influences functional yield in an E. coli chaperone-assisted expression system. <u>PLoS One. 13 (4): e0195868.</u>
- 29. Ascione, A. *et al.* (2019) Development of a novel human phage display-derived anti-LAG3 scFv antibody targeting CD8<sup>+</sup> T lymphocyte exhaustion. <u>BMC Biotechnol. 19</u> (1): 67.

- 30. Zoccola, D. *et al.* (2017) Structural and functional analysis of coral Hypoxia Inducible Factor. PLoS One. 12 (11): e0186262.
- 31. Kimura, K. *et al.* (2021) Overexpression of human BAG3<sup>P209L</sup> in mice causes restrictive cardiomyopathy. <u>Nat Commun. 12 (1): 3575.</u>
- 32. Dongdem, J.T. *et al.* (2021) Modification of small ubiquitin-related modifier 2 (SUMO2) by phosphoubiquitin in HEK293T cells. <u>Proteomics. 21 (15): e2000234.</u>
- 33. Chuang, H.C. *et al.* (2021) Effect of cell-permeable grouper Manganese Superoxide Dismutase on environmental stress in fish. <u>Protein Expr Purif. 187: 105951.</u>
- 34. Cheng, C.M. *et al.* (2021) Heterologous expression of bacterial CotA-laccase, characterization and its application for biodegradation of malachite green. <u>Bioresour Technol.</u> 340: 125708.
- 35. De Vooght, L. *et al.* (2022) Targeting the tsetse-trypanosome interplay using genetically engineered *Sodalis glossinidius*.. <u>PLoS Pathog. 18 (3): e1010376.</u>
- 36. Minami, S.A. *et al.* (2022) Production of novel SARS-CoV-2 Spike truncations in Chinese hamster ovary cells leads to high expression and binding to antibodies <u>Biotechnol</u> J. Jun 3 [Epub ahead of print]

Storage	Store at +4°C. DO NOT FREEZE.
	This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	12 months from date of despatch
Acknowledgements	His-tag is a registered trademark of EMD Biosciences.
Health And Safety	Material Safety Datasheet documentation #10040 available at:
Information	https://www.bio-rad-antibodies.com/SDS/MCA1396A
	10040
Regulatory	For research purposes only

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M365245:200529'

### Printed on 15 Apr 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint