

Datasheet: MCA1749A647

BATCH NUMBER 152580

Description:	MOUSE ANTI PIG CD4 ALPHA:Alexa Fluor® 647			
Specificity:	CD4 ALPHA			
Other names:	CD4			
Format:	ALEXA FLUOR® 647			
Product Type:	Monoclonal Antibody			
Clone:	MIL17			
Isotype:	lgG2b			
Quantity:	100 TESTS/1ml			

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	•			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Pig		
Product Form	Purified IgG conjug	ated to Alexa Fluor 647	- liquid
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®647	650	665
Preparation	Purified IgG prepar supernatant	ed by affinity chromatog	raphy on Protein A f
Buffer Solution	Phosphate buffered	d saline	
Preservative	0.09% Sodium Azio	de (NaN ₃)	
Stabilisers	1% Bovine Serum	Albumin	
Approx. Protein	IgG concentration (0.05 mg/ml	

Immunogen

Leucocytes isolated from porcine gut lamina propria.

Specificity

Mouse anti Porcine CD4 alpha, clone MIL17 recognizes a ~55 kDa porcine homologue to the human CD4 antigen found on the surface of helper-T cells. MIL-17 was confirmed as a member of the CD4 alpha cluster at the 'Third International Workshop on Swine Leukocyte Differentiation Antigens' (<u>Haverson et al. 2001</u>). Porcine CD4 is a type 1 transmembrane member of the immunoglobulin superfamily.

Pigs appear unusual amongst mammalian species as they appear to have four populations of resting T lymphocytes. In addition to the two populations of mutually exclusive CD4+/CD8- and CD4-/CD8+ lymphocytes, they also appear to have significant populations of CD4-/CD8- and CD4+/CD8+ cells. Lymphoblasts with a double positive phenotype have been described in other species but this is not the case for mature T lymphocytic calls (<u>Saalmuller et al. 1987</u>)

Mouse anti Pig CD4 alpha, clone MIL17 stains a population of cells with characteristic lymphocyte morphology in immunohistochemistry (Inman *et al.* 2010).

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

- 1. Saalmüller A *et al.* (2001) Summary of workshop findings for porcine T-lymphocyte-specific monoclonal antibodies. <u>Vet Immunol Immunopathol.</u> 80 (1-2): 35-52.
- 2. Castellano, G. *et al.* (2010) Therapeutic targeting of classical and lectin pathways of complement protects from ischemia-reperfusion-induced renal damage. <u>Am J Pathol. 176:</u> 1648-59.
- 3. Inman, C.F. *et al.* (2010) Dendritic cells interact with CD4 T cells in intestinal mucosa. <u>J</u> Leukoc Biol. 88 (3): 571-8.
- 4. Kick AR *et al.* (2011) Evaluation of peripheral lymphocytes after weaning and vaccination for *Mycoplasma hyopneumoniae*. Res Vet Sci. 91 (3): e68-72.
- 5. Kick, A.R. *et al.* (2012) Effects of stress associated with weaning on the adaptive immune system in pigs. J Anim Sci. 90: 649-56.
- 6. Goujon, J.M. *et al.* (2000) Influence of cold-storage conditions on renal function of autotransplanted large pig kidneys. Kidney Int. 58: 838-50.
- 7. Tambuyzer BR *et al.* (2012) Osteopontin alters the functional profile of porcine microglia *in vitro*. Cell Biol Int. 36 (12): 1233-8.
- 8. Tuchscherer, M. *et al.* (2012) Effects of inadequate maternal dietary protein:carbohydrate ratios during pregnancy on offspring immunity in pigs. <u>BMC Vet Res.</u> 8: 232.
- 9. Cao, D. *et al.* (2010) Synthetic innate defence regulator peptide enhances in vivo immunostimulatory effects of CpG-ODN in newborn piglets. <u>Vaccine</u>. 28: 6006-13.
- 10. Clapperton, M. *et al.* (2005) Associations of weight gain and food intake with leukocyte sub-sets in Large White pigs <u>Livestock Production Science</u> 96: 249-60
- 11. Clapperton, M. *et al.* (2005) Innate immune traits differ between Meishan and Large White pigs. Vet Immunol Immunopathol. 104: 131-44.
- 12. Clapperton, M. *et al.* (2008) Pig peripheral blood mononuclear leucocyte subsets are heritable and genetically correlated with performance. <u>Animal. 2: 1575-84.</u>

- 13. Faure, J.P. *et al.* (2002) Polyethylene glycol reduces early and long-term cold ischemia-reperfusion and renal medulla injury. <u>J Pharmacol Exp Ther. 2002</u> Sep;302(3):861-70.
- 14. Faure, J.P. *et al.* (2004) Evidence for protective roles of polyethylene glycol plus high sodium solution and trimetazidine against consequences of renal medulla ischaemia during cold preservation and reperfusion in a pig kidney model. <u>Nephrol Dial Transplant.</u> 19: 1742-51.
- 15. Inman, C.F. *et al.* (2012) Neonatal colonisation expands a specific intestinal antigen-presenting cell subset prior to CD4 T-cell expansion, without altering T-cell repertoire. PLoS One. 7(3): e33707.
- 16. Kick, A.R. *et al.* (2012) Effects of stress associated with weaning on the adaptive immune system in pigs. <u>J Anim Sci. 90: 649-56.</u>
- 17. Langerhuus, S.N. *et al.* (2010) Brief report: biomarkers of aortic vascular prosthetic graft infection in a porcine model with Staphylococcus aureus. <u>Eur J Clin Microbiol Infect Dis.</u> 29: 1453-6.
- 18. Lu, X. *et al.* (2012) Genome-wide association study for T lymphocyte subpopulations in swine. BMC Genomics. 13: 488.
- 19. Monroy-Salazar, H.G. *et al.* (2012) Effects of a live yeast dietary supplement on fecal coliform counts and on peripheral blood CD4+ and CD8+ lymphocyte subpopulations in nursery pigs. J Swine Health Prod 20: 276-282.
- 20. Shi, K. *et al.* (2008) Changes in peripheral blood leukocyte subpopulations in piglets co-infected experimentally with porcine reproductive and respiratory syndrome virus and porcine circovirus type 2. <u>Vet Microbiol. 129: 367-77.</u>
- 21. Spreeuwenberg, M.A. *et al.* (2001) Small intestine epithelial barrier function is compromised in pigs with low feed intake at weaning. J Nutr. 131: 1520-7.
- 22. Tambuyzer, B.R. *et al.* (2012) Osteopontin alters the functional profile of porcine microglia in vitro. Cell Biol Int. 36: 1233-8.
- 23. Zelnickova, P. *et al.* (2007) Intracellular cytokine detection by flow cytometry in pigs: fixation, permeabilization and cell surface staining. J Immunol Methods. 327: 18-29.
- 24. Kvist, P.H. *et al.* (2010) Effect of subcutaneous glucose sensor implantation on skin mRNA expression in pigs. <u>Diabetes Technol Ther. 12: 791-9.</u>
- 25. Lefevre, E.A. *et al.* (2012) Immune responses in pigs vaccinated with adjuvanted and non-adjuvanted A(H1N1)pdm/09 influenza vaccines used in human immunization programmes. PLoS One. 7(3): e32400.
- 26. Akershoek, J.J. *et al.* (2016) Cell therapy for full-thickness wounds: are fetal dermal cells a potential source? Cell Tissue Res. 364 (1): 83-94.
- 27. Liu J *et al.* (2016) The Role of Porcine Monocyte Derived Dendritic Cells (MoDC) in the Inflammation Storm Caused by *Streptococcus suis* Serotype 2 Infection. <u>PLoS One.</u> 11 (3): e0151256.
- 28. Liermann, W. *et al.* (2017) Effects of two commercial diets and technical feed treatment on stomach lesions and immune system of fattening pigs. <u>J Anim Physiol Anim Nutr (Berl)</u>. 101 (5): e414-26.
- 29. Gardner, D.S. *et al.* (2016) Remote effects of acute kidney injury in a porcine model. Am J Physiol Renal Physiol. 310 (4): F259-71.
- 30. Hemmink, J.D. *et al.* (2016) Distinct immune responses and virus shedding in pigs following aerosol, intra-nasal and contact infection with pandemic swine influenza A virus, A(H1N1)09. <u>Vet Res. 47 (1): 103.</u>

31. Dąbrowski, M. et al. (2017) The Effect of Deoxynivalenol on Selected Populations of
Immunocompetent Cells in Porcine Blood-A Preliminary Study. Molecules. 22 (5)Apr 26
[Epub ahead of print].
32. Hsu, W.T. et al. (2013) Prostaglandin E2 potentiates mesenchymal stem cell-induced
IL-10+IFN-γ+CD4+ regulatory T cells to control transplant arteriosclerosis. <u>J Immunol. 190</u>
<u>(5): 2372-80.</u>
33. López, E. et al. (2019) Identification of very early inflammatory markers in a porcine
myocardial infarction model. BMC Vet Res. 15 (1): 91.
34. Hu, Z. et al. (2019) Genomic variant in porcine TNFRSF1A gene and its effects on
TNF signaling pathway in vitro. Gene. 700: 105-9.

35. Fogle, J.E. *et al.* (2019) Antibiotic Therapy Does Not Alter the Humoral Response to Vaccination for Porcine Circovirus 2 in Weaned Pigs. <u>Vet Sci. 6 (2)May 30 [Epub ahead of print].</u>

Further Reading

1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. Vet Res. 39: 54.

Storage

Store at +4°C or at -20°C if preferred.

Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information

Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA1749A647 10041

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG2b NEGATIVE CONTROL:Alexa Fluor® 647 (MCA691A647)

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

North & South Tel: +1 800 265 7376 America Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739

Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

'M379506:210325'

Email: antibody_sales_de@bio-rad.com

Printed on 29 Feb 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint