

Datasheet: MCA1969F

Description:	MOUSE ANTI HUMAN CD62E:FITC
Specificity:	CD62E
Other names:	E-SELECTIN
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	CL2/6
Isotype:	lgG2a
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.biorad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	•			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjug	gated to Fluorescein Isoth	niocyanate Isomer 1
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation Buffer Solution	Antibody purified fr	rom tissue culture supern	atant
Preservative	0.02% sodium azio	de (NaN ₃)	
Stabilisers	1% bovine serum albumin		
External Database	IIn:Duct:		
Links	UniProt:		

Related reagents

P16581

Entrez Gene:

6401 SELE Related reagents

Synonyms	ELAM1					
RRID	AB_323219					
Specificity	Mouse anti Human CD62E antibody, clone CL2/6 recognizes human CD62E, a member of the selectin family which is expressed on endothelial cells at sites of inflammation.					
	Mouse anti Human CD62E antibody, clone CL2/6 has been shown to inhibit the binding of human neutrophils to stimulated HUVECS.					
Flow Cytometry	Use 10μl of the suggested working dilution to label 10 ⁶ cells in 100μl					
References	 Kishimoto, T.K. <i>et al.</i> (1991) Antibodies against human neutrophil LECAM-1 (LAM-1/Leu-8/DREG-56 antigen) and endothelial cell ELAM-1 inhibit a common CD18-independent adhesion pathway <i>in vitro</i>. Blood. 78: 805-11. Abbassi, O. <i>et al.</i> (1993) E-selectin supports neutrophil rolling <i>in vitro</i> under conditions of flow. J Clin Invest. 92 (6): 2719-30. Nieuwland, R. <i>et al.</i> (2000) Cellular origin and procoagulant properties of microparticles in meningococcal sepsis. Blood. 95: 930-5. Rius, C. <i>et al.</i> (2010) Trans- but not cis-resveratrol impairs angiotensin-II-mediated vascular inflammation through inhibition of NF-κB activation and peroxisome proliferator-activated receptor-gamma upregulation. J Immunol. 185: 3718-27. Marin, C. <i>et al.</i> (2011) Mediterranean diet reduces endothelial damage and improves the regenerative capacity of endothelium. Am J Clin Nutr. 93: 267-74. Wassmer, S.C. <i>et al.</i> (2011) Vascular endothelial cells cultured from patients with cerebral or uncomplicated malaria exhibit differential reactivity to TNF. Cell Microbiol. 13: 198-209. Lo. C. Y. <i>et al.</i> (2013) The use of surface immobilization of P-selectin glycoprotein ligand-1 on mesenchymal stem cells to facilitate selectin mediated cell tethering and rolling Biomaterials. 34 (33): 8213-22. 					
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be					
	protected from light.					
Guarantee	12 months from date of despatch					
Health And Safety Information	Material Safety Datasheet documentation #10058 available at: https://www.bio-rad-antibodies.com/SDS/MCA1969F 10058					

Related Products

Recommended Negative Controls

MOUSE IgG2a NEGATIVE CONTROL:FITC (MCA929F)

Email: antibody_sales_us@bio-rad.com

Recommended Useful Reagents

HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

North & South Tel: +1 800 265 7376 America Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +44 (0)1865 852 739 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M423280:231011'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint