

Datasheet: MCA2111EL

Description:	MOUSE ANTI BOVINE INTERLEUKIN-10:Low Endotoxin	
Specificity:	IL-10	
Format:	Low Endotoxin	
Product Type:	Monoclonal Antibody	
Clone:	CC320	
Isotype:	lgG1	
Quantity:	0.5 mg	

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	•			
Immunohistology - Frozen				
Immunohistology - Paraffin				
ELISA	•			5ug/ml - 10ug/ml
Immunoprecipitation				
Western Blotting				
Functional Assays	-			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) Membrane permeabilization is required for this application. Bio-Rad recommends the use of Leucoperm (BUF09) for this purpose.

Target Species	Bovine
Species Cross Reactivity	Reacts with: Horse, Sheep, Goat N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture

supernatant

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 0.01 EU/µg
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Plasmid cDNA encoding bovine IL-10.
External Database Links	UniProt: P43480 Related reagents
	Entrez Gene: 281246 IL10 Related reagents
RRID	AB_2264769
Fusion Partners	Spleen cells from immunised Balb/c mice were fused with cells of the mouse sp2/0 myeloma cell line.
Specificity	Mouse anti Bovine Interleukin-10 antibody, clone CC320 recognizes bovine IL-10.
	Mouse anti Bovine Interleukin-10 antibody, clone CC320 has been shown to neutralize the activity of bovine IL-10 as measured by the inhibition of the inhibitory activity of IL-10 on IFN gamma synthesis (Buza et al. 2004).
References	 Kwong, L.S. <i>et al.</i> (2002) Development of an ELISA for bovine IL-10. <u>Vet Immunol Immunopathol.</u> 85 (3-4): 213-23. Buza JJ <i>et al.</i> (2004) Neutralization of interleukin-10 significantly enhances gamma interferon expression in peripheral blood by stimulation with Johnin purified protein derivative and by infection with <i>Mycobacterium avium</i> subsp. <i>paratuberculosis</i> in experimentally infected cattle with paratuberculosis. <u>Infect Immun.</u> 72 (4): 2425-8. Bannerman, D.D. <i>et al.</i> (2004) Characterization of the bovine innate immune response to intramammary infection with <i>Klebsiella pneumoniae</i>. <u>J Dairy Sci.</u> 87: 2420-32. Bannerman, D.D. <i>et al.</i> (2004) <i>Escherichia coli</i> and <i>Staphylococcus aureus</i> elicit differential innate immune responses following intramammary infection. <u>Clin Diagn Lab Immunol.</u> 11: 463-72. Abbott, J.R. <i>et al.</i> (2005) Rapid and long-term disappearance of CD4+ T lymphocyte responses specific for <i>Anaplasma marginale</i> major surface protein-2 (MSP2) in MSP2 vaccinates following challenge with live <i>A. marginale</i>. <u>J Immunol.</u> 174: 6702-15. Berger, S.T. and Griffin, F.T. (2006) A comparison of ovine monocyte-derived

- macrophage function following infection with *Mycobacterium avium* ssp. *avium* and *Mycobacterium avium* ssp. *paratuberculosis*. Immunol Cell Biol. 84: 349-56.
- 7. Lei, L. and Hostetter, J.M. (2007) Limited phenotypic and functional maturation of bovine monocyte-derived dendritic cells following *Mycobacterium avium* subspecies *paratuberculosis* infection *in vitro*. <u>Vet Immunol Immunopathol. 120: 177-86.</u>
- 8. Denis, M. *et al.* (2007) Enhancement of the sensitivity of the whole-blood gamma interferon assay for diagnosis of *Mycobacterium bovis* infections in cattle. <u>Clin Vaccine</u> Immunol. 14: 1483-9.
- 9. Davis, T.L. and Pate, J.L. (2007) Bovine luteal cells stimulate proliferation of major histocompatibility nonrestricted gamma delta T cells. Biol Reprod. 77: 914-22.
- 10. Pascale, F. *et al.* (2008) Plasmacytoid dendritic cells migrate in afferent skin lymph. <u>J. Immunol.</u> 180: 5963-72.
- 11. Karcher, E.L. *et al.* (2008) Osteopontin: a novel cytokine involved in the regulation of *Mycobacterium avium* subspecies *paratuberculosis* infection in periparturient dairy cattle. <u>J</u> <u>Dairy Sci. 91: 3079-91.</u>
- 12. Souza, M. *et al.* (2008) Pathogenesis and immune responses in gnotobiotic calves after infection with the genogroup II.4-HS66 strain of human norovirus. <u>J Virol. 82:</u> 1777-86.
- 13. Flynn, R.J. and Mulcahy, G. (2008) Possible role for Toll-like receptors in interaction of Fasciola hepatica excretory/secretory products with bovine macrophages. <u>Infect Immun.</u> 76: 678-84.
- 14. Dean, G.S. *et al.* (2008) Isoniazid treatment of *Mycobacterium bovis* in cattle as a model for human tuberculosis. Tuberculosis (Edinb). 88: 586-94.
- 15. Weiss DJ *et al.* (2008) Bovine monocyte TLR2 receptors differentially regulate the intracellular fate of *Mycobacterium avium* subsp. *paratuberculosis* and *Mycobacterium avium* subsp. *avium*. J Leukoc Biol. 83 (1): 48-55.
- 16. Lybeck, K.R. *et al.* (2009) Neutralization of interleukin-10 from CD14(+) monocytes enhances gamma interferon production in peripheral blood mononuclear cells from *Mycobacterium avium* subsp. *paratuberculosis*-infected goats. Clin Vaccine Immunol. 16 (7): 1003-11.
- 17. Parker, D.G. *et al.* (2009) A steroid-inducible promoter for the cornea. <u>Br J Ophthalmol</u>. 93: 1255-9.
- 18. Hoek, A. *et al.* (2009) Subpopulations of bovine WC1(+) gammadelta T cells rather than CD4(+)CD25(high) Foxp3(+) T cells act as immune regulatory cells *ex vivo*. <u>Vet Res. 40:6.</u>
- 19. Ferret-Bernard, S. *et al.* (2010) Cellular and molecular mechanisms underlying the strong neonatal IL-12 response of lamb mesenteric lymph node cells to R-848. <u>PLoS One.</u> 5: e13705.
- 20. Parker, D.G. *et al.* (2010) Lentivirus-mediated gene transfer of interleukin 10 to the ovine and human cornea. <u>Clin Experiment Ophthalmol. 38: 405-13.</u>
- 21. Rinaldi, M. *et al* (2010) A sentinel function for teat tissues in dairy cows: dominant innate immune response elements define early response to *E. coli* mastitis. <u>Funct Integr</u> Genomics. 10: 21-38.
- 22. Scandurra, G.M. *et al.* (2010) Assessment of live candidate vaccines for paratuberculosis in animal models and macrophages. Infect Immun. 78 (3): 1383-9.
- 23. Yamaguchi, T. *et al.* (2010) Generation of IFN-gamma-producing cells that recognize the major piroplasm surface protein in *Theileria orientalis*-infected bovines. <u>Vet Parasitol.</u>

171: 207-15.

- 24. Coad, M. *et al.* (2010) Repeat tuberculin skin testing leads to desensitisation in naturally infected tuberculous cattle which is associated with elevated interleukin-10 and decreased interleukin-1 beta responses. Vet Res. 41: 14.
- 25. Wenz, J.R. *et al.* (2010) Factors associated with concentrations of select cytokine and acute phase proteins in dairy cows with naturally occurring clinical mastitis. <u>J Dairy Sci.</u> 93: 2458-70.
- 26. Ferret-Bernard, S. *et al.* (2011) Mesenteric lymph node cells from neonates present a prominent IL-12 response to CpG oligodeoxynucleotide via an IL-15 feedback loop of amplification. <u>Vet Res. 42:19.</u>
- 27. Stabel, J.R. and Robbe-Austerman, S. (2011) Early immune markers associated with *Mycobacterium avium* subsp. *paratuberculosis* infection in a neonatal calf model. <u>Clin</u> Vaccine Immunol. 18: 393-405.
- 28. Ma&347;lanka T *et al.* (2012) The presence of CD25 on bovine WC1+ gammadelta T cells is positively correlated with their production of IL-10 and TGF-beta, but not IFN-gamma. Pol J Vet Sci. 15 (1): 11-20.
- 29. Lybeck, K.R. *et al.* (2012) Intestinal Strictures, Fibrous Adhesions and High Local Interleukin-10 Levels in Goats Infected Naturally with *Mycobacterium avium* subsp. *paratuberculosis*. J Comp Pathol. 148: 157-72.
- 30. Dooley, L.M. *et al.* (2015) Effect of mesenchymal precursor cells on the systemic inflammatory response and endothelial dysfunction in an ovine model of collagen-induced arthritis. <u>PLoS One. 10 (5): e0124144.</u>
- 31. Pomeroy, B. *et al.* (2015) Monocyte-derived dendritic cells from late gestation cows have an impaired ability to mature in response to *E. coli* stimulation in a receptor and cytokine-mediated fashion. <u>Vet Immunol Immunopathol. 167 (1-2): 22-9.</u>
- 32. Pomeroy, B. *et al.* (2016) Impact of *in vitro* treatments of physiological levels of estradiol and progesterone observed in pregnancy on bovine monocyte-derived dendritic cell differentiation and maturation. Vet Immunol Immunopathol. 182: 37-42.
- 33. Rainard, P. *et al.* (2016) Innate and Adaptive Immunity Synergize to Trigger Inflammation in the Mammary Gland. <u>PLoS One. 11 (4): e0154172.</u>
- 34. Rodrigues, V. *et al.* (2017) Development of a bead-based multiplexed assay for simultaneous quantification of five bovine cytokines by flow cytometry. <u>Cytometry A. 91</u> (9): 901-7.
- 35. Westover, A. *et al.* (2017) Effect of Human Amnion Epithelial Cells on the Acute Inflammatory Response in Fetal Sheep. <u>Front Physiol. 8: 871.</u>
- 36. Cassady-Cain, R.L. *et al.* (2017) Inhibition of Antigen-Specific and Nonspecific Stimulation of Bovine T and B Cells by Lymphostatin from Attaching and Effacing *Escherichia coli*. Infect Immun. 85 (2): e00845-16.
- 37. Jimbo, S. *et al.* (2019) Natural and inducible regulatory B cells are widely distributed in ovine lymphoid tissues. <u>Vet Immunol Immunopathol. 211: 44-8.</u>
- 38. Stabel, J.R. & Bannantine, J.P. (2019) Divergent Antigen-Specific Cellular Immune Responses during Asymptomatic Subclinical and Clinical States of Disease in Cows Naturally Infected with *Mycobacterium avium*. subsp. *paratuberculosis*. <u>Infect Immun.</u> 88(1):e00650-19.
- 39. Galinsky, R. *et al.* (2020) Tumor necrosis factor inhibition attenuates white matter gliosis after systemic inflammation in preterm fetal sheep. <u>J Neuroinflammation</u>. 17 (1): 92.
- 40. Borghesi, J. et al. (2020) Effects of doxorubicin associated with amniotic membrane

stem cells in the treatment of canine inflammatory breast carcinoma (IPC-366) cells. <u>BMC</u> Vet Res. 16 (1): 353.

- 41. Ciliberti, M.G. *et al.* (2020) Nexus Between Immune Responses and Oxidative Stress: The Role of Dietary Hydrolyzed Lignin in *ex vivo* Bovine Peripheral Blood Mononuclear Cell Response. Front Vet Sci. 7: 9.
- 42. Davidson, J.O. *et al.* (2021) Window of opportunity for human amnion epithelial stem cells to attenuate astrogliosis after umbilical cord occlusion in preterm fetal sheep. <u>Stem Cells Transl Med. 10 (3): 427-40.</u>
- 43. Stabel, J.R. *et al.* (2021) Comparative cellular immune responses in calves after infection with *Mycobacterium avium.* subsp. *paratuberculosis.*, *M. avium.* subsp. *avium.*, *M. kansasii.* and *M. bovis.*. Vet Immunol Immunopathol. 237: 110268.
- 44. Ciliberti, M.G. *et al.* (2022) Green extraction of bioactive compounds from wine lees and their bio-responses on immune modulation using *in vitro* sheep model. <u>J Dairy Sci.</u> 105 (5): 4335-53.
- 45. Santillo, A. *et al.* (2022) Feeding tannins to dairy cows in different seasons improves the oxidative status of blood plasma and the antioxidant capacity of cheese. <u>J Dairy Sci.</u> 105 (11): 8609-20.
- 46. Bouroutzika, E. *et al.* (2023) Melatonin Administration to Pregnant Ewes for Coccidiosis Control in Their Offspring. <u>Animals (Basel)</u>. 13 (14): 2381.
- 47. de Silva, K.. *et al.* (2018) Defining resilience to mycobacterial disease: Characteristics of survivors of ovine paratuberculosis. <u>Vet Immunol Immunopathol. 195: 56-64.</u>
- 48. Stabel, J.R. *et al.* (2020) Comparison of Sheep, Goats, and Calves as Infection Models for *Mycobacterium avium* subsp. *paratuberculosis*. <u>Vet Immunol Immunopathol</u>. 225: 110060.
- 49. Ciliberti, M.G. *et al.* (2024) Role of hazelnut skin supplementation on plasma antioxidant status and cytokine profile in growing lambs. <u>Front Vet Sci. 11: 1340141.</u> 50. Hong, S. *et al.* (2024) Impact of an Injectable Trace Mineral Supplement on the Immune Response and Outcome of *Mannheimia haemolytica* Infection in Feedlot Cattle. Biol Trace Elem Res. Jun 10 [Epub ahead of print].

Further Reading

1. Hamza, E. *et al.* (2007) Modulation of allergy incidence in icelandic horses is associated with a change in IL-4-producing T cells. Int Arch Allergy Immunol. 144: 325-37.

Storage

Store at -20°C only.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Health And Safety Information

Material Safety Datasheet documentation #10162 available at:

https://www.bio-rad-antibodies.com/SDS/MCA2111EL

10162

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR12...) RPE

Goat Anti Mouse IgG IgA IgM (STAR87...) HRP

Goat Anti Mouse IgG (STAR76...) RPE

Goat Anti Mouse IgG (STAR70...) FITC

Rabbit Anti Mouse IgG (STAR13...) HRP

Goat Anti Mouse IgG (Fc) (STAR120...) FITC, HRP

Rabbit Anti Mouse IgG (STAR9...) <u>FITC</u>

Goat Anti Mouse IgG (STAR77...) HRP

Goat Anti Mouse IgG (H/L) (STAR117...) Alk. Phos., DyLight®488, DyLight®550,

DyLight®650, DyLight®680, DyLight®800,

FITC, HRP

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL:Low Endotoxin (MCA928EL)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M437823:250319'

Printed on 19 Mar 2025

© 2025 Bio-Rad Laboratories Inc | Legal | Imprint