

Datasheet: MCA2257A647

| Description:  | MOUSE ANTI HUMAN CD226:Alexa Fluor® 647 |
|---------------|---|
| Specificity:  | CD226                                   |
| Other names:  | DNAM-1                                  |
| Format:       | ALEXA FLUOR® 647                        |
| Product Type: | Monoclonal Antibody                     |
| Clone:        | DX11                                    |
| Isotype:      | IgG1                                    |
| Quantity:     | 100 TESTS/1ml                           |
|               |   |

# **Product Details**

# **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

|                | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry |     |    |                | Neat               |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| Target Species                    | Human  |                       |                                |          |
|-----------------------------------|--|-----------------------|--------------------------------|----------|
| Product Form                      | Purified IgG conjugate                                   | ed to Alexa Fluor® 64 | 7 - liquid                     |          |
| Max Ex/Em                         | Fluorophore  | Excitation Max (nm)   | Emission Max (nm)              |          |
|                                   | Alexa Fluor®647  | 650                   | 665                            |          |
| Preparation  Buffer Solution      | Purified IgG prepared supernatant  Phosphate buffered sa |                       | raphy on Protein A from tissue | e cultur |
| Preservative<br>Stabilisers       | 0.09% sodium azide (<br>1% bovine serum albu             | . 0,                  |                                |          |
| Approx. Protein<br>Concentrations | IgG concentration 0.0                                    | 5 mg/ml               |                                |          |

| ı | m   | m   | u | n  | o | a | e | n |
|---|-----|-----|---|----|---|---|---|---|
| • | ••• | ••• | • | •• | • | 3 | • |   |

Human cytotoxic T lymphocyte clone

# External Database Links

#### **UniProt:**

Q15762 Related reagents

#### **Entrez Gene:**

10666 CD226 Related reagents

#### **Synonyms**

DNAM1

#### **Fusion Partners**

Spleen cells from immunized BALB/c mouse were fused with cells of the Sp2/0 myeloma cell line

### **Specificity**

Mouse anti Human CD226 antibody, clone DX11 recognizes human CD226, a ~65 kDa glycoprotein, also known as DNAM1 (DNAX accessory molecule-1). CD226 is broadly expressed on T-cells, NK cells, platelets, monocytes and a subset of B cells. CD226 is also expressed by a subset of CD3 positive thymocytes.

Mouse anti Human CD226 antibody, clone DX11 is reported to inhibit T- and NK cell mediated cytotoxicity against tumor cell targets and to block TNF alpha and IFN gamma secretion by alloantigen-specific T-cells (Kojima *et al.* 2003).

# Flow Cytometry

Use 10µl of the suggested working dilution to label 1 x 10<sup>6</sup> cells in 100µl

# References

- 1. Shibuya, A. *et al.* (1996) DNAM-1, a novel adhesion molecule involved in the cytolytic function of T lymphocytes. Immunity. 4 (6): 573-81.
- 2. Shibuya, A. *et al.* (1998) Protein kinase C is involved in the regulation of both signaling and adhesion mediated by DNAX accessory molecule-1 receptor. <u>J Immunol. 161 (4):</u> 1671-6.
- 3. Kojima, H. *et al.* (2003) CD226 mediates platelet and megakaryocytic cell adhesion to vascular endothelial cells. <u>J Biol Chem. 278 (38): 36748-53.</u>
- 4. Manes, T.D. and Pober, J.S. (2011) Identification of Endothelial Cell Junctional Proteins and Lymphocyte Receptors Involved in Transendothelial Migration of Human Effector Memory CD4+ T Cells. <u>J Immunol. 186: 1763-8.</u>
- 5. Ardolino, M. *et al.* (2011) DNAM-1 ligand expression on Ag-stimulated T lymphocytes is mediated by ROS-dependent activation of DNA-damage response: relevance for NK-T cell interaction. <u>Blood</u>. 117: 4778-86.
- 6. Soriani, A. *et al.* (2009) ATM-ATR-dependent up-regulation of DNAM-1 and NKG2D ligands on multiple myeloma cells by therapeutic agents results in enhanced NK-cell susceptibility and is associated with a senescent phenotype. <u>Blood. 113: 3503-11.</u>
- 7. Fionda, C. *et al.* (2009) Heat shock protein-90 inhibitors increase MHC class I-related chain A and B ligand expression on multiple myeloma cells and their ability to trigger NK cell degranulation. <u>J Immunol.</u> 183 (7): 4385-94.
- 8. Matusali, G. *et al.* (2012) The Human Immunodeficiency Virus Type 1 Nef and Vpu Proteins Downregulate the Natural Killer Cell-Activating Ligand PVR. <u>J Virol. 86:</u> 4496-504.
- 9. Fionda, C. et al. (2015) Nitric oxide donors increase PVR/CD155 DNAM-1 ligand

expression in multiple myeloma cells: role of DNA damage response activation. <u>BMC</u> Cancer. 15 (1): 17.

- 10. Molfetta, R. *et al.* (2019) The Ubiquitin-proteasome pathway regulates Nectin2/CD112 expression and impairs NK cell recognition and killing. <u>Eur J Immunol. 49 (6): 873-83.</u>
- 11. Vulpis, E. *et al.* (2022) Impact on NK cell functions of acute versus chronic exposure to extracellular vesicle-associated MICA: Dual role in cancer immunosurveillance. <u>J Extracell Vesicles</u>. 11 (1): e12176.
- 12. Molfetta, R. *et al.* (2020) CD155: A Multi-Functional Molecule in Tumor Progression. Int J Mol Sci. 21 (3): 922.
- 13. Zitti, B. *et al.* (2017) Innate immune activating ligand SUMOylation affects tumor cell recognition by NK cells. Sci Rep. 7 (1): 10445.
- 14. Mekhloufi, A. *et al.* (2020) Bone Marrow Stromal Cell-Derived IL-8 Upregulates PVR Expression on Multiple Myeloma Cells via NF-kB Transcription Factor. <u>Cancers (Basel).</u> 12 (2): 440.

#### **Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

#### Guarantee

12 months from date of despatch

# Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

# Health And Safety Information

Material Safety Datasheet documentation #10041 available at: <a href="https://www.bio-rad-antibodies.com/SDS/MCA2257A647">https://www.bio-rad-antibodies.com/SDS/MCA2257A647</a> 10041

# Regulatory

For research purposes only

# Related Products

# **Recommended Negative Controls**

MOUSE IgG1 NEGATIVE CONTROL: Alexa Fluor® 647 (MCA928A647)

# **Recommended Useful Reagents**

# HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M413618:221123'

# Printed on 19 Oct 2023

© 2023 Bio-Rad Laboratories Inc | Legal | Imprint