

Datasheet: MCA2334A488

BATCH NUMBER 162449

Description:	MOUSE ANTI BOVINE TNF ALPHA:Alexa Fluor® 488			
Specificity:	TNF ALPHA			
Format:	ALEXA FLUOR® 488			
Product Type:	Monoclonal Antibody			
Clone:	CC327			
Isotype:	lgG2b			
Quantity:	100 TESTS/1ml			

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	-			1/10 - 1/50

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)Membrane permeabilisation is required for this application. Bio-Rad recommends the use of Leucoperm $^{\text{TM}}$ (Product Code <u>BUF09</u>) for this purpose.

Target Species	Bovine			
Species Cross	Reacts with: Fallow	v deer		
Reactivity	reactivity is derived	tivity and working conditi I from testing within our lications from the originate	aboratories, peer-reviev	wed publica
Product Form	Purified IgG conjug	gated to Alexa Fluor® 48	8 - liquid	
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)	
	Alexa Fluor®488	495	519	
Preparation	Purified IgG prepar	red by affinity chromatog	raphy on Protein A from	າ tissue cult

Buffer Solution	Phosphate buffered saline	
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum Albumin	
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml	
Immunogen	Recombinant bovine TNF alpha.	
External Database Links	UniProt: Q06599 Related reagents Entrez Gene: 280943 TNF Related reagents	
Synonyms	TNFA, TNFSF2	
RRID	AB_567242	
Specificity	Mouse anti Bovine TNF alpha antibody, clone CC327 recogni 17.5 kDa cytokine, expressed by many different stimulated cell ty macrophages, endothelial cells, fibroblasts and both T and B-lym. The production of TNF alpha is induced by a variety of factors, d and includes bacterial toxins, IL-1, PDGF, IFN-beta, NGF, Oncos infections. The presence of TNF alpha is responsible for diverse anti-tumour and toxic effects and under certain conditions is also self-stimulation and inhibition.	ypes including monocytes, aphocytes. Ilependant upon cell type statin M and viral immunomodulatory,
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in	n 100ul.
References	 Hope, J.C. <i>et al.</i> (2003) Maturation of bovine dendritic cells by Immunol Immunopathol. 95 (1-2): 21-31. Whelan, A.O. <i>et al.</i> (2003) Modulation of the bovine delayed-ty responses to defined mycobacterial antigens by a synthetic bacter Immun. 71 (11): 6420-5. Guergnon J <i>et al.</i> (2003) A tumour necrosis factor alpha autocomproliferation and nuclear factor-kappaβ activation of <i>Theileria para Cell Microbiol.</i> 5 (10): 709-16. Kwong, L.S. <i>et al.</i> (2010) Production and characterization of two bovine tumournecrosis factor alpha (TNF-alpha) and their cross TNF-alpha. Vet Immunol Immunopathol. 135: 320-4. Wenz, J.R. <i>et al.</i> (2010) Factors associated with concentration acute phase proteins in dairy cows with naturally occurring clinical 93: 2458-70. 	ype hypersensitivity erial lipopeptide. Infect crine loop contributes to rva-transformed B cells. wo monoclonal antibodies es-reactivity with ovine as of select cytokine and al mastitis. J Dairy Sci.

6. Rinaldi, M. *et al* (2010) A sentinel function for teat tissues in dairy cows: dominant

innate immune response elements define early response to *E. coli* mastitis. <u>Funct Integr</u> Genomics. 10: 21-38.

- 7. Sow, F.B. *et al.* (2011) Respiratory syncytial virus is associated with an inflammatory response in lungs and architectural remodeling of lung-draining lymph nodes of newborn lambs. Am J Physiol Lung Cell Mol Physiol. 300 (1): L12-24.
- 8. Simojoki, H. *et al.* (2011) Innate immune response in experimentally induced bovine intramammary infection with *Staphylococcus simulans* and *S. epidermidis*. <u>Vet Res. 42:</u> 49.
- 9. Whelan, A.O. *et al.* (2011) Development of an Antibody to Bovine IL-2 Reveals Multifunctional CD4 T(EM) Cells in Cattle Naturally Infected with Bovine Tuberculosis. PLoS One. 6: e29194.
- 10. García-Jiménez, W.L. (2012) Histological and immunohistochemical characterisation of Mycobacterium bovis induced granulomas in naturally infected fallow deer (*Dama dama*). <u>Vet Immunol Immunopathol. 149: 66-75.</u>
- 11. Redondo, E. *et al.* (2014) Induction of interleukin-8 and interleukin-12 in neonatal ovine lung following experimental inoculation of bovine respiratory syncytial virus. <u>J Comp Pathol. 150 (4): 434-48.</u>
- 12. Luisa, C. *et al.* (2016) Evaluation of serum markers of blood redox homeostasis and inflammation in PCB naturally contaminated heifers undergoing decontamination <u>Science</u> of The Total Environment. 542: 653-64.
- 13. Maggioli, M.F. *et al.* (2016) Increased TNF-α/IFN-γ/IL-2 and Decreased TNF-α/IFN-γ Production by Central Memory T Cells Are Associated with Protective Responses against Bovine Tuberculosis Following BCG Vaccination. Front Immunol. 7: 421.
- 14. Rutigliano, H.M. *et al.* (2016) Trophoblast Major Histocompatibility Complex Class I Expression Is Associated with Immune-Mediated Rejection of Bovine Fetuses Produced by Cloning. <u>Biol Reprod. 95 (2): 39.</u>
- 15. Camejo, M.I. *et al.* (2014) TNF-alpha in bulls experimentally infected with *Trypanosoma vivax*: a pilot study. <u>Vet Immunol Immunopathol.</u> 162 (3-4): 192-7.
- 16. Jolly A *et al.* (2016) Evidence of a pro-apoptotic effect of specific antibodies in a bovine macrophage model of infection with *Mycobacterium avium* subsp. paratuberculosis. Vet Immunol Immunopathol. 169: 47-53.
- 17. Rodrigues, V. *et al.* (2017) Development of a bead-based multiplexed assay for simultaneous quantification of five bovine cytokines by flow cytometry. <u>Cytometry A. 91</u> (9): 901-7.
- 18. Sirak, A. *et al.* (2021) Cellular and Cytokine Responses in Lymph Node Granulomas of Bacillus Calmette Guérin (BCG)-Vaccinated and Non-vaccinated Cross-Breed Calves Naturally Infected With *Mycobacterium bovis*. <u>Front Vet Sci. 8: 698800</u>.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee

12 months from date of despatch

Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information

Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA2334A488

10041

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG2b NEGATIVE CONTROL: Alexa Fluor® 488 (MCA691A488)

North & South Tel: +1 800 265 7376 America

Fax: +1 919 878 3751

Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_uk@bio-rad.com Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M385382:210513'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint