

Datasheet: MCA2334A488

BATCH NUMBER 166701

| Description: | MOUSE ANTI BOVINE TNF ALPHA:Alexa Fluor® 488 | | | |
|---------------|--|--|--|--|
| Specificity: | TNF ALPHA | | | |
| Format: | ALEXA FLUOR® 488 | | | |
| Product Type: | Monoclonal Antibody | | | |
| Clone: | CC327 | | | |
| Isotype: | lgG2b | | | |
| Quantity: | 100 TESTS/1ml | | | |
| | | | | |

Product Details

| Applications | This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-</u> | | | | | | |
|-----------------------------|--|------------------------|-------------------------|----------------------|--|--|--|
| | rad-antibodies.com/protocols. | | | | | | |
| | | Yes No | Not Determined | Suggested Dilution | | | |
| | Flow Cytometry (1) | • | | 1/10 - 1/50 | | | |
| | Where this antibody has not been tested for use in a particular technique this does not | | | | | | |
| | necessarily exclude its use in such procedures. It is recommended that the user titrates | | | | | | |
| | the antibody for use ir | n their own system usi | ng appropriate negative | e/positive controls. | | | |
| | (1) Membrane permeabilization is required for this application. The use of | | | | | | |
| | Leucoperm (Product Code <u>BUF09</u>) is recommended for this purpose. | | | | | | |
| Target Species | Bovine | | | | | | |
| Species Cross Reactivity | Reacts with: Fallow deer N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. | | | | | | |
| Product Form | Purified IgG conjugated to Alexa Fluor® 488 - liquid | | | | | | |
| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) | | | | |
| | Alexa Fluor®488 | 495 | 519 | | | | |
| Preparation | Purified IgG prepared supernatant | by affinity chromatog | raphy on Protein A fron | n tissue culture | | | |

| Buffer Solution | Phosphate buffered saline |
|-----------------------------------|--|
| Preservative Stabilisers | 0.09% Sodium Azide 1% Bovine Serum Albumin |
| Approx. Protein Concentrations | IgG concentration 0.05 mg/ml |
| Immunogen | Recombinant bovine TNF alpha. |
| External Database Links | UniProt: Q06599 Related reagents Entrez Gene: 280943 TNF Related reagents |
| Synonyms | TNFA, TNFSF2 |
| RRID | AB_567242 |
| Specificity | Mouse anti Bovine TNF alpha antibody, clone CC327 recognizes bovine TNF alpha, a 17.5 kDa cytokine, expressed by many different stimulated cell types including monocytes, macrophages, endothelial cells, fibroblasts and both T and B-lymphocytes. The production of TNF alpha is induced by a variety of factors, dependant upon cell type and includes bacterial toxins, IL-1, PDGF, IFN-beta, NGF, Oncostatin M and viral infections. The presence of TNF alpha is responsible for diverse immunomodulatory, anti-tumour and toxic effects and under certain conditions is also capable of self-stimulation and inhibition. |
| Flow Cytometry | Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul. |
| References | Hope, J.C. <i>et al.</i> (2003) Maturation of bovine dendritic cells by lipopeptides. <u>Vet</u> <u>Immunol Immunopathol. 95 (1-2): 21-31.</u> Whelan, A.O. <i>et al.</i> (2003) Modulation of the bovine delayed-type hypersensitivity responses to defined mycobacterial antigens by a synthetic bacterial lipopeptide. <u>Infect</u> <u>Immun. 71 (11): 6420-5.</u> Guergnon J <i>et al.</i> (2003) A tumour necrosis factor alpha autocrine loop contributes to proliferation and nuclear factor-kappaβ activation of <i>Theileria parva</i>-transformed B cells. <u>Cell Microbiol. 5 (10): 709-16.</u> Kwong, L.S. <i>et al.</i> (2010) Production and characterization of two monoclonal antibodies to bovine tumournecrosis factor alpha (TNF-alpha) and their cross-reactivity with ovine TNF-alpha. <u>Vet Immunol Immunopathol. 135: 320-4.</u> Wenz, J.R. <i>et al.</i> (2010) Factors associated with concentrations of select cytokine and acute phase proteins in dairy cows with naturally occurring clinical mastitis. J Dairy Sci. |

innate immune response elements define early response to *E. coli* mastitis. <u>Funct Integr</u> <u>Genomics. 10: 21-38.</u>

7. Sow, F.B. *et al.* (2011) Respiratory syncytial virus is associated with an inflammatory response in lungs and architectural remodeling of lung-draining lymph nodes of newborn lambs. <u>Am J Physiol Lung Cell Mol Physiol. 300 (1): L12-24.</u>

8. Simojoki, H. *et al.* (2011) Innate immune response in experimentally induced bovine intramammary infection with *Staphylococcus simulans* and *S. epidermidis*. <u>Vet Res. 42:</u> 49.

9. Whelan, A.O. *et al.* (2011) Development of an Antibody to Bovine IL-2 Reveals Multifunctional CD4 T(EM) Cells in Cattle Naturally Infected with Bovine Tuberculosis. <u>PLoS One. 6: e29194.</u>

10. García-Jiménez, W.L. (2012) Histological and immunohistochemical characterisation of Mycobacterium bovis induced granulomas in naturally infected fallow deer (*Dama dama*). Vet Immunol Immunopathol. 149: 66-75.

11. Redondo, E. *et al.* (2014) Induction of interleukin-8 and interleukin-12 in neonatal ovine lung following experimental inoculation of bovine respiratory syncytial virus. <u>J Comp</u> Pathol. 150 (4): 434-48.

12. Luisa, C. *et al.* (2016) Evaluation of serum markers of blood redox homeostasis and inflammation in PCB naturally contaminated heifers undergoing decontamination <u>Science</u> of The Total Environment. 542: 653-64.

13. Maggioli, M.F. *et al.* (2016) Increased TNF- α /IFN- γ /IL-2 and Decreased TNF- α /IFN- γ Production by Central Memory T Cells Are Associated with Protective Responses against Bovine Tuberculosis Following BCG Vaccination. <u>Front Immunol. 7: 421.</u>

14. Rutigliano, H.M. *et al.* (2016) Trophoblast Major Histocompatibility Complex Class I Expression Is Associated with Immune-Mediated Rejection of Bovine Fetuses Produced by Cloning. <u>Biol Reprod. 95 (2): 39.</u>

15. Camejo, M.I. *et al.* (2014) TNF-alpha in bulls experimentally infected with *Trypanosoma vivax*: a pilot study. <u>Vet Immunol Immunopathol. 162 (3-4): 192-7.</u>

16. Jolly A *et al.* (2016) Evidence of a pro-apoptotic effect of specific antibodies in a bovine macrophage model of infection with *Mycobacterium avium* subsp. paratuberculosis. Vet Immunol Immunopathol. 169: 47-53.

17. Rodrigues, V. *et al.* (2017) Development of a bead-based multiplexed assay for simultaneous quantification of five bovine cytokines by flow cytometry. <u>Cytometry A. 91</u> (<u>9): 901-7.</u>

18. Sirak, A. *et al.* (2021) Cellular and Cytokine Responses in Lymph Node Granulomas of Bacillus Calmette Guérin (BCG)-Vaccinated and Non-vaccinated Cross-Breed Calves Naturally Infected With *Mycobacterium bovis*. <u>Front Vet Sci. 8: 698800</u>.

StorageThis product is shipped at ambient temperature. It is recommended to aliquot and store at
-20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for
short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee 12 months from date of despatch

| Acknowledgements | This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com |
|----------------------------------|--|
| Health And Safety Information | Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA2334A488 10041 |

Regulatory For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG2b NEGATIVE CONTROL:Alexa Fluor® 488 (MCA691A488)

| North & South | Tel: +1 800 265 7376 | Worldwide | Tel: +44 (0)1865 852 700 | Europe | Tel: +49 (0) 89 8090 95 21 |
|---------------|---------------------------------|-----------|---------------------------------|--------|--------------------------------------|
| America | Fax: +1 919 878 3751 | | Fax: +44 (0)1865 852 739 | | Fax: +49 (0) 89 8090 95 50 |
| | Email: antibody_sales_us@bio-ra | id.com | Email: antibody_sales_uk@bio-ra | ad.com | Email: antibody_sales_de@bio-rad.com |

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M404136:220820'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint