

Datasheet: MCA2411A647

BATCH NUMBER 167451

Description:	MOUSE ANTI DOG CD34:Alexa Fluor® 647
Specificity:	CD34
Format:	ALEXA FLUOR® 647
Product Type:	Monoclonal Antibody
Clone:	1H6
Isotype:	lgG1
Quantity:	100 TESTS/1ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry				Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Purified IgG conjugated to Alexa Fluor®647- liquid Ex/Em Fluorophore Alexa Fluor®647 Fluorophore Bexcitation Max (nm) Alexa Fluor®647 Fluorophore Bexcitation Max (nm) Fluorophore Bexcitation Max				
Fluorophore Excitation Max (nm) Emission Max (nm) Alexa Fluor®647 650 665 Furation Purified IgG prepared by affinity chromatography on Protein A from supernatant For Solution Phosphate buffered saline Fraction Phosphate buffered saline Fraction Phosphate buffered saline Fraction Phosphate buffered saline For Solution Phosphate buffe	get Species	Dog		
Alexa Fluor®647 650 665 Purified IgG prepared by affinity chromatography on Protein A frosupernatant Phosphate buffered saline Protein 0.09% Sodium Azide isers 1% Bovine Serum Albumin	oduct Form	Purified IgG conjugate	ed to Alexa Fluor®647	'- liquid
Purified IgG prepared by affinity chromatography on Protein A frosupernatant Phosphate buffered saline Protein O.09% Sodium Azide Isers IgG concentration 0.05 mg/ml	x Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm
supernatant r Solution Phosphate buffered saline rvative 0.09% Sodium Azide isers 1% Bovine Serum Albumin ox. Protein InG concentration 0.05 mg/ml		Alexa Fluor®647	650	665
rvative 0.09% Sodium Azide isers 1% Bovine Serum Albumin ox. Protein	eparation	•	by affinity chromatog	raphy on Protein A
isers 1% Bovine Serum Albumin ox. Protein IgG concentration 0.05 mg/ml	ffer Solution	Phosphate buffered sa	aline	
ox. Protein InG concentration 0.05 mg/ml	eservative	0.09% Sodium Azide		
IdG concentration () ()5 mg/ml	ıbilisers	1% Bovine Serum	Albumin	
	prox. Protein ncentrations	IgG concentration 0.08	5 mg/ml	

Immuno	gen
--------	-----

Canine CD34 fusion protein.

External Database

Links

UniProt:

Q28270 Related reagents

Entrez Gene:

415130 CD34 Related reagents

RRID

AB 2074486

Fusion Partners

Spleen cells from immunized BALB/c mice were fused with cells of the mouse NS-1/FOX-NY myeloma cell line.

Specificity

Mouse anti dog CD34 antibody, clone 1H6 recognizes the canine homologue of CD34, a glycosylated type 1 transmembrane protein of approximately 110 kDa (<u>McSweeney et al. 1998</u>) expressed on the cell suface of endothelial cells and haematopoietic stem cells.

Mouse anti dog CD34 antibody, clone 1H6 is a key marker of canine hematopoietic progenitor cells and is reported for use in CD34+ enrichment assays, (Goerner et al. 2001) and (Horn et al. 2004).

Flow Cytometry

Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul.

References

- 1. Goerner, M. *et al.* (1999) The use of granulocyte colony-stimulating factor during retroviral transduction on fibronectin fragment CH-296 enhances gene transfer into hematopoietic repopulating cells in dogs. <u>Blood. 94 (7): 2287-92.</u>
- 2. Bhattacharya, V. *et al.* (2000) Enhanced endothelialization and microvessel formation in polyester grafts seeded with CD34(+) bone marrow cells. Blood. 95 (2): 581-5.
- 3. Goerner, M. *et al.* (2001) Sustained multilineage gene persistence and expression in dogs transplanted with CD34(+) marrow cells transduced by RD114-pseudotype oncoretrovirus vectors. Blood. 98 (7): 2065-70.
- 4. Georges, G. *et al.* (2001) Engraftment of DLA-haploidentical marrow with ex vivo expanded, retrovirally transduced cytotoxic T lymphocytes. <u>Blood. 98:3447-55.</u>
- 5. Horn, P.A. *et al.* (2004) Efficient lentiviral gene transfer to canine repopulating cells using an overnight transduction protocol. <u>Blood</u>. 103 (10): 3710-6.
- 6. Avallone, G. *et al.* (2007) The spectrum of canine cutaneous perivascular wall tumors: morphologic, phenotypic and clinical characterization. <u>Vet Pathol. 44 (5): 607-20.</u>
- 7. Palmieri, C. *et al.* (2013) Use of electron microscopy to classify canine perivascular wall tumors. Vet Pathol. 50 (2): 226-33.
- 8. Bearden, R.N. *et al.* (2017) *In-vitro* characterization of canine multipotent stromal cells isolated from synovium, bone marrow, and adipose tissue: a donor-matched comparative study. <u>Stem Cell Res Ther. 8 (1): 218.</u>
- 9. Trindade, A.B. *et al.* (2017) Mesenchymal-like stem cells in canine ovary show high differentiation potential. Cell Prolif. Oct 08 [Epub ahead of print].
- 10. Lee, S.H. *et al.* (2016) Impact of local injection of brain-derived neurotrophic factor-expressing mesenchymal stromal cells (MSCs) combined with intravenous MSC delivery in a canine model of chronic spinal cord injury. Cytotherapy. Oct 28 [Epub ahead of print].

- 11. Muir, P. *et al.* (2016) Autologous Bone Marrow-Derived Mesenchymal Stem Cells Modulate Molecular Markers of Inflammation in Dogs with Cruciate Ligament Rupture. PLoS One. 11 (8): e0159095.
- 12. Rajawat, Y.S. *et al.* (2021) *In Vivo* Gene Therapy for Canine SCID-X1 Using Cocal-Pseudotyped Lentiviral Vector. <u>Hum Gene Ther. 32 (1-2): 113-27.</u>
- 13. Grudzien, M. *et al.* (2021) A newly established canine NK-type cell line and its cytotoxic properties. <u>Vet Comp Oncol. 19 (3): 567-77.</u>
- 14. Tongu, E.A.O. *et al.* (2021) Allogenic mesenchymal stem cell-conditioned medium does not affect sperm parameters and mitigates early endometrial inflammatory responses in mares. Theriogenology. 169: 1-8.
- 15. Jaensch, S. *et al.* (2022) Clinicopathologic and immunophenotypic features in dogs with presumptive large granular lymphocyte leukaemia <u>Australian Veterinary Journal.</u> [Epub ahead of print].
- 16. Salari Sedigh, H. *et al.* (2023) *In vitro* investigation of canine periodontal ligament-derived mesenchymal stem cells: A possibility of promising tool for periodontal regeneration. J Oral Biol Craniofac Res. 13 (3): 403-11.
- 17. Papa, P.M. *et al.* (2023) Intratesticular transplantation of allogenic mesenchymal stem cells mitigates testicular destruction after induced heat stress in Miniature-horse stallions. J Equine Vet Sci. 132: 104961.
- 18. Rezaei, M. *et al.* (2019) Transplantation of Bone Marrow-Derived Mesenchymal Stem Cells, Platelet-Rich Plasma, and Fibrin Glue for Periodontal Regeneration. <u>Int J Periodontics Restorative Dent. 39 (1): e32-e45.</u>
- 19. Yang, V.K. *et al.* (2021) Intravenous administration of allogeneic Wharton jelly-derived mesenchymal stem cells for treatment of dogs with congestive heart failure secondary to myxomatous mitral valve disease. Am J Vet Res. 82 (6): 487-93.
- 20. Crain, S.K. *et al.* (2019) Extracellular Vesicles from Wharton's Jelly Mesenchymal Stem Cells Suppress CD4 Expressing T Cells Through Transforming Growth Factor Beta and Adenosine Signaling in a Canine Model. <u>Stem Cells Dev. 28 (3): 212-26.</u>
- 21. Sheng, R. *et al.* (2023) Prognostic significance of CD25 expression in dogs with a noninvasive diagnosis of B-cell lymphoma treated with CHOP chemotherapy. <u>Vet Comp</u> Oncol. 21 (1): 28-35.
- 22. Millanta, F. *et al.* (2020) Cytologic grading of canine and feline spindle-cell sarcomas of soft tissues and its correlation with histologic grading. <u>Top Companion Anim Med. 41:</u> 100458.

Further Reading

1. McSweeney, P. *et al.* (1996) Canine CD34: cloning of the cDNA and evaluation of an antiserum to recombinant protein. <u>Blood</u>. 88:1992-2003.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee

12 months from date of despatch

Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information

Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA2411A647

10041

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL: Alexa Fluor® 647 (MCA928A647)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M385407:210513'

Printed on 25 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint

Email: antibody_sales_us@bio-rad.com