

Datasheet: MCA2419A488

BATCH NUMBER 157014

Description:	MOUSE ANTI HUMAN CD62P:Alexa Fluor® 488		
Specificity:	CD62P		
Other names:	P-SELECTIN		
Format:	ALEXA FLUOR® 488		
Product Type:	Monoclonal Antibody		
Clone:	Psel.KO.2.7		
Isotype:	lgG1		
Quantity:	100 TESTS/1ml		

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	-			Neat - 1/5

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Species Cross	Reacts
Reactivity	N.B. Ar

Human

Target Species

Reacts with: Mouse, Horse, Bovine, Rat, Goat, Cat, Sheep

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml		
Immunogen	P-selectin transfected 300.19 cells.		
External Database Links	UniProt: P16109 Related reagents P42201 Related reagents		
	Entrez Gene: 6403 SELP Related reagents 281486 SELP Related reagents 25651 Selp Related reagents		
Synonyms	GMRP, GRMP		
RRID	AB_566869		
Fusion Partners	Spleen cells from immunised CD62P knock-out mice (strain C57/B6) were fused with cells of the NS-1 myeloma cell line.		
Specificity	Mouse anti human CD62P, clone Psel.KO.2.7, recognizes human P-Selectin. CD62P is a ~140 kDa surface antigen expressed by activated platelets and endothelial cells, and plays an important role in adhesive processes between leucocytes and endothelial cells.		
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.		
References	 Massaguer, A. <i>et al.</i> (2000) Production and characterization of monoclonal antibodies against conserved epitopes of P-selectin (CD62P). <u>Tissue Antigens. 56 (2): 117-28.</u> Massaguer, A. <i>et al.</i> (2003) Characterization of platelet and soluble-porcine P-selectin (CD62P). <u>Vet Immunol Immunopathol. 96 (3-4): 169-81.</u> Massaguer, A. <i>et al.</i> (2002) Reactivity of CD62P workshop mAbs with resting and activated platelets from different animal species. In: Leucocyte Typing VII. Edited by Mason, D. <i>et al.</i> Oxford University Press, pp 342-3. Shirasuna, K. <i>et al.</i> (2012) Rapid accumulation of polymorphonuclear neutrophils in the Corpus luteum during prostaglandin F(2α)-induced luteolysis in the cow. <u>PLoS One. 7: e29054.</u> Johnson, C.A. Jr. <i>et al.</i> (2011) Platelet activation in ovines undergoing sham surgery or implant of the second generation PediaFlow pediatric ventricular assist device. <u>Artif Organs. 35: 602-13.</u> Iwaszko-Simonik, A. <i>et al.</i> (2015) Expression of surface platelet receptors (CD62P and 		

CD41/61) in horses with recurrent airway obstruction (RAO). <u>Vet Immunol Immunopathol.</u> 164 (1-2): 87-92.

- 7. Johnson, C.A. Jr *et al.* (2008) Flow cytometric assays for quantifying activated ovine platelets. Artif Organs. 32 (2): 136-45.
- 8. Johnson, C.A. Jr *et al.* (2011) Biocompatibility assessment of the first generation PediaFlow pediatric ventricular assist device. <u>Artif Organs. 35 (1): 9-21.</u>
- 9. Johnson, C.A. Jr *et al.* (2011) Platelet activation after implantation of the Levitronix PediVAS in the ovine model. ASAIO J. 57 (6): 516-21.
- 10. Shankarraman, V. *et al.* (2014) Biocompatibility Assessment of the CentriMag-Novalung Adult ECMO Circuit in a Model of Acute Pulmonary Hypertension. <u>ASAIO J. 60</u> (4): 429-35.
- 11. Trichler, S.A. *et al.* (2013) Ultra-pure platelet isolation from canine whole blood. <u>BMC</u> Vet Res. 9: 144.
- 12. Iwaszko-Simonik, A. *et al.* (2015) Expression of surface platelet receptors (CD62P and CD41/61) in horses with recurrent airway obstruction (RAO). <u>Vet Immunol Immunopathol.</u> 164 (1-2): 87-92.
- 13. Johnson, C.A. Jr. *et al.* (2011) Platelet activation in ovines undergoing sham surgery or implant of the second generation PediaFlow pediatric ventricular assist device. <u>Artif Organs.</u> 35 (6): 602-13.
- 14. ChanC, H.H. *et al.* (2017) Shear Stress-Induced Total Blood Trauma in Multiple Species. Artif Organs. 41 (10): 934-947.

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

12 months from date of despatch

Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information

Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA2419A488 10041

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL: Alexa Fluor® 488 (MCA928A488)

Recommended Useful Reagents

HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M376345:210128'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint

Email: antibody_sales_us@bio-rad.com