Datasheet: MCA2538A700 BATCH NUMBER 166932

Description:	MOUSE ANTI HUMAN CD79a:Alexa Fluor® 700
Specificity:	CD79a
Other names:	MB-1
Format:	ALEXA FLUOR® 700
Product Type:	Monoclonal Antibody
Clone:	HM57
lsotype:	lgG1
Quantity:	100 TESTS/1ml

Product Details

Applications	This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-</u> rad-antibodies com/protocols					
		Yes	No	Not Determined	Suggested Dilution	
	Flow Cytometry (1)	•			Neat	
	Where this product has not been tested for use in a particular technique this does not					
	necessarily exclude its use in such procedures. Suggested working dilutions are given as					
	a guide only. It is recon	nmended the	at the use	er titrates the product	for use in their own	
	system using appropria	ate negative/	positive c	ontrols.		
	(1) Membrane permea	abilization is	s require	d for this application	n. The use of	
	Leucoperm (Product	Code <u>BUF0</u>	<u>9</u>) is reco	ommended for this p	ourpose.	
Taurat On a sing						
larget Species	Human					
Species Cross	Human Reacts with: Mouse, R	abbit, Horse	, Pig, Moi	nkey, Rat, Bovine, Gu	uinea Pig, Fallow deer,	
Species Cross Reactivity	Reacts with: Mouse, R American Bison, Red d	abbit, Horse leer, Ferret,	, Pig, Moi Goat, Do	nkey, Rat, Bovine, Gu g	uinea Pig, Fallow deer,	
Species Cross Reactivity	Reacts with: Mouse, R American Bison, Red d N.B. Antibody reactivity	abbit, Horse leer, Ferret, y and workin	, Pig, Moi Goat, Do g conditio	nkey, Rat, Bovine, Gu g ons may vary betweer	uinea Pig, Fallow deer, n species. Cross	
Species Cross Reactivity	Reacts with: Mouse, Re American Bison, Red d N.B. Antibody reactivity reactivity is derived from	abbit, Horse leer, Ferret, y and workin m testing wit	, Pig, Moi Goat, Do g conditio hin our la	nkey, Rat, Bovine, Gu g ons may vary betweer boratories, peer-revie	uinea Pig, Fallow deer, n species. Cross ewed publications or	
Species Cross Reactivity	Reacts with: Mouse, Ra American Bison, Red d N.B. Antibody reactivity reactivity is derived from personal communication	abbit, Horse leer, Ferret, y and workin m testing wit ons from the	, Pig, Mor Goat, Do g conditio hin our la originato	nkey, Rat, Bovine, Gu g ons may vary betweer boratories, peer-revie rs. Please refer to refe	uinea Pig, Fallow deer, n species. Cross ewed publications or erences indicated for	
Species Cross Reactivity	Reacts with: Mouse, Ra American Bison, Red d N.B. Antibody reactivity reactivity is derived from personal communication further information.	abbit, Horse leer, Ferret, y and workin m testing wit ons from the	, Pig, Moi Goat, Do g conditio hin our la originato	nkey, Rat, Bovine, Gu g ons may vary betweer boratories, peer-revie rs. Please refer to refe	uinea Pig, Fallow deer, n species. Cross ewed publications or erences indicated for	
Species Cross Reactivity Product Form	Human Reacts with: Mouse, Ra American Bison, Red d N.B. Antibody reactivity reactivity is derived from personal communication further information.	abbit, Horse leer, Ferret, y and workin m testing wit ons from the d to Alexa Fl	, Pig, Moi Goat, Do g conditio hin our la originato uor 700 -	nkey, Rat, Bovine, Gu g ons may vary betweer boratories, peer-revie rs. Please refer to refe	uinea Pig, Fallow deer, n species. Cross ewed publications or erences indicated for	
Species Cross Reactivity Product Form Max Ex/Em	Human Reacts with: Mouse, Ra American Bison, Red d N.B. Antibody reactivity reactivity is derived from personal communication further information. Purified IgG conjugated Fluorophore	abbit, Horse leer, Ferret, y and workin m testing wit ons from the d to Alexa Fl Excitation N	, Pig, Mor Goat, Do g conditio hin our la originato uor 700 - lax (nm)	nkey, Rat, Bovine, Gu g ons may vary between boratories, peer-revie rs. Please refer to refe liquid Emission Max (nm)	uinea Pig, Fallow deer, n species. Cross ewed publications or erences indicated for	

Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml
Immunogen	Synthetic peptide corresponding to 202-216 amino acid sequence of human mb-1
External Database Links	UniProt: P11912 Related reagents Entrez Gene: 973 CD79A Related reagents
Synonyms	IGA, MB1
Fusion Partners	Spleen cells from immunized Balb/c mice were fused with cells of the Sp2/0 myeloma cell line
Specificity	Mouse anti Human CD79a antibody, clone HM57 recognizes an epitope within the cytoplasmic domain of CD79a. CD79a, also known as mb-1, is a 45 kDa protein that is expressed by B lymphocytes during differentiation from early pre-B cell stage through to plasma cells.
	The CD79a molecule associates with CD79b (B29) to form a heterodimer that is non-covalently linked to surface immunoglobulin, forming the B-cell receptor (BCR) complex. The CD79a/CD79b heterodimers are also necessary for intracellular signaling following antigen-binding to surface immunoglobulin.
Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^{6} cells in 100ul.
References	 Jones, M. <i>et al.</i> (1993) Detection of T and B cells in many animal species using cross-reactive anti-peptide antibodies. J Immunol. 150 (12): 5429-35. Nelson, D.D. <i>et al.</i> (2010) CD8(+)/perforin(+)/WC1(-) gammadelta T cells, not CD8(+) alphabeta T cells, infiltrate vasculitis lesions of American bison (<i>Bison bison</i>) with experimental sheep-associated malignant catarrhal fever. Vet Immunol Immunopathol. 136: 284-91. De Schauwer, C. <i>et al.</i> (2012) In search for cross-reactivity to immunophenotype equine mesenchymal stromal cells by multicolor flow cytometry. Cytometry A. 81 (4): 312-23. Spaas, J.H. <i>et al.</i> (2013) Culture and characterisation of equine peripheral blood mesenchymal stromal cells. Vet J. 195 (1): 107-13. Moore, P.F. <i>et al.</i> (2013) Canine inflamed nonepitheliotropic cutaneous T-cell lymphoma: a diagnostic conundrum. Vet Dermatol. 24 (1): 204-11.e44-5.

6. De Schauwer, C. *et al.* (2014) Characterization and profiling of immunomodulatory genes of equine mesenchymal stromal cells from non-invasive sources. <u>Stem Cell Res</u> <u>Ther. 5 (1): 6.</u>

7. Bozkurt, Y.A., *et al.* (2014) Histological and immunohistological studies of the structure of lymph nodes in Kilis goats. <u>Biotech Histochem. 89(6):440-5.</u>

8. Paebst, F. *et al.* (2014) Comparative immunophenotyping of equine multipotent mesenchymal stromal cells: an approach toward a standardized definition. <u>Cytometry A.</u> <u>85 (8): 678-87.</u>

9. Gelain ME *et al.* (2014) CD44 in canine leukemia: analysis of mRNA and protein expression in peripheral blood. <u>Vet Immunol Immunopathol. 159 (1-2): 91-6.</u>

10. Claessen, C. *et al.* (2015) Equid herpesvirus 1 (EHV1) infection of equine mesenchymal stem cells induces a pUL56-dependent downregulation of select cell surface markers. <u>Vet Microbiol. 176 (1-2): 32-9.</u>

11. Aresu, L. *et al.* (2015) Canine indolent and aggressive lymphoma: clinical spectrum with histologic correlation. <u>Vet Comp Oncol. 13 (4): 348-62.</u>

12. Poggi, A. *et al.* (2015) Flow cytometric evaluation of ki67 for the determination of malignancy grade in canine lymphoma. <u>Vet Comp Oncol. 13 (4): 475-80.</u>

13. Froment, R. & Bédard, C. (2016) Marked hyperphosphatasemia associated with an acute leukemia in a Great Dane. <u>Vet Clin Pathol. 45 (3): 459-65.</u>

14. Schinköthe J >*et al.* (2016) Characterization of tuberculous granulomas in different stages of progression and associated tertiary lymphoid tissue in goats experimentally infected with *Mycobacterium avium* subsp. *hominissuis*. <u>Comp Immunol Microbiol Infect</u> Dis. 47: 41-51.

15. Novacco, M. *et al.* (2016) Prognostic factors in canine acute leukaemias: a retrospective study. <u>Vet Comp Oncol. 14 (4): 409-16.</u>

16. Hillmann, A. *et al.* (2016) Comparative Characterization of Human and Equine Mesenchymal Stromal Cells: A Basis for Translational Studies in the Equine Model. <u>Cell</u> <u>Transplant. 25 (1): 109-24.</u>

17. Long, H. *et al.* (2016) Polyostotic Lymphoma in a Ferret (*Mustela putorius furo*). J Comp Pathol. 154 (4): 341-4.

18. Wessels, M. *et al.* (2017) Systemic necrotizing polyarteritis in three weaned lambs from one flock. J Vet Diagn Invest 29 (5) :733-37.

19. Uitterdijk, A. *et al.* (2017) Time course of VCAM-1 expression in reperfused myocardial infarction in swine and its relation to retention of intracoronary administered bone marrow-derived mononuclear cells. <u>PLoS One. 12 (6): e0178779.</u>

20. Nagata, K. *et al.* (2017) Epstein-Barr Virus Lytic Reactivation Activates B Cells Polyclonally and Induces Activation-Induced Cytidine Deaminase Expression: A Mechanism Underlying Autoimmunity and Its Contribution to Graves' Disease. <u>Viral</u> <u>Immunol. 30 (3): 240-9.</u>

21. Collins, J.J.P. *et al.* (2018) Impaired Angiogenic Supportive Capacity and Altered Gene Expression Profile of Resident CD146⁺ Mesenchymal Stromal Cells Isolated from Hyperoxia-Injured Neonatal Rat Lungs. <u>Stem Cells Dev. 27 (16): 1109-24.</u>

22. Murphy, E.G. *et al.* (2019) First detection of Hepatitis E virus (Orthohepevirus C) in wild brown rats (*Rattus norvegicus.*) from Great Britain. <u>Zoonoses Public Health. 66 (6)</u>: 686-94.

23. Skovdal, S.M. *et al.* (2019) Inhaled nebulized glatiramer acetate against Gram-negative bacteria is not associated with adverse pulmonary reactions in healthy,

	young adult female pigs. <u>PLoS One. 14 (10): e0223647.</u>
	24. Mu&mtildeoz-Silvestre, A. <i>et al.</i> (2020) Pathogenesis of Intradermal Staphylococcal
	Infections: Rabbit Experimental Approach to Natural <i>Staphylococcus aureus</i> Skin
	Infections. Am J Pathol. 190 (6): 1188-210.
	25. Forner, R. et al. (2021) Distribution difference of colostrum-derived B and T cells
	Subsets in gills and sows. <u>PLoS One. 16 (5): e0249366.</u>
	rat renal fibrosis followed by repeated injection of cisplatin. <u>J Vet Med Sci. 83 (9): 1435-42.</u>
	27. Carroll, C.S.E. et al. (2021) Simple and effective bacterial-based intratumoral cancer
	immunotherapy. <u>J Immunother Cancer.9(9):e002688.</u>
	28. Maciag, S.S. et al. (2022) The influence of source of porcine colostrum in
	development of early immune ontogeny in the piglet <u>Res Sq. Mar 24 [Epub ahead of print].</u>
	29. do Prado Duzanski, A. <i>et al.</i> (2022) Cell-mediated immunity and expression of MHC
	class I and class II molecules in dogs naturally infected by canine transmissible venereal
	tumor: Is there complete spontaneous regression outside the experimental CTVT? Res
	<u>Vet Sci. 145: 193-204.</u>
	30. Korbonits, L. et al. (2022) Mycobacterium avium subsp. paratuberculosis Infected
	Cows Reveal Divergent Immune Response in Bovine Peripheral Blood Derived
	Lymphocyte Proteome. <u>Metabolites. 12 (10): 924.</u>
	31. Schilloks, M.C. et al. (2023) Effects of GHR Deficiency and Juvenile Hypoglycemia on
	Immune Cells of a Porcine Model for Laron Syndrome. <u>Biomolecules. 13 (4): 597.</u>
	sz. dos Santos, M.C. et al. (2023) Effect of yeast extracted p-glucans of the infiniture
	periods Livestock Science, 275: 105280
	33 Haach V et al. (2023) A polyvalent virosomal influenza vaccine induces broad cellular.
	and humoral immunity in pigs Virol J 20 (1): 181
	34. Viitanen, S.J. <i>et al.</i> (2023) <i>Escherichia coli</i> -associated follicular cystitis in dogs:
	Clinical and pathologic characterization. J Vet Intern Med. 37 (3): 1059-66.
	35. Martini, V. <i>et al.</i> (2018) A retrospective study of flow cytometric characterization of
	suspected extranodal lymphomas in dogs. J Vet Diagn Invest, 30 (6): 830-6.
	36. Cha, S. <i>et al.</i> (2023) Non-B, Non-T Acute Lymphoblastic Leukemia in a Cat Journal of
	Veterinary Clinics. 40 (4): 298-302.
	37. Giese, I.M. et al. (2020) Chronic Hyperglycemia Drives Functional Impairment of
	Lymphocytes in Diabetic INS(C94Y) Transgenic Pigs. Front Immunol. 11: 607473.
	38. Keller, A.K. et al. (2019) POST-TRAUMATIC OCULAR LYMPHOMA IN THREE
	RABBITS (ORYCTOLAGUS CUNICULUS) Journal of Exotic Pet Medicine. 28: 154-61.
Further Reading	1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. <u>Vet Res. 39: 54.</u>
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.
	Avoid repeated freezing and thawing as this may denature the antibody. Storage in
	frost-free freezers is not recommended. This product is photosensitive and should be

	protected from light.	
Guarantee	12 months from date of despatch	
Acknowledgements	This product is provided under an intellectual property licence Corporation. The transfer of this product is contingent on the b product solely in research, excluding contract research or any and the buyer must not sell or otherwise transfer this product of diagnostic, therapeutic or prophylactic purposes; (b) testing, an services, or information in return for compensation on a per-tes or quality assurance or quality control, or (d) resale, whether or research. For information on purchasing a license to this product as described above, contact Life Technologies Corporation, 57 CA 92008 USA or outlicensing@thermofisher.com	from Life Technologies ouyer using the purchased fee for service research, or its components for (a) nalysis or screening st basis; (c) manufacturing or not resold for use in uct for purposes other than 791 Van Allen Way, Carlsbad
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA2538A700 10041	
Regulatory	For research purposes only	
Dolotod Droduc	to	

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL:Alexa Fluor® 700 (MCA928A700)

Recommended Useful Reagents

HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

North & South	Tel: +1 800 265 7376	Norldwide	Tel: +44 (0)1865 852 700	Europe	Tel: +49 (0) 89 8090 95 21
America	Fax: +1 919 878 3751		Fax: +44 (0)1865 852 739		Fax: +49 (0) 89 8090 95 50
	Email: antibody_sales_us@bio-rad.c	om	Email: antibody_sales_uk@bio-rad	.com	Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M416282:230222'

Printed on 08 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint