

# Datasheet: MCA336P BATCH NUMBER 154655

Description:	RAT ANTI MOUSE IgG1 HEAVY CHAIN:HRP
Specificity:	IgG1 HEAVY CHAIN
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	LO-MG1-2
Isotype:	IgG1
Quantity:	0.5 mg

# **Product Details**

## **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			•	
Immunohistology - Paraffin				
ELISA	•			500ng/ml
Western Blotting			•	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using approriate negative/positive controls.

Target Species	Mouse	
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid	
Preparation	Purified IgG prepared by affinity chromatography from tissue culture supernatant	
Buffer Solution	Phosphate buffered saline	
Preservative Stabilisers	0.01% Thiomersal 50% Glycerol	
Approx. Protein Concentrations	IgG concentration 1 mg/ml	

lmmunogen	Purified mouse IgG1 from BALB/c mice	
External Database	UniProt:	
	P01869 Related reagents	
	P01868 Related reagents	
	Entrez Gene:	
	16017 Ighg1 Related reagents	
	16017 Ighg1 Related reagents	
Synonyms	lgh-4	
RRID	AB_321957	
Fusion Partners	Spleen cells from immunised LOU/c rats were fused with cells of the rat IR983F myeloma cell line.	
Specificity	Rat anti Mouse IgG1 Heavy Chain antibody, clone LO-MG1-2 recognizes murine IgG1, and does not bind other mouse immunoglobulin classes or subclasses.	
ELISA	This product may be used in a direct ELISA or as a detection reagent in a sandwich ELISA together with MCA1289 as the capture reagent.	
References	<ol> <li>Song, J. et al. (2000) Heterogeneous distribution of isoactins in cultured vascular smooth muscle cells does not reflect segregation of contractile and cytoskeletal domains. J Histochem Cytochem. 48 (11): 1441-52.</li> <li>Denis, O. et al. (1993) Resting B cells can act as antigen presenting cells in vivo and induce antibody responses. Int Immunol. 5 (1): 71-8.</li> <li>Nakanishi, S. et al. (2010) Sequence analysis of a bacteriocinogenic plasmid of Clostridium butyricum and expression of the bacteriocin gene in Escherichia coli.</li> </ol>	
	Anaerobe. 16: 253-7.	

- 4. Echeverria, P.C. *et al.* (2006) Potent antigen-specific immunity to Toxoplasma gondii in adjuvant-free vaccination system using Rop2-Leishmania infantum Hsp83 fusion protein. Vaccine. 24: 4102-10.
- 5. Huang, C.H. *et al.* (2011) Airway inflammation and IgE production induced by dust mite allergen-specific memory/effector Th2 cell line can be effectively attenuated by IL-35. <u>J. Immunol.</u> 187: 462-71.
- 6. Agallou, M. *et al.* (2014) Experimental Validation of Multi-Epitope Peptides Including Promising MHC Class I- and II-Restricted Epitopes of Four Known Leishmania infantum Proteins. <u>Front Immunol. 5: 268.</u>
- 7. Doerfler, P.A. *et al.* (2015) BAFF Blockade Prevents Anti-Drug Antibody Formation in a Mouse Model of Pompe Disease. <u>Clin Immunol. pii: S1521-6616(15)00125-4.</u>
- 8. Ramos, J.D.A. *et al.* (2009) Characterization of Blo t 11 Monoclonal Antibodies with Constant Region Mutations Phil Sci Lett. 2(1): 38-48
- 9. Blackwell, N.M. & Else, K.J. (2002) A comparison of local and peripheral parasite-specific antibody production in different strains of mice infected with Trichuris muris. Parasite Immunol. 24 (4): 203-11.

- 10. Hall, G. *et al.* (2003) Suppression of allergen reactive Th2 mediated responses and pulmonary eosinophilia by intranasal administration of an immunodominant peptide is linked to IL-10 production. Vaccine. 21 (5-6): 549-61.
- 11. Hjerpe, C. *et al.* (2010) Dendritic cells pulsed with malondialdehyde modified low density lipoprotein aggravate atherosclerosis in Apoe(-/-) mice. <u>Atherosclerosis</u>. 209 (2): 436-41.
- 12. Kretschmer, B. *et al.* (2015) Anti-CD83 promotes IgG1 isotype switch in marginal zone B cells in response to TI-2 antigen. <u>Immunobiology</u>. 220 (8): 964-75.
- 13. Doerfler, P.A. *et al.* (2016) Copackaged AAV9 Vectors Promote Simultaneous Immune Tolerance and Phenotypic Correction of Pompe Disease. Hum Gene Ther. 27 (1): 43-59.
- 14. Kato, G. *et al.* (2014) β2 adrenergic agonist attenuates house dust mite-induced allergic airway inflammation through dendritic cells. BMC Immunol. 15: 39.
- 15. Margaroni, M. *et al.* (2017) Vaccination with poly(D,L-lactide-co-glycolide) nanoparticles loaded with soluble *Leishmania* antigens and modified with a TNFα-mimicking peptide or monophosphoryl lipid A confers protection against experimental visceral leishmaniasis. Int J Nanomedicine. 12: 6169-84.
- 16. DeGiovanni, C. *et al.* (2019) Cancer Vaccines Co-Targeting HER2/Neu and IGF1R. Cancers (Basel). 11 (4) Apr 11 [Epub ahead of print].
- 17. Doshi, B.S. *et al.* (2021) B cell-activating factor modulates the factor VIII immune response in hemophilia A. <u>J Clin Invest.</u> 131(8):142906.

#### Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Health And Safety
Information

12 months from date of despatch

Material Safety Datasheet documentation #10097 available at: https://www.bio-rad-antibodies.com/SDS/MCA336P

10097

Regulatory

Guarantee

For research purposes only

# Related Products

## **Recommended Useful Reagents**

AbGUARD® HRP STABILIZER PLUS (BUF052A)

AbGUARD® HRP STABILIZER PLUS (BUF052B)

AbGUARD® HRP STABILIZER PLUS (BUF052C)

TMB CORE (BUF056A)

TMB CORE+ (BUF062A)

TMB SIGNAL+ (BUF054A)

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

North & South Tel: +1 800 265 7376 America Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739

Europe

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Email: antibody\_sales\_us@bio-rad.com

Email: antibody\_sales\_uk@bio-rad.com

'M367587:200529'

Email: antibody\_sales\_de@bio-rad.com

## Printed on 29 Mar 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint