

Datasheet: MCA643F

Description:	MOUSE ANTI RAT CD44:FITC		
Specificity:	CD44		
Other names:	H-CAM, PGP-1		
Format:	FITC		
Product Type:	Monoclonal Antibody		
Clone:	OX-50		
Isotype:	lgG1		
Quantity:	0.1 mg		

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry				Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species Rat Product Form Purified IgG conjugated to Fluorescein Isothiocyanate Isomer
Product Form Purified IgG conjugated to Fluorescein Isothiocvanate Isomer
3 - , 3
ax Ex/Em Fluorophore Excitation Max (nm) Emission Max (n
FITC 490 525
paration Purified IgG prepared by affinity chromatography on Protein A
supernatant fer Solution Phosphate buffered saline

Immunogen	Rat T cell blasts.
External Database Links	UniProt: P26051 Related reagents Entrez Gene:
	25406 Cd44 Related reagents
RRID	AB_321722
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS0/1 myeloma cell line.
Specificity	Mouse anti Rat CD44 antibody, clone OX-50 recognizes the rat CD44 cell surface antigen, also known as Extracellular matrix receptor III, P90 lymphocyte homing/adhesion receptor, HUTCH-I, Hermes antigen, Hyaluronate receptor, Phagocytic glycoprotein 1, PGP-1 or Phagocytic glycoprotein I.
	CD44 is a 482 amino acid ~85 kDa single pass type I transmembrane glycoprotein, expressed by T cells, B cells, macrophages and thymocytes, with expression being increased following activation.
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	 Stevenson, K.S. <i>et al.</i> (2009) Isolation, characterization, and differentiation of thy1.1-sorted pancreatic adult progenitor cell populations. <u>Stem Cells Dev. 18:1389-98.</u> Jiang, T.S. <i>et al.</i> (2010) Reconstruction of the corneal epithelium with induced marrow mesenchymal stem cells in rats. <u>Mol Vis. 16: 1304-16.</u> Kanellis, J. <i>et al.</i> (2010) JNK signalling in human and experimental renal
	ischaemia/reperfusion injury. Nephrol Dial Transplant. 25: 2898-908.
	4. Li, S. <i>et al.</i> (2010) Upregulation of CXCR4 favoring neural-like cells migration via AKT activation. Neurosci Res. 67: 293-9.
	5. Stephens, L.A. et al. (2004) Phenotypic characterization of regulatory CD4+CD25+ T
	cells in rats. Int Immunol. 16: 365-75. 6. Rice, C.M. et al. (2010) Multipotent adult progenitor cell isolation and proliferation in cytokine and serum-free medium conditioned by rat B104 cells. Br J Haematol. 148: 441-4.
	7. Carulli, D. <i>et al</i> (2006) Composition of perineuronal nets in the adult rat cerebellum and
	the cellular origin of their components. <u>J Comp Neurol</u> . 494: 559-77.
	8. Galtrey, C.M. <i>et al.</i> (2008) Distribution and synthesis of extracellular matrix proteoglycans, hyaluronan, link proteins and tenascin-R in the rat spinal cord. <u>Eur J</u>
	Neurosci. 27: 1373-90.
	9. Hofmann, N. et al. (2002) Increased expression of ICAM-1, VCAM-1, MCP-1, and
	MIP-1 alpha by spinal perivascular macrophages during experimental allergic
	encephalomyelitis in rats. <u>BMC Immunol. 3: 11.</u> 10. Walther, M. <i>et al.</i> (2001) Exogenous antigen containing perivascular phagocytes
	10. Training, W. Ot an (2001) Enogonous uningent containing perivascular phagocytes

induce a non-encephalitogenic extravasation of primed lymphocytes. <u>J Neuroimmunol.</u>

117: 30-42.

- 11. Suzuki, A. *et al.* (2006) Localization of CD44 and hyaluronan in the synovial membrane of the rat temporomandibular joint. <u>Anat Rec A Discov Mol Cell Evol Biol. 288:</u> 646-52.
- 12. Goransson, V. *et al.* (2004) Renal hyaluronan accumulation and hyaluronan synthase expression after ischaemia-reperfusion injury in the rat. <u>Nephrol Dial Transplant.19:</u> 823-30.
- 13. Campbell, N.G. *et al.* (2016) Cell Size Critically Determines Initial Retention of Bone Marrow Mononuclear Cells in the Heart after Intracoronary Injection: Evidence from a Rat Model. <u>PLoS One. 11 (7): e0158232.</u>
- 14. Bejar, M.T. *et al.* (2016) Inhibition of Notch rescues the angiogenic potential impaired by cardiovascular risk factors in epicardial adipose stem cells. FASEB J. 30 (8): 2849-59.
- 15. Rochefort, G.Y. *et al.* (2006) Multipotential mesenchymal stem cells are mobilized into peripheral blood by hypoxia. Stem Cells. 24 (10): 2202-8.
- 16. Redondo, J. *et al.* (2015) Reductions in kinesin expression are associated with nitric oxide-induced axonal damage. J Neurosci Res. 93 (6): 882-92.
- 17. Huang, X. *et al.* (2019) MRI Tracking of SPIO- and *Fth1*-Labeled Bone Marrow Mesenchymal Stromal Cell Transplantation for Treatment of Stroke. <u>Contrast Media Mol Imaging</u>. 2019: 5184105.
- 18. Aminzadeh, A. *et al.* (2020) Investigating The Alterations of Oxidative Stress Status, Antioxidant Defense Mechanisms, MAP Kinase and Mitochondrial Apoptotic Pathway in Adipose-Derived Mesenchymal Stem Cells from STZ Diabetic Rats. <u>Cell J. 22 (Suppl 1):</u> 38-48.
- 19. Paiva, R.G. *et al.* (2020) Stem cells in end-to-side neurorrhaphy. Experimental study in rats Acta Cirúrgica Brasileira. 35 (12) [Epub ahead of print].
- 20. Hou, B. *et al.* (2018) Comparison of the Effects of BMSC-derived Schwann Cells and Autologous Schwann Cells on Remyelination Using a Rat Sciatic Nerve Defect Model. <u>Int J Biol Sci. 14 (13): 1910-22.</u>
- 21. Porwal, K. *et al.*kfz172 (2019) Increased bone marrow-specific adipogenesis by clofazimine causes impaired fracture healing, osteopenia and osteonecrosis without extraskeletal effects in rats. <u>Toxicol Sci. kfz172.</u>
- 22. Chang, H.H. *et al.* (2019) Intrarenal Transplantation of Hypoxic Preconditioned Mesenchymal Stem Cells Improves Glomerulonephritis through Anti-Oxidation, Anti-ER Stress, Anti-Inflammation, Anti-Apoptosis, and Anti-Autophagy. <u>Antioxidants (Basel). 9 (1): 2.</u>
- 23. Wu, J. *et al.* (2022) Reinforcing the function of bone graft via the Ca-P ceramics dynamic behavior-enhanced osteogenic microenvironment for optimal bone regeneration and reconstruction <u>Applied Materials Today</u>. 27: 101465.
- 24. Cheng, Y.H. *et al.* (2022) Intrarenal Arterial Transplantation of Dexmedetomidine Preconditioning Adipose Stem-Cell-Derived Microvesicles Confers Further Therapeutic Potential to Attenuate Renal Ischemia/Reperfusion Injury through miR-122-5p/Erythropoietin/Apoptosis Axis. <u>Antioxidants (Basel)</u>. 11(9): 1702.

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA643F 10041
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL:FITC (MCA1209F)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M385312:210513'

Printed on 14 Apr 2025

© 2025 Bio-Rad Laboratories Inc | Legal | Imprint