

# Datasheet: MCA87B BATCH NUMBER 162168

Description:	MOUSE ANTI HUMAN CD45:Biotin		
Specificity:	CD45		
Other names:	LCA		
Format:	Biotin		
<b>Product Type:</b>	Monoclonal Antibody		
Clone:	F10-89-4		
Isotype:	lgG2a		
Quantity:	0.1 mg		

#### **Product Details**

#### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	-			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml
lmmunogen	Human T lymphocytes.

## **External Database UniProt:** Links P08575 Related reagents **Entrez Gene:** 5788 PTPRC Related reagents **Synonyms CD45 RRID** AB\_323920 **Fusion Partners** Spleen cells from immunized BALB/c mice were fused with cells of the mouse NS-1 myeloma cell line. **Specificity**

Mouse anti Human CD45 antibody, clone F10-89-4 recognizes the human CD45 cell

surface antigen, also known as leucocyte common antigen (LCA). CD45 is a complex molecule existing in a number of isoforms.

Antibodies recognizing a common epitope on all of these isoforms are termed CD45 whilst those recognizing only individual isoforms are termed CD45RA or CD45RO etc.

Mouse anti Human CD45 antibody, clone F10-89-4 reacts with all forms of CD45 expressed by all haematopoietic cells, except erythrocytes, having a higher level of expression on lymphocytes than on granulocytes. It is routinely tested in flow cytometry on human peripheral blood leucocytes.

Mouse anti Human CD45 antibody, clone F10-89-4, has been validated for use on the Genesis Cell Isolation System with the CelSelect Slide<sup>TM</sup> technology.

#### **Flow Cytometry**

Use 10ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

#### References

- 1. Quenby, S et al. (1999) Pre-implantation endometrial leukocytes in women with recurrent miscarriage. Human Reprod. 14(9):2386-2391.
- 2. Hauser, P.V. et al. (2010) Stem cells derived from human amniotic fluid contribute to acute kidney injury recovery. Am J Pathol. 177: 2011-21.
- 3. Mallam, E. et al. (2010) Characterization of in vitro expanded bone marrow-derived mesenchymal stem cells from patients with multiple sclerosis. Mult Scler. 16: 909-18.
- 4. Marrinucci, D. et al. (2010) Cytomorphology of circulating colorectal tumor cells:a small case series. J Oncol. 2010: 861341.
- 5. Paul, G. et al. (2012) The adult human brain harbors multipotent perivascular mesenchymal stem cells. PLoS One. 7: e35577.
- 6. De Schauwer, C. et al. (2012) In search for cross-reactivity to immunophenotype equine mesenchymal stromal cells by multicolor flow cytometry. Cytometry A. 81 (4): 312-23.
- 7. Kazane, S.A. et al. (2012) Site-specific DNA-antibody conjugates for specific and sensitive immuno-PCR. Proc Natl Acad Sci U S A. 109: 3731-6.
- 8. Spaas, J.H. et al. (2013) Culture and characterisation of equine peripheral blood mesenchymal stromal cells. Vet J. 195 (1): 107-13.
- 9. Sadarangani, A. et al. (2015) GLI2 inhibition abrogates human leukemia stem cell

- dormancy. J Transl Med. 13: 98.
- 10. Gunawardene, P. *et al.* (2015) Association Between Circulating Osteogenic Progenitor Cells and Disability and Frailty in Older Persons: The Nepean Osteoporosis and Frailty Study. <u>J Gerontol A Biol Sci Med Sci. pii: glv190.</u>
- 11. Mohamed Suhaimi, N.A. *et al.* (2015) Non-invasive sensitive detection of KRAS and BRAF mutation in circulating tumor cells of colorectal cancer patients. <u>Mol Oncol. 9 (4):</u> 850-60.
- 12. Ruiz, C. *et al.* (2015) Limited genomic heterogeneity of circulating melanoma cells in advanced stage patients. <u>Phys Biol. 12 (1): 016008</u>.
- 13. Gogoi P *et al.* (2016) Development of an Automated and Sensitive Microfluidic Device for Capturing and Characterizing Circulating Tumor Cells (CTCs) from Clinical Blood Samples. PLoS One. 11 (1): e0147400.
- 14. Gomiero, C. *et al.* (2016) Tenogenic induction of equine mesenchymal stem cells by means of growth factors and low-level laser technology. Vet Res Commun. 40 (1): 39-48.
- 15. Bianchessi, M. *et al.* (2016) Effect of Fibroblast Growth Factor 2 on Equine Synovial Fluid Chondroprogenitor Expansion and Chondrogenesis. Stem Cells Int. 2016: 9364974.
- 16. Branly, T. *et al.* (2017) Characterization and use of Equine Bone Marrow Mesenchymal Stem Cells in Equine Cartilage Engineering. Study of their Hyaline Cartilage Forming Potential when Cultured under Hypoxia within a Biomaterial in the Presence of BMP-2 and TGF-β1. <u>Stem Cell Rev Rep. 13 (5): 611-30.</u>
- 17. GarikipatiV, N.S. *et al.* (2018) Isolation and characterization of mesenchymal stem cells from human fetus heart. PLoS One. 13 (2): e0192244.
- 18. Shishido, S.N. *et al.* (2019) Circulating tumor cells as a response monitor in stage IV non-small cell lung cancer. J Transl Med. 17 (1): 294.
- 19. Welter, L. *et al.* (2020) Treatment response and tumor evolution: lessons from an extended series of multianalyte liquid biopsies in a metastatic breast cancer patient. <u>Cold Spring Harb Mol Case Stud.</u> 6 (6): a005819.
- 20. Ndacayisaba, L.J. *et al.* (2022) Enrichment-Free Single-Cell Detection and Morphogenomic Profiling of Myeloma Patient Samples to Delineate Circulating Rare Plasma Cell Clones <u>Curr Oncol. 29 (5): 2954-72.</u>
- 21. Shishido, S.N. *et al.* (2022) Liquid Biopsy Landscape in Patients with Primary Upper Tract Urothelial Carcinoma. Cancers (Basel). 14 (12): 3007.
- 22. Chai, S. *et al.* (2022) Identification of epithelial and mesenchymal circulating tumor cells in clonal lineage of an aggressive prostate cancer case. NPJ Precis Oncol. 6 (1): 41.
- 23. Zhu, J. *et al.* (2022) Sequential Method for Analysis of CTCs and Exosomes from the Same Sample of Patient Blood. <u>ACS Omega. 7 (42): 37581-88.</u>
- 24. Setayesh, S.M. *et al.* (2022) Multianalyte liquid biopsy to aid the diagnostic workup of breast cancer. NPJ Breast Cancer. 8 (1): 112.
- 25. Ndacayisaba, L.J. *et al.* (2022) Characterization of BCMA Expression in Circulating Rare Single Cells of Patients with Plasma Cell Neoplasms. <u>Int J Mol Sci. 23 (21): 13427.</u>
- 26. Qi, E. *et al.* (2023) Investigation of liquid biopsy analytes in peripheral blood of individuals after SARS-CoV-2 infection. <u>EBioMedicine</u>. 90: 104519.
- 27. Seo, J. *et al.* (2023) Plasticity of circulating tumor cells in small cell lung cancer. <u>Sci</u> Rep. 13 (1): 11775.
- 28. Setayesh, S.M. *et al.* (2023) Targeted single-cell proteomic analysis identifies new liquid biopsy biomarkers associated with multiple myeloma. <u>NPJ Precis Oncol. 7 (1): 95.</u>
- 29. Welter, L. et al. (2023) Cell State and Cell Type: Deconvoluting Circulating Tumor Cell

Populations in Liquid Biopsies by Multi-Omics. Cancers (Basel). 15 (15): 3949. 30. Shishido, S.N. et al. (2024) Cancer-related cells and oncosomes in the liquid biopsy of pancreatic cancer patients undergoing surgery. NPJ Precis Oncol. 8 (1): 36. Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. Guarantee 12 months from date of despatch **Health And Safety** Material Safety Datasheet documentation #10041 available at: Information https://www.bio-rad-antibodies.com/SDS/MCA87B 10041

### Related Products

#### **Recommended Useful Reagents**

**HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)** 

North & South Tel: +1 800 265 7376 America

Regulatory

Fax: +1 919 878 3751

Worldwide

For research purposes only

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody\_sales\_uk@bio-rad.com

Tel: +49 (0) 89 8090 95 21 Europe Fax: +49 (0) 89 8090 95 50

Email: antibody\_sales\_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M381501:210512'

#### Printed on 17 Apr 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint

Email: antibody\_sales\_us@bio-rad.com