

## Datasheet: MCA87PECY5.5

|                      |                                 |
|----------------------|---------------------------------|
| <b>Description:</b>  | MOUSE ANTI HUMAN CD45:RPE-Cy5.5 |
| <b>Specificity:</b>  | CD45                            |
| <b>Other names:</b>  | LCA                             |
| <b>Format:</b>       | RPE-CY5.5                       |
| <b>Product Type:</b> | Monoclonal Antibody             |
| <b>Clone:</b>        | F10-89-4                        |
| <b>Isotype:</b>      | IgG2a                           |
| <b>Quantity:</b>     | 100 TESTS                       |

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪   |    |                | Neat               |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

|                                 |  |                            |                          |
|---------------------------------|--|----------------------------|--------------------------|
| <b>Target Species</b>           | Human  |                            |                          |
| <b>Product Form</b>             | Purified IgG conjugated to R. Phycoerythrin-Cy5.5 (RPE-Cy5.5) - lyophilized  |                            |                          |
| <b>Reconstitution</b>           | Reconstitute with 1.0 ml distilled water<br>Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution |                            |                          |
| <b>Max Ex/Em</b>                | <b>Fluorophore</b>   | <b>Excitation Max (nm)</b> | <b>Emission Max (nm)</b> |
|                                 | RPE-Cy5.5 488nm laser  | 496                        | 695                      |
| <b>Preparation</b>              | Purified IgG prepared by affinity chromatography on Protein A  |                            |                          |
| <b>Buffer Solution</b>          | Phosphate buffered saline  |                            |                          |
| <b>Preservative Stabilisers</b> | 0.09% Sodium Azide<br>1% Bovine Serum Albumin  |                            |                          |

5% Sucrose

---

|                  |                      |
|------------------|----------------------|
| <b>Immunogen</b> | Human T lymphocytes. |
|------------------|----------------------|

---

|                                |  |
|--------------------------------|--|
| <b>External Database Links</b> | <b>UniProt:</b><br><a href="#">P08575</a> <a href="#">Related reagents</a>         |
|                                | <b>Entrez Gene:</b><br><a href="#">5788</a> PTPRC <a href="#">Related reagents</a> |

---

|                 |      |
|-----------------|------|
| <b>Synonyms</b> | CD45 |
|-----------------|------|

---

|             |             |
|-------------|-------------|
| <b>RRID</b> | AB_10844351 |
|-------------|-------------|

---

|                        |  |
|------------------------|--|
| <b>Fusion Partners</b> | Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS-1 myeloma cell line. |
|------------------------|--|

---

|                    |   |
|--------------------|---|
| <b>Specificity</b> | <p><b>Mouse anti Human CD45 antibody, clone F10-89-4</b> recognizes the human CD45 cell surface antigen, also known as leucocyte common antigen (LCA). CD45 is a complex molecule existing in a number of isoforms.</p> <p>Antibodies recognizing a common epitope on all of these isoforms are termed CD45 whilst those recognizing only individual isoforms are termed CD45RA or CD45RO etc.</p> <p>Mouse anti Human CD45 antibody, clone F10-89-4 reacts with all forms of CD45 expressed by all haematopoietic cells, except erythrocytes, having a higher level of expression on lymphocytes than on granulocytes.</p> <p>Mouse anti Human CD45 antibody, clone F10-89-4 is routinely tested in flow cytometry on human peripheral blood leucocytes.</p> |
|--------------------|---|

---

|                       |  |
|-----------------------|--|
| <b>Flow Cytometry</b> | Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells or 100ul whole blood |
|-----------------------|--|

---

|                   |   |
|-------------------|---|
| <b>References</b> | <ol style="list-style-type: none"><li>1. Dalchau, R. <i>et al.</i> (1980) Monoclonal antibody to a human leukocyte-specific membrane glycoprotein probably homologous to the leukocyte-common (L-C) antigen of the rat. <a href="#">Eur J Immunol. 10 (10): 737-44.</a></li><li>2. Quenby, S <i>et al.</i> (1999) Pre-implantation endometrial leukocytes in women with recurrent miscarriage. <a href="#">Human Reprod. 14(9):2386-2391.</a></li><li>3. Hauser, P.V. <i>et al.</i> (2010) Stem cells derived from human amniotic fluid contribute to acute kidney injury recovery. <a href="#">Am J Pathol. 177: 2011-21.</a></li><li>4. Mallam, E. <i>et al.</i> (2010) Characterization of <i>in vitro</i> expanded bone marrow-derived mesenchymal stem cells from patients with multiple sclerosis. <a href="#">Mult Scler. 16: 909-18.</a></li><li>5. Marrinucci, D. <i>et al.</i> (2010) Cytomorphology of circulating colorectal tumor cells:a small case series. <a href="#">J Oncol. 2010: 861341.</a></li><li>6. Kazane, S.A. <i>et al.</i> (2012) Site-specific DNA-antibody conjugates for specific and sensitive immuno-PCR. <a href="#">Proc Natl Acad Sci U S A. 109: 3731-6.</a></li><li>7. Paul, G. <i>et al.</i> (2012) The adult human brain harbors multipotent perivascular</li></ol> |
|-------------------|---|

- mesenchymal stem cells. [PLoS One. 7: e35577.](#)
8. Sadarangani, A. *et al.* (2015) GLI2 inhibition abrogates human leukemia stem cell dormancy. [J Transl Med. 13: 98.](#)
  9. Gunawardene, P. *et al.* (2015) Association Between Circulating Osteogenic Progenitor Cells and Disability and Frailty in Older Persons: The Nepean Osteoporosis and Frailty Study. [J Gerontol A Biol Sci Med Sci. pii: glv190.](#)
  10. Gogoi P *et al.* (2016) Development of an Automated and Sensitive Microfluidic Device for Capturing and Characterizing Circulating Tumor Cells (CTCs) from Clinical Blood Samples. [PLoS One. 11 \(1\): e0147400.](#)
  11. Spaas, J.H. *et al.* (2013) Culture and characterisation of equine peripheral blood mesenchymal stromal cells. [Vet J. 195 \(1\): 107-13.](#)
  12. Gomiero, C. *et al.* (2016) Tenogenic induction of equine mesenchymal stem cells by means of growth factors and low-level laser technology. [Vet Res Commun. 40 \(1\): 39-48.](#)
  13. De Schauwer, C. *et al.* (2012) In search for cross-reactivity to immunophenotype equine mesenchymal stromal cells by multicolor flow cytometry. [Cytometry A. 81 \(4\): 312-23.](#)
  14. Bianchessi, M. *et al.* (2016) Effect of Fibroblast Growth Factor 2 on Equine Synovial Fluid Chondroprogenitor Expansion and Chondrogenesis. [Stem Cells Int. 2016: 9364974.](#)
  15. Mohamed Suhaimi, N.A. *et al.* (2015) Non-invasive sensitive detection of KRAS and BRAF mutation in circulating tumor cells of colorectal cancer patients. [Mol Oncol. 9 \(4\): 850-60.](#)
  16. Ruiz, C. *et al.* (2015) Limited genomic heterogeneity of circulating melanoma cells in advanced stage patients. [Phys Biol. 12 \(1\): 016008.](#)
  17. Branly, T. *et al.* (2017) Characterization and use of Equine Bone Marrow Mesenchymal Stem Cells in Equine Cartilage Engineering. Study of their Hyaline Cartilage Forming Potential when Cultured under Hypoxia within a Biomaterial in the Presence of BMP-2 and TGF- $\beta$ 1. [Stem Cell Rev. Jun 09 \[Epub ahead of print\].](#)
  18. GarikipatiV, N.S. *et al.* (2018) Isolation and characterization of mesenchymal stem cells from human fetus heart. [PLoS One. 13 \(2\): e0192244.](#)
  19. Shishido, S.N. *et al.* (2019) Circulating tumor cells as a response monitor in stage IV non-small cell lung cancer. [J Transl Med. 17 \(1\): 294.](#)
  20. Welter, L. *et al.* (2020) Treatment response and tumor evolution: Lessons from an extended series of multi-analyte liquid biopsies in a metastatic breast cancer patient. [Cold Spring Harb Mol Case Stud. Nov 17 \[Epub ahead of print\].](#)

---

|                                      |  |
|--------------------------------------|--|
| <b>Storage</b>                       | <p>Prior to reconstitution store at +4°C.</p> <p>After reconstitution store at +4°C.</p> <p>DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use</p> |
| <b>Guarantee</b>                     | 12 months from date of despatch  |
| <b>Acknowledgements</b>              | Cy® and CyDye® are registered trademarks of GE Healthcare  |
| <b>Health And Safety Information</b> | Material Safety Datasheet documentation #20487 available at: 20487: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/20487.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/20487.pdf</a>  |

---

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M375685:210104'

**Printed on 22 Oct 2021**