

Datasheet: MCA890F BATCH NUMBER 166134

Description:	MOUSE ANTI HUMAN CYTOKERATIN 14:FITC		
Specificity:	CYTOKERATIN 14		
Format:	FITC		
Product Type:	Monoclonal Antibody		
Clone:	LL002		
Isotype:	IgG3		
Quantity:	0.1 mg		

Product Details

Applications	This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <u>www.bio-rad-antibodies.com/protocols</u> .						
		Yes No	Not Determined	Suggested Dilution			
	Flow Cytometry (1)	•		Neat - 1/10			
	Immunofluorescence	-					
Target Species	necessarily exclude it a guide only. It is reco system using appropr (1) Membrane perme	s use in such proce mmended that the iate negative/positive abilization is requ	user titrates the antibody	ng dilutions are given as y for use in their own n. The use of			
Species Cross Reactivity	Reacts with: Elephant, Dog, Pig, Lion N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.						
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid						
Max Ex/Em	Fluorophore FITC	Excitation Max (nr 490	n) Emission Max (nm) 525				
Preparation	Purified IgG prepared	by affinity chromat	ography on Protein A fro	m tissue culture			

	supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1mg/ml
Immunogen	Last 15 C-terminal residues of human cytokeratin 14 conjugated to thyroglobulin.
External Database Links	UniProt: <u>P02533</u> <u>Related reagents</u> Entrez Gene: <u>3861</u> KRT14 <u>Related reagents</u>
RRID	AB_872024
Specificity	Mouse anti Human Cytokeratin 14 antibody, clone LL002 recognizes cytokeratin 14, a type I intermediate filament, expressed by stratifying epithelial cells and can be used to distinguish these cell types from simple epithelial cells, which do not express cytokeratin 14. Mouse anti cytokeratin 14, clone LL002 has been reported to be suitable for use in Western blotting (<u>Alam <i>et al.</i> 2011</u>)
Flow Cytometry	Use 10ul of the suggested working dilution to label 1×10^6 cells in 100ul.
References	 Wetzels, R.H. <i>et al.</i> (1989) Detection of basement membrane components and basal cell keratin 14 in noninvasive and invasive carcinomas of the breast. <u>Am J Pathol. 134 (3)</u>: <u>571-9.</u> Lane, E.B. & Alexander, C.M. (1990) Use of keratin antibodies in tumor diagnosis. <u>Semin Cancer Biol. 1 (3)</u>: 165-79. Munz, B. <i>et al.</i> (1999) Overexpression of activin A in the skin of transgenic mice reveals new activities of activin in epidermal morphogenesis, dermal fibrosis and wound repair. <u>EMBO J. 18: 5205-15.</u> Stumpf, P. & Welsch, U. (2002) Cutaneous eccrine glands of the foot pads of the rock hyrax (<i>Procavia capensis.</i>, <i>Hyracoidea, Mammalia.</i>). <u>Cells Tissues Organs. 171 (2-3)</u>: <u>215-26.</u> Stumpf, P. <i>et al.</i> (2004) Cutaneous eccrine glands of the foot pads of the small Madagascan tenrec (<i>Echinops telfairi, Insectivora, Tenrecidae.</i>): skin glands in a primitive mammal. <u>Cell Tissue Res. 315 (1): 59-70.</u> Varley, C.L. <i>et al.</i> (2004) Activation of peroxisome proliferator-activated receptor-gamma reverses squamous metaplasia and induces transitional differentiation in normal human urothelial cells. <u>Am J Pathol. 164: 1789-98.</u>

7. Richardson, G.D. *et al.* (2004) CD133, a novel marker for human prostatic epithelial stem cells. J Cell Sci. 117 (Pt 16): 3539-45.

8. Stumpf, P. & Welsch, U. (2004) Secretory and defensive functions of the duct system of the lactating mammary gland of the African elephant (*Loxodonta africana*, Proboscidea) <u>Zoomorphology. 123 (3): 155–67</u>

9. Collins, A.T. *et al*et al. (2005) Prospective identification of tumorigenic prostate cancer stem cells. <u>Cancer Res. 65: 10946-51.</u>

10. Matos, A.J. *et al.* (2006) Detection of lymph node micrometastases in malignant mammary tumours in dogs by cytokeratin immunostaining. <u>Vet Rec. 158: 626-30.</u>

11. Faustino, A.M. & Dias Pereira, P. (2007) A salivary malignant myoepithelioma in a dog. <u>Vet J. 173 (1): 223-6.</u>

12. Holliday, D. *et al.* (2009) Novel multicellular organotypic models of normal and malignant breast: tools for dissecting the role of the microenvironment in breast cancer progression. <u>Breast Cancer Res. 11: R3</u>

13. Eastman, R. .Jr. *et al.* (2010) Fibroblast growth factor-10 signals development of von Brunn's nests in the exstrophic bladder. <u>Am J Physiol Renal Physiol.299:F1094-110.</u>

14. Takahashi, C. *et al.* (2010) Newly established cell lines from mouse oral epithelium regenerate teeth when combined with dental mesenchyme. <u>In Vitro Cell Dev Biol Anim.</u> 46: 457-68.

15. Alam H *et al.* (2011) Loss of keratins 8 and 18 leads to alterations in α 6 β 4-integrinmediated signalling and decreased neoplastic progression in an oral-tumour-derived cell line. J Cell Sci. 124 (Pt 12): 2096-106.

16. Abdeen, S.K. *et al.* (2011) Wwox inactivation enhances mammary tumorigenesis. Oncogene. 30: 3900-6.

17. Clark, S.E. *et al.* (2011) Molecular subtyping of DCIS: heterogeneity of breast cancer reflected in pre-invasive disease. <u>Br J Cancer. 104: 120-7.</u>

18. Turner, A. *et al.* (2011) Transplantation of autologous differentiated urothelium in an experimental model of composite cystoplasty. <u>Eur Urol. 59: 447-54.</u>

19. Mwase, M. *et al.* (2013) Cutaneous Squamous Cell Carcinoma presenting as a Wound with Discharging Sinus Tracts in a Wild African Lion (*Panthera leo*). <u>J Comp Pathol. pii:</u> <u>S0021-9975(13)00106-0</u>.

20. Caceres, S. *et al.* (2015) Establishment and Characterization of a New Cell Line of Canine Inflammatory Mammary Cancer: IPC-366. <u>PLoS One. 10 (3): e0122277.</u>

21. Honda, K. & Tomooka, Y. (2016) Nerve-independent and ectopically additional induction of taste buds in organ culture of fetal tongues. *In Vitro* Cell Dev Biol Anim. 52 (9): 911-9.

22. Ogihara, K. & Madarame, H. (2020) Pancreatic adenosquamous carcinoma with invasion to the spleen in a cat. <u>J Vet Med Sci. 82 (9): 1395-9.</u>

23. Hustler, A. *et al.* (2018) Differential transcription factor expression by human epithelial cells of buccal and urothelial derivation. <u>Exp Cell Res. 369 (2): 284-294.</u>

24. Sakai, Y. *et al.* (2019) A clonal stem cell line established from a mouse mammary placode with ability to generate functional mammary glands. <u>In Vitro Cell Dev Biol Anim.</u> <u>55 (10): 861-71.</u>

25. Yasuno, K. *et al.* (2013) Atypical canine mammary adenoma characterized by cystic ducts comprising a single layer of basaloid cells with myoepithelial differentiation. <u>J Vet Med Sci. 75 (8): 1095-9.</u>

Further Reading	1. Moll, R. <i>et al.</i> (1982) The catalog of human cytokeratins: patternormal epithelia, tumors and cultured cells. <u>Cell. 31 (1): 11-24.</u>	erns of expression in
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.	
	Avoid repeated freezing and thawing as this may denature the a frost-free freezers is not recommended. This product is photose protected from light.	, ,
Guarantee	12 months from date of despatch	
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: https://www.bio-rad-antibodies.com/SDS/MCA890F 10041	
Regulatory	For research purposes only	

Related Products

Recommended Useful Reagents

LEUCOPERM (BUF09) HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

North & South	Tel: +1 800 265 7376	Worldwide	Tel: +44 (0)1865 852 700	Europe	Tel: +49 (0) 89 8090 95 21
America	Fax: +1 919 878 3751		Fax: +44 (0)1865 852 739		Fax: +49 (0) 89 8090 95 50
	Email: antibody_sales_us@b	io-rad.com	Email: antibody_sales_uk@bio	o-rad.com	Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M423276:231010'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint