

Datasheet: MCA923FT

Description:	MOUSE ANTI HUMAN CD49d:FITC
Specificity:	CD49d
Other names:	INTEGRIN ALPHA 4 CHAIN, VLA-4
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	44H6
lsotype:	lgG1
Quantity:	25 μg

Product Details

RRID	AB_2130044				
Applications	from testing within our	laboratories, pee refer to reference ase visit <u>www.bio</u>	r-reviewe es indica -rad-antik		al communications from n. For general protocol
		Yes	No	Not Determined	Suggested Dilution
	Flow Cytometry				Neat
	exclude its use in such	procedures. Sug user titrates the	ggested v	In a particular technique vorking dilutions are give for use in their own syste	
Target Species	Human				
Product Form	Purified IgG conjugate	d to Fluorescein	lsothiocy	anate Isomer 1 (FITC) - I	iquid
Max Ex/Em	Fluorophore FITC	Excitation Max (n 490	m) Emi	ssion Max (nm) 525	
Preparation	Purified IgG prepared I	by affinity chroma	atography	on Protein G from tissue	e culture supernatant
Buffer Solution	Phosphate buffered sa	line			
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum A	Albumin			
Approx. Protein Concentrations	IgG concentration 0.1	mg/ml			
Immunogen	HOON pre-B leukaemi	a cell line.			
External Database Links	UniProt:				

P13612 Related reagents

	Entrez Gene: <u>3676</u> ITGA4 <u>Related reagents</u>
Synonyms	CD49D
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the SP2.0-Ag14 mouse myeloma cell line.
Specificity	Mouse anti Human CD49d monoclonal antibody, clone 44H6 recognizes human CD49d also known as integrin alpha-4 or VLA-4 subunit alpha. CD49d is a ~150kDa single pass type 1 transmembrane glycoprotein with seven FG-GAP repeats, characteristic of alpha integrins, in its extracellular domain. CD49d can be proteolytically cleaved to yield fragments of 80 and 70 kDa (Hemler <i>et al.</i> 1987). CD49d associates with either CD29 to form VLA-4 or with Integrin beta-7 to form The Peyer's patches-specific homing receptor LPAM-1, involved in the lymphocyte migration and homing to gut-associated lymphoid tissue (Sackstein 2006) through its interaction with MadCam-1, preferentially expressed on Peyer's patch high endothelial venules and postcapillary venules in lamina propria (Briskin <i>et al.</i> 1997)
	Mouse anti human CD49d, clone 44H6 is reported to induce homotypic aggregation of pre-B cell lines including <u>HOON</u> and <u>NALM-6</u> (<u>Letarte <i>et al.</i> 1993</u>).
Flow Cytometry	Use 10ul of the antibody to label 100ul whole blood or 10 ⁶ cells in 100ul.
References	 Quackenbush, E.J. & Letarte, M. (1985) Identification of several cell surface proteins of non-T, non-B acute lymphoblastic leukemia by using monoclonal antibodies. J Immunol. 134 (2): 1276-85. Letarte, M. <i>et al.</i> (1993) Homotypic aggregation of pre-B leukemic cell lines by antibodies to VLA integrins correlates with their expression of CD9. Leukemia 7(1): 93-103. Hartz B <i>et al.</i> (2011) α 4 integrin levels on mobilized peripheral blood stem cells predict rapidity of engraftment in patients receiving autologous stem cell transplantation. Blood. 118 (8): 2362-5. Ling, K.L. <i>et al.</i> (2007) Modulation of CD103 expression on human colon carcinoma-specific CTL. JImmunol. 178: 2908-15. Eigenmann, PA. <i>et al.</i> (1999) The mucosal adhesion receptor alpha4beta7 integrin is selectively increased in lymphocytes stimulated with beta-lactoglobulin in children allergic to cow's milk. J Allergy Clin Immunol. 103: 931-6. Rabe, H. <i>et al.</i> (2011) Higher proportions of circulating FOXP3+ and CTLA-4+ regulatory T cells are associated with lower fractions of memory CD4+ T cells in infants. J Leukoc Biol. 90: 1133-40. Seoh, J.Y. <i>et al.</i> (2001) Cell cycling status of human cord blood CD34+ cells during ex vivo expansion is related to the level of very late antigen expression. J Korean Med Sci. 16: 20-4. Seoh, J.Y. (2001) VLA-4 expression and cell cycling status during ex vivo expansion of human cord blood CD34+ cells. Acta Haematol. 105: 111-5. Dogan, A. <i>et al.</i> (1997) Expression of lymphocyte homing receptors and vascular addressins in low-grade gastric B-cell lymphomas of mucosa-associated lymphoid tissue. Am J Pathol. 151: 1361-9. Dogan, A. <i>et al.</i> (1998) Adhesion molecules in iris biopsy specimens from patients with uveitis. Br J Ophthalmol. 82: 432-7. Humphries, J.D. and Humphries, M.J. (2007) CD14 is a ligand for the integrin alpha4beta1. FEBS Lett. 581: 757-63. Grindeb

	responses by combinations of monoclonal antibodies specific to CD2 and NKG2D and extracellular
	signal-regulated kinase kinase inhibitor. <u>Immunology. 130: 545-55.</u>
	14. Stadlmann, S. et al. (2003) Disruption of the integrity of human peritoneal mesothelium by
	interleukin-1beta and tumor necrosis factor-alpha. Virchows Arch. 443: 678-85.
	15. Paul, G. <i>et al.</i> (2012) The adult human brain harbors multipotent perivascular mesenchymal stem cells. <u>PLoS One. 7: e35577.</u>
	16. Spring, F.A. <i>et al.</i> (2001) Intercellular adhesion molecule-4 binds alpha(4)beta(1) and alpha(V)-
	family integrins through novel integrin-binding mechanisms. <u>Blood. 98: 458-66.</u>
	17. Woolhouse, I.S. <i>et al.</i> (2005) Endothelial interactions of neutrophils under flow in chronic
	obstructive pulmonary disease. Eur Respir J. 25: 612-7.
	18. McGilvray, I.D. <i>et al.</i> (1997) VLA-4 integrin cross-linking on human monocytic THP-1 cells
	induces tissue factor expression by a mechanism involving mitogen-activated protein kinase. J Biol
	Chem. 272: 10287-94.
	19. Sakamoto, T.M. et al. (2013) Altered red cell and platelet adhesion in hemolytic diseases:
	Hereditary spherocytosis, paroxysmal nocturnal hemoglobinuria and sickle cell disease. Clin
	<u>Biochem. pii: S0009-9120(13)00442-6.</u>
	20. Zucchetto, A. et al. (2012) The CD49d/CD29 complex is physically and functionally associated
	with CD38 in B-cell chronic lymphocytic leukemia cells. Leukemia. 26 (6): 1301-12.
Storage	Store at +4°C or at -20°C if preferred.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted.
Storage	
Storage	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.
Storage	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be
Storage Guarantee	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product
Guarantee	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. 18 months from date of despatch.
Guarantee Health And Safety	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. 18 months from date of despatch. Material Safety Datasheet documentation #10041 available at:
Guarantee	This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. 18 months from date of despatch.

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL:FITC (MCA928F)

Recommended Useful Reagents

HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)

America

North & South Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700 Europe Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com

Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com

'M343349:190110'

Printed on 11 Oct 2019

© 2019 Bio-Rad Laboratories Inc | Legal | Imprint