

Datasheet: MCA929A647

BATCH NUMBER 1014

| Description: | MOUSE IgG2a NEGATIVE CONTROL:Alexa Fluor® 647 |
|---------------|---|
| Specificity: | MOUSE IgG2a NEGATIVE CONTROL |
| Format: | ALEXA FLUOR® 647 |
| Product Type: | Negative/Isotype Control |
| Isotype: | lgG2a |
| Quantity: | 100 TESTS/1ml |
| | |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | • | | | * |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures.

^{*} It is recommended that the user dilutes the antibody for use in their own system to a concentration equivalent to their test reagents.

| Target Species | Negative Control | |
|-----------------------------------|---|--|
| Product Form | Purified IgG conjugated to Alexa Fluor® 647 - liquid | |
| Preparation | Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant | |
| Buffer Solution | Phosphate buffered saline | |
| Preservative Stabilisers | 0.09% Sodium Azide 1% Bovine Serum Albumin | |
| Approx. Protein Concentrations | IgG concentration 0.05 mg/ml | |
| Immunogen | Activated rat T-helper cells. | |
| Fusion Partners | Spleen cells from immunised BALB/c mice were fused with cells of the NS1 mouse | |

myeloma cell line.

Specificity

Mouse IgG2a negative control antibody, clone OX-34 is suitable for use as a negative control reagent for the measurement of non-specific binding of mouse monoclonal antibodies of isotype IgG2a to human tissue.

Clone MRC OX-34 recognises a rat cell surface marker, and therefore cannot be used as a negative control in this species.

This product is routinely tested in flow cytometry on rat splenocytes to confirm antibody activity and on human whole blood to test for suitability as a negative control.

Test results have shown that MCA929 is also suitable for use as a negative control with bovine, ovine, porcine, equine, canine, lapine and guinea-pig tissues.

This antibody may not be suitable for intracellular staining on some cell types.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10^6 cells in 100ul.

References

- 1. Avigdor, A. *et al.* (2004) CD44 and hyaluronic acid cooperate with SDF-1 in the trafficking of human CD34+ stem/progenitor cells to bone marrow. <u>Blood. 103: 2981-2989.</u>
- 2. Kamble, N.M. *et al.* (2016) Interaction of a live attenuated *Salmonella gallinarum* vaccine candidate with chicken bone marrow-derived dendritic cells. <u>Avian Pathol. Jan</u> 26:1-24. [Epub ahead of print]

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee

18 months from date of despatch.

Acknowledgements

This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com

Health And Safety Information

Material Safety Datasheet documentation available at: https://www.bio-rad-antibodies.com/SDS/MCA929A647

Material Safety Datasheet Documentation #10041 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376 America Fax: +1 919 878 3751 **Worldwide** Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +1 919 878 3751 Fax: +44 (0)1865 852 739
Email: antibody_sales_us@bio-rad.com Email: antibody_sales_uk@bio-rad.com

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M281135:160211'

Printed on 27 Mar 2025

© 2025 Bio-Rad Laboratories Inc | Legal | Imprint