

Datasheet: MCA929P647

BATCH NUMBER 157436

| | |
|----------------------|---|
| Description: | MOUSE IgG2a NEGATIVE CONTROL:RPE-Alexa Fluor® 647 |
| Specificity: | MOUSE IgG2a NEGATIVE CONTROL |
| Format: | RPE-ALEXA FLUOR® 647 |
| Product Type: | Negative/Isotype Control |
| Isotype: | IgG2a |
| Quantity: | 100 TESTS/1ml |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------|-----|----|----------------|--------------------|
| Flow Cytometry | ■ | | | * |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. * It is recommended that the user dilutes the antibody for use in their own system to a concentration equivalent to their test reagent.

| | | | |
|------------------------|---|----------------------------|--------------------------|
| Target Species | Negative Control | | |
| Product Form | Purified IgG conjugated to R. Phycoerythrin (RPE) - Alexa Fluor® 647 - lyophilized | | |
| Reconstitution | Reconstitute with 1.0 ml distilled water | | |
| Max Ex/Em | Fluorophore | Excitation Max (nm) | Emission Max (nm) |
| | RPE-Alexa Fluor®647 488nm laser | 496 | 667 |
| | RPE-Alexa Fluor®647 561nm laser | 546 | 667 |
| Preparation | Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant | | |
| Buffer Solution | Phosphate buffered saline | | |
| Preservative | 0.09% Sodium Azide | | |
| Stabilisers | 1% Bovine Serum Albumin | | |

5% Sucrose

| | |
|------------------|-------------------------------|
| Immunogen | Activated rat T-helper cells. |
|------------------|-------------------------------|

| | |
|-------------|-----------|
| RRID | AB_567345 |
|-------------|-----------|

| | |
|------------------------|---|
| Fusion Partners | Spleen cells from immunised BALB/c mice were fused with cells of the NS1 mouse myeloma cell line. |
|------------------------|---|

| | |
|--------------------|--|
| Specificity | <p>Mouse IgG2a negative control antibody, clone OX-34 is suitable for use as a negative control reagent for the measurement of non-specific binding of mouse monoclonal antibodies of isotype IgG2a to human tissue.</p> <p>Clone MRC OX-34 recognises a rat cell surface marker, and therefore cannot be used as a negative control in this species.</p> <p>This product is routinely tested in flow cytometry on rat splenocytes to confirm antibody activity and on human whole blood to test for suitability as a negative control.</p> <p>Test results have shown that MCA929 is also suitable for use as a negative control with bovine, ovine, porcine, equine, canine, lapine and guinea-pig tissues.</p> <p>This antibody may not be suitable for intracellular staining on some cell types.</p> |
|--------------------|--|

| | |
|-----------------------|---|
| Flow Cytometry | Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul. |
|-----------------------|---|

| | |
|-------------------|---|
| References | <ol style="list-style-type: none">1. Avigdor, A. <i>et al.</i> (2004) CD44 and hyaluronic acid cooperate with SDF-1 in the trafficking of human CD34+ stem/progenitor cells to bone marrow. Blood. 103 (8): 2981-9.2. Kamble, N.M. <i>et al.</i> (2016) Interaction of a live attenuated <i>Salmonella gallinarum</i> vaccine candidate with chicken bone marrow-derived dendritic cells. Avian Pathol. Jan 26:1-24. [Epub ahead of print]3. Wattedegera, S.R. <i>et al.</i> (2017) Enhancing the toolbox to study IL-17A in cattle and sheep. Vet Res. 48 (1): 20. |
|-------------------|---|

| | |
|----------------|---|
| Storage | <p>Prior to reconstitution store at +4°C. Following reconstitution store at +4°C.</p> <p>DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p> |
|----------------|---|

| | |
|------------------|---------------------------------|
| Guarantee | 12 months from date of despatch |
|------------------|---------------------------------|

| | |
|-------------------------|---|
| Acknowledgements | The Alexa Fluor® dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, except for use in combination with microarrays, and are covered by pending and issued patents. |
|-------------------------|---|

Health And Safety Information Material Safety Datasheet documentation #20487 available at:
<https://www.bio-rad-antibodies.com/SDS/MCA929P647>
20487

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL:RPE-Alexa Fluor® 647 \(MCA1210P647\)](#)

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M378909:210311'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)