

Datasheet: OBT0009 BATCH NUMBER 168814

Description:	RABBIT ANTI HUMAN TdT KIT
Name:	TdT
Other names:	TERMINAL DEOXYNUCLEOTIDYL TRANSFERASE
Format:	Kit
Product Type:	Kits
Isotype:	Polyclonal IgG
Quantity:	20 TESTS

Product Details

Applications	This product has been re derived from testing within communications from the information. For general p rad-antibodies.com/protoc	n our labo originator protocol re	ratories, s. Please	peer-reviewed public refer to references i	ations or personal indicated for further
		Yes	No	Not Determined	Suggested Dilution
	Immunofluorescence	•			
	Where this product has n			•	•
	necessarily exclude its us a guide only. It is recomm system using appropriate	ended that	at the use	er titrates the product	•
Target Species	Human				
Product Form	Kit				
Buffer Solution	Phosphate buffered saline	9			
Preservative Stabilisers	All components contain 0 Components 2, 2C and 3 Component M contains 8	contain 1	% Bovine	(⁰)	
External Database Links	UniProt: <u>P04053</u> <u>Related re</u> Entrez Gene: 1791 DNTT Relate	eagents d reagents	S		

Synonyms	TDT
Product Information	Rabbit anti Human TdT kit is a set of reagents intended to allow the identification of Terminal deoxynucleotidyl transferase positive (TdT ⁺) cells by combined fluorescence and phase contrast microscopy.
Reagents In The Kit	Kit contains ready-to-use reagents sufficient for twenty slide tests:
	1 Normal goat serum (10% normal goat serum in PBS/azide) - 1.0ml
	2 Rabbit anti- human -TdT antibody mixture - 0.5ml
	2C Negative control antibody - 0.5ml
	3 FITC-goat anti- rabbit IgG F(ab') ₂ - 1.0ml.
	M Mounting medium (80% glycerol in PBS/azide) - 1.0ml
Additional Reagents Required	Cold absolute methanol Phosphate buffered saline (PBS)
Instructions For Use	Standard bone marrow or peripheral blood smears may be used for immunofluorescence staining. Cytocentrifuge preparations of cell suspensions made up in RPMI-1640 tissue culture medium containing 10% foetal calf serum are also suitable. Frozen sections should be cut as thin as possible. Store unfixed slide preparations in a dessicator at room temperature.
	1. Circle two areas of interest (1 cm diameter) with diamond pencil on the slide to be stained. One will be used as a control area and the other is for the test. Fix smears for 30 min. in absolute methanol at 4°C immediately before staining. Air dry.
	2. Block by adding a drop of normal goat serum to both circled areas. Let stand for 5 min. Rinse off goat serum with PBS.
	3. Apply 25 ul rabbit anti human TdT antibody mixture to the test area and 25 ul of negative control to control area and allow to react for 30 min in a humid chamber.
	4. Rinse avoiding cross-contamination of reagents. Now wash for 10 min. in a Coplin jar with several changes of PBS. Agitate occasionally for thorough washing.
	5. Wipe the slide, avoiding cells, and shake excess wash solution off circled area. Spread 25 ul of FITC goat anti rabbit IgG F(ab') ₂ and incubate in humid chamber for 30 min in the dark.
	6. Rinse off and wash for 10 min. in a Coplin jar with several changes of PBS.
	7. Apply small drops of glycerol/PBS (mounting medium) to two glass cover slips and invert circled areas on slide over dropped cover slips.

Storage		This produ	uct is shippe	ed at ambient temperatu	re. It is recor	nmended to aliquot and store
			•	en thawed, aliquot the s 4 weeks) and store the	•	eded. Keep aliquots at 2-8°C f quots at -20°C.
		•		ng and thawing as this r ot recommended.	nay denature	e the antibody. Storage in
Guarante	e	Guarantee	ed until date	of expiry. Please see p	roduct label.	
Health A	nd Safety	Material S	afety Datasl	heet documentation #10)504 #10147	#10143 available at:
1						
Informati	on	https://ww		ntibodies.com/SDS/OBT	0009	
Informati		<u>https://ww</u> 1 - Norma	l goat serum	ntibodies.com/SDS/OBT n (10504)	0009	
Informati		<u>https://ww</u> 1 - Norma M - Mount	l goat serum ting Medium	ntibodies.com/SDS/OBT n (10504)	0009	
Informati Regulato		https://ww 1 - Norma M - Mount 2, 2C & 3	l goat serum ting Medium	ntibodies.com/SDS/OBT n (10504) (10147) Preparations (10143)	0009	
	ry	https://ww 1 - Norma M - Mount 2, 2C & 3 For reseat	l goat serum ting Medium - Antibody F	ntibodies.com/SDS/OBT n (10504) (10147) Preparations (10143)	0009 Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50

Printed on 29 Aug 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint