

Datasheet: PUR006

BATCH NUMBER 121110

Description:	PROTEUS PROTEIN A MINI PURIFICATION STARTER KIT
Name:	ANTIBODY PURIFICATION KIT PROTEIN A
Format:	MINI - Spin Columns, Buffers, U/F Spinners
Product Type:	Purification Kit
Quantity:	2 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers	Buffers contain 0.09% Sodium Azide (NaN ₃) Resin contains 0.01% Thiomersal
---------------------------------	---

Product Information	See Binding affinities for Protein A and G
----------------------------	--

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit	2 x Protein A Spin Column Plugs 2 x Spin Columns 4 x 2.2 ml Centrifuge Tubes 2 x 10 kDa MWCO Ultrafiltration Spinners 1 x 60 ml of 1 x Buffer A 1 x 20 ml of 1 x Buffer B1 1 x 20 ml of 1 x Buffer B2 1 x 10 ml of 1 x Buffer C Plug Insertion Tool Instructions
----------------------------	---

References	1. Raeisi, H. <i>et al.</i> (2019) Applying the pthA effector protein of <i>Xanthomonas citri</i> . subsp. <i>citri</i> . for production of specific antibodies and its application for detection of infected plants J Plant Pathol. 102 (1): 79-87.
Storage	Store at +4°C. DO NOT FREEZE.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10251 #10253 #10254 #10252 available at: https://www.bio-rad-antibodies.com/SDS/PUR006 Protein A Binding Buffer A (10251) Elution Buffer B2 (10253) Neutralisation Buffer C (10254) Elution Buffer B1 (10252)
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN G MINI PURIFICATION STARTER KIT \(PUR014\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350457:190307'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)