

Datasheet: PUR011 BATCH NUMBER 230507

| Description: PROTEUS PROTEIN G MIDI PURIFICATION SPIN COLUMN PACE Name: ANTIBODY PURIFICATION KIT PROTEIN G Format: MIDI - Spin Columns Product Type: Purification Kit Quantity: 12 UNITS | | |
|---|----------------------|--|
| Format: MIDI - Spin Columns Product Type: Purification Kit | Description: | PROTEUS PROTEIN G MIDI PURIFICATION SPIN COLUMN PACK |
| Product Type: Purification Kit | Name: | ANTIBODY PURIFICATION KIT PROTEIN G |
| | Format: | MIDI - Spin Columns |
| Quantity: 12 UNITS | Product Type: | Purification Kit |
| | Quantity: | 12 UNITS |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------|-----|----|-----------------------|---------------------------|
| Protein Purification | | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| Preservative |
|--------------|
| Stabilisers |

20% Ethanol

Product Information

See Binding affinities for Protein A and G

Recommended Reading - Proteus Protein A Purification Handbook

For more information visit Proteus: Protein Purification Spin Kits

Reagents In The Kit

12 x Protein G Spin Column Plugs

12 x Spin Columns

24 x 50 ml Centrifuge Tubes

Plug Insertion Tool

Instructions

References

1. Montaudon, E. *et al.* (2015) Cardiac effects of long-term active immunization with the second extracellular loop of human β 1- and/or β 3-adrenoceptors in Lewis rats. Pharmacol

Res. 100: 210-9.

Storage

Store at +4°C. DO NOT FREEZE

| Guarantee | Guaranteed until date of expiry. Please see product label. |
|----------------------------------|---|
| Health And Safety Information | Material Safety Datasheet documentation #qtec44-5 available at: https://www.bio-rad-antibodies.com/SDS/PUR011 qtec44-5 Not required |
| Regulatory | For research purposes only |

Related Products

Recommended Useful Reagents

MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT (MMT1)
RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT (RMT1)

MYCOPLASMA REMOVAL AGENT (BUF035)

PROTEUS PROTEIN A MIDI PURIFICATION SPIN COLUMN PACK (PUR002)

 North & South
 Tel: +1 800 265 7376
 Worldwide
 Tel: +44 (0)1865 852 700
 Europe
 Tel: +49 (0) 89 8090 95 21

 America
 Fax: +1 919 878 3751
 Fax: +44 (0)1865 852 739
 Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets 'M350460:190307'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | Legal | Imprint