

Datasheet: PUR012

BATCH NUMBER 204901

| | |
|----------------------|--|
| Description: | PROTEUS PROTEIN G MIDI PURIFICATION KIT |
| Name: | ANTIBODY PURIFICATION KIT PROTEIN G |
| Format: | MIDI - Spin Columns, Buffers, U/F Spinners |
| Product Type: | Purification Kit |
| Quantity: | 4 UNITS |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------|-----|----|----------------|--------------------|
| Protein Purification | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Preservative Stabilisers | Buffers contain 0.09% Sodium Azide (NaN ₃) Resin contains 20% Ethanol |
|---------------------------------|--|

Product Information See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

| | |
|----------------------------|--|
| Reagents In The Kit | 4 x Protein G Spin Column Plugs 4 x Spin Columns 8 x 50 ml Centrifuge Tubes 4 x 10 kDa MWCO Ultrafiltration Spinners 1 x 250 ml of 1 x Buffer G 1 x 125 ml of 1 x Buffer B2 1 x 30 ml of 1 x Buffer C Plug Insertion Tool Instructions |
|----------------------------|--|

References 1. Abdelkrim, M.A. *et al.* (2012) Autoantibodies against cardiac $\beta(1)$ -adrenoceptor do not

affect the low-affinity state $\beta(1)$ -adrenoceptor-mediated inotropy in rat cardiomyocytes.

[Can J Physiol Pharmacol. 90: 407-14.](#)

2. Soler, L. *et al.* (2011) Development and validation of a novel competitive ELISA for the detection of serum amyloid A in pigs. [Vet J. 190: e7-11.](#)

3. Kimura, K. *et al.* (2012) Pyrroloquinoline quinone stimulates epithelial cell proliferation by activating epidermal growth factor receptor through redox cycling. [Free Radic Biol Med. 53: 1239-51.](#)

4. Montaudon, E. *et al.* (2014) Cardiovascular effects of beta1 and beta3-adrenergic receptor autoantibodies in Lewis rat.

[Article in French] [Ann Cardiol Angeiol \(Paris\). 63: 128-34.](#)

5. Yuan, W. *et al.* (2015) Quantitative analysis of immunoglobulin subclasses and subclass specific glycosylation by LC-MS-MRM in liver disease. [J Proteomics. 116: 24-33.](#)

6. Broecker, F. *et al.* (2016) Multivalent display of minimal *Clostridium difficile* glycan epitopes mimics antigenic properties of larger glycans. [Nat Commun. 7: 11224.](#)

7. Montaudon E *et al.* (2015) Cardiac effects of long-term active immunization with the second extracellular loop of human $\beta 1$ - and/or $\beta 3$ -adrenoceptors in Lewis rats. [Pharmacol Res. 100: 210-9.](#)

Storage

Store at +4°C. DO NOT FREEZE.

Guarantee

Guaranteed until date of expiry. Please see product label.

Health And Safety Information

Material Safety Datasheet documentation #10125 #10253 #10254 available at: <https://www.bio-rad-antibodies.com/SDS/PUR012>

Protein G Binding Buffer G (10125)

Elution Buffer B2 (10253)

Neutralisation Buffer C (10254)

Regulatory

For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN A MIDI PURIFICATION KIT \(PUR003\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M350461:190307'

Printed on 25 Mar 2023