

## Datasheet: PUR032

**BATCH NUMBER 23311**

<b>Description:</b>	PROTEUS MINI ENDOTOXIN REMOVAL SPIN COLUMN PACK
<b>Name:</b>	ENDOTOXIN REMOVAL MINI KIT
<b>Format:</b>	MINI - Spin Columns
<b>Product Type:</b>	Purification Kit
<b>Quantity:</b>	12 UNITS

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures.

### Preservative Stabilisers

94% Preservative base of resin containing:  
20% Aqueous Ethanol  
80% 0.1M Sodium Chloride

### Endotoxin Level

This kit can be used for samples with endotoxin loads < 3,000 EU. For endotoxin loads <30,000 EU both the Midi Standard and High Capacity Endotoxin Removal kits can be used. Endotoxin loads up to 1,000,000 EU require the use of a Midi High Capacity Endotoxin Removal kit.

### Product Information

This kit offers rapid, one step, high grade (>95%) removal of endotoxin from a wide range of protein solutions, with >90% protein recovery.

We recommend reading the [Proteus Endotoxin Removal Kit Handbook](#) prior to use.

For more information visit [Proteus: Protein Purification Spin Kits](#).

### Reagents In The Kit

12 x vials containing 0.25ml Endotoxin Binding Resin  
12 x Mini Spin Column barrels (20ml capacity)  
12 x 50ml Centrifuge Tubes  
12 x Clear Spin Push Caps  
12 x Yellow Incubation Screw Caps  
Instructions for use

<b>Additional Reagents Required</b>	0.2µm syringe filter units Low endotoxin pre-equilibration buffer e.g PBS pH7.4
<b>Instructions For Use</b>	Instructions for use can be found at <a href="http://www.bio-rad-antibodies.com/uploads/IFU/PUR032.pdf">www.bio-rad-antibodies.com/uploads/IFU/PUR032.pdf</a>
<b>References</b>	1. Negus D <i>et al.</i> (2015) Parenteral Administration of Capsule Depolymerase EnvD Prevents Lethal Inhalation Anthrax Infection. <a href="#">Antimicrob Agents Chemother. pii: AAC.01547-15.</a>
<b>Storage</b>	Upon receipt of kit, remove the vials of resin and store them at +4°C. DO NOT FREEZE RESIN VIALS OR STORE AT ROOM TEMPERATURE. Store remainder of kit at room temperature.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #20510 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PUR032">https://www.bio-rad-antibodies.com/SDS/PUR032</a> 20510
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)  
'M417827:230417'

Printed on 18 Jan 2024