

Datasheet: STAR129F

BATCH NUMBER 164543

| | |
|----------------------|---|
| Description: | GOAT ANTI HUMAN LAMBDA LIGHT CHAIN:FITC |
| Specificity: | LAMBDA LIGHT CHAIN |
| Format: | FITC |
| Product Type: | Polyclonal Antibody |
| Isotype: | Polyclonal IgG |
| Quantity: | 1 mg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry | ▪ | | | Neat - 1/10 |
| Immunohistology - Frozen | ▪ | | | |
| Immunohistology - Paraffin | | | ▪ | |
| Immunofluorescence | ▪ | | | Neat - 1/10 |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

| Target Species | Human | | | | | | |
|-----------------------|--|-------------------|---------------------|-------------------|------|-----|-----|
| Product Form | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid | | | | | | |
| Max Ex/Em | <table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>FITC</td> <td>490</td> <td>525</td> </tr> </tbody> </table> | Fluorophore | Excitation Max (nm) | Emission Max (nm) | FITC | 490 | 525 |
| Fluorophore | Excitation Max (nm) | Emission Max (nm) | | | | | |
| FITC | 490 | 525 | | | | | |

Antiserum Preparation Antisera to human lambda chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

| | |
|---------------------------------|--|
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | <0.1% Sodium Azide (NaN ₃) |
| Approx. Protein | IgG concentration 1.0mg/ml |

Concentrations

Immunogen Human lambda light chains.

External Database Links

UniProt:

[P0CG04](#) [Related reagents](#)

Entrez Gene:

[28815](#) IGLV2-14 [Related reagents](#)

RRID AB_1102725

Specificity **Goat anti Human Lambda Light Chain antibody** recognizes human lambda light chains. This antibody has been cross absorbed with pooled human myeloma proteins bearing kappa light chains.

Flow Cytometry Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Ben-horin, S. *et al.* (2016) Cross-immunogenicity: antibodies to infliximab in Remicade-treated patients with IBD similarly recognise the biosimilar Remsima. [Gut. 65 \(7\): 1132-8.](#)
2. Magro, F. *et al.* (2018) The performance of Remicade[®]- optimized quantification assays in the assessment of Flixabi[®] levels. [Therap Adv Gastroenterol. 11: 1756284818796956.](#)

Storage This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light.

Guarantee Guaranteed for 12 months from the date of despatch or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: <https://www.bio-rad-antibodies.com/SDS/STAR129F10040>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M402052:220718'

Printed on 01 Mar 2024
